



EuroElone®

Genetic Sequencer

DNBSEQ-G400*

To be used with HotMPS exclusively



Welcome HotMPS**: Bringing choice to sequencing

Strengthen your daily sequencing capabilities with the unique combination of our new sequencing chemistry and the proven performance of our DNBSEQ-G400*.

Leading Life Science Innovation

© About MGI

MGI Tech Co., Ltd. (referred to as MGI) is committed to building core tools and technologies to lead life science through intelligent innovation. MGI focuses on R&D, production and sales of DNA sequencing instruments, reagents, and related products to support life science research, agriculture, precision medicine and healthcare. MGI is a leading producer of clinical high-throughput gene sequencers, and its multi-omics platforms include genetic sequencing, medical imaging, and laboratory automation.

MGI, founded in 2016 in Shenzhen, China, is one of the few companies worldwide that can develop and mass-produce clinical genetic sequencers. Providing real-time, comprehensive, life-long solutions, its vision is to lead life science innovation.

At present, MGI has more than 1,800 employees and operates in more than 70 countries and regions serving more than 1,100 customers around the world, with 8 international subsidiaries including Wuhan, Qingdao, Changchun, Kunshan and Hongkong, China, Japan, UAE, Latvia, and USA.

© Disclaimer

DNBSEQ-G400RS (Cat.No.900-000493-00) and DNBSEQ-G400 (Cat.No.900-000492-00) whose software has been configured for HotMPS MUST be used in conjunction with MGI's HotMPS sequencing reagent, and MUST NOT be used with MGI's CoolMPS or StandardMPS reagents (or with any reagents containing a 3'O-azidomethyl blocking group). This is important because (i) such sequencers will not work effectively with MGI's CoolMPS or StandardMPS reagents (or with any reagents containing a 3'O-azidomethyl blocking group); and (ii) the use of such reagents will give rise to the risk of patent infringement proceedings. Use of MGI's CoolMPS or StandardMPS reagents (or with any reagents containing a 3'O-azidomethyl blocking group) with such sequencers will invalidate any warranty which may have been provided by MGI, and any liability for intellectual property infringement arising from the use of such reagents is excluded from any IP indemnity.

The HotMPS sequencing reagent can only be used with DNBSEQ-G400RS (Cat.No.900-000493-00) and DNBSEQ-G400 (Cat.No.900-000492-00) whose software has been configured for HotMPS, and will not work effectively with MGI's sequencer whose software has not been properly configured for HotMPS. Use of these reagents with sequencers whose software has not been configured for HotMPS will invalidate any warranty which may have been provided by MGI, and any liability for intellectual property infringement arising from the use of such reagents is excluded from any IP indemnity.

© About

DNBSEQ-G400*

01 Product Introduction

Technical Principle

Hardware Platform

Performance Parameters

Adapted Applications

02 Solutions

Versatile library preparation, sequencing & analysis solution

Flexible sequencing & analysis solution

Fluent library preparation & sequencing solution

03 Applications

WGS

WES

RNA-Seq

10X Genomics 3'RNA-Seq

Pathogen Sequencing – COVID-19

04 Appendix

Hardware Parameters

Ordering Information


MGI Global Presence

Flexible, high quality
Activate your daily sequencing capability

Product Introduction

DNBSEQ-G400* is a versatile benchtop sequencer, providing users with comprehensive, flexible, and efficient sequencing options.

DNBSEQ-G400* sequencer supports a wide range of applications including scientific research, clinical research, disease prevention, environment studies and agriculture, etc., increasing the popularity of high-throughput sequencing systems in medical and scientific research fields.



Multiple choices

- **Dual flow cell system**

One or two flow cells:
1500M-3600M reads/run

High quality

- **Proven DNBSEQ™ technology**

No error accumulation
Low amplification bias
Low index hopping

High speed

- **Short sequencing & time**

FCL SE50 from DNB to
FASTQ takes only 14 hours


Versatile

- **Wide range of applications**

Suitable for research and clinical applications

© Technical Principle

MGI'S PROPRIETARY 「DNBSEQ™」 TECHNOLOGY




- ↑ INCREASED ACCURACY
- ↓ DECREASED DUPLICATES
- ↓ REDUCED INDEX HOPPING

DNA Nanoball sequencing technology - No accumulation of amplification errors

© Hardware Platform

DNBSEQ-G400* sequencer utilizes an innovative flow cell system which can support various sequencing modes and an optimized optical and biochemical system that enables the whole sequencing process to be completed within a short period of time, offering the user a simplified and streamlined sequencing experience.



© Performance Parameters

| Effective Reads/Flow Cell ⁽¹⁾ | Reads Lengths | Data Output/Flow Cell ¹ | Run Time ² | Q30 ³ |
|--|---------------|------------------------------------|-----------------------|------------------|
| 1500-1800M | SE50 | 75-90G | 15.5H | >90% |
| | SE100 | 150-180G | 25.5H | >85% |
| | PE50 | 150-180G | 28.5H | >90% |
| | PE100 | 300-360G | 50.5H | >85% |

1. The maximum number of effective reads are based on the sequencing of an internal standard library. Actual output may vary depending on sample type and library preparation method.

2. Run time was calculated based on dual-slides mode and takes sample loading, sequencing, base calling and data processing in account.

3. The percentage of base above Q30 is the average of an internal standard library over the entire run. Actual performance is affected by factors such as sample type, library quality, and insert fragment length.

© Adapted Applications

DNBSEQ-G400* supports four read lengths, from fertility testing and rapid pathogen detection (SE50, SE100), to tumor detection, transcriptome, WGBS, WES, WGS (PE100), comprehensively cover the application needs of scientific research and clinical fields.



| Application type | Recommended read length | Recommended data per sample | Recommended sample numbers ⁴ for a single run on DNBSEQ-G400* | | Recommended max sample numbers for annual run on DNBSEQ-G400* |
|------------------------------|-------------------------|--|--|-------------|---|
| | | | 1 FCL | 2 FCL | 2 FCL |
| | | | 1800M reads | 3600M reads | 3600M reads |
| NIPT | SE50 | 5M reads/sample | 275 samples | 550 samples | 137,500 samples |
| Pathogen Fast Identification | SE50/SE100 | 25M reads/sample | 55 samples | 110 samples | 18,810/27,500 samples |
| Single cell RNA-Seq | PE100 | 5000cells, 100K reads/-cell, 100 Gb/sample | 4 samples | 8 samples | 800 samples |
| Cancer small panel | | 1 Gb/sample | 275 samples | 550 samples | 550,000 samples |
| Cancer large panel | | 5 Gb/sample | 55 samples | 110 samples | 11,000 samples |
| Transcriptome | | 40M reads/sample | 35 samples | 70 samples | 7,000samples |
| WES | | 15 Gb/sample (100×) | 20 samples | 40 samples | 4,000 samples |
| WGS | | 90 Gb/sample (30×) | 4 samples | 8 samples | 800 samples |

4. Theoretical sample numbers are calculated with consideration to pooling variation and applications. For reference only.

© Total Package

Versatile Library Prep, Sequencing & Analysis Solution

Genetic Sequencer DNBSEQ-G400* fully automated workflow (& all scenarios applicable)



© Highlights

Efficient

- High efficiency with fully automated comprehensive workflow from sample to report
- Less than 30 mins of manual operations from sample to sequencing

Intelligent

- ZLIMS controls fully automated workflow from sample to report
- Compatible with multiple MGI and third-party software
- MegaBOLT is compatible with MGI and third-party data

Flexible

- Multiple sequencing and analysis options to expand application scenarios
- The above equipments can be selected according to needs

Applications

Whole Genome Sequencing (WGS)

Case 1

Background

Whole genome sequencing (WGS) is designed to fully decode the genetic information of individuals and populations. With the increasing number of species with known genomic sequences, WGS is widely used in animal and plant breeding, population studies, disease research, clinical diagnosis and drug development, and has become one of the most rapid and effective sequencing methods.

In this case study the DNBSEQ-G400* was used together with the HotMPS High-throughput Sequencing Set** for whole genome sequencing.



Figure 1-1 Whole Genome Sequencing (WGS)

Samples: Human Cell Line NA12878

Library kit: MGIEasy PCR-Free DNA Library Prep Set¹, MGIEasy FS PCR-Free DNA Library Prep Kit¹

Sequencing Parameter: Genetic Sequencer DNBSEQ-G400RS* with HotMPS High-throughput Sequencing Set (G400 HM FCL PE100)**

Result

- The Q30 ratios of the bases in all samples are over 93%.
- Low duplication rate (<3%) and high unique mapping rate (>94%) indicate higher data efficiency
- The coverage in high GC content region is even




Figure 1-2 WGS GC distribution curve

Conclusion

DNBSEQ™ Technology can generate accurate and reliable whole-genome sequencing data, with high Q30, high unique mapping rate, and low duplication rate.




Figure 1-3 WGS coverage

1. MGIEasy FS PCR-Free DNA Library Prep Kit uses enzymatic fragmentation while MGIEasy PCR-Free DNA Library Prep Set uses mechanical fragmentation.

Table 1 WGS 30X data quality and analysis

| Parameters | MGI FS PCR-Free | MGI PCR-Free | Competitor N PCR-Free |
|----------------------------------|-----------------|--------------|-----------------------|
| Clean Q30 (%) | 93.04% | 95.73% | 92.19% |
| GC Content (%) | 41.15% | 41.60% | 41.56% |
| Mapping Rate (%) | 98.94% | 98.92% | 99.29% |
| Unique mapping Rate (%) | 94.67% | 94.45% | 87.87% |
| Duplicate Rate (%) | 2.08% | 2.24% | 9.46% |
| Mismatch Rate (%) | 0.48% | 0.47% | 0.76% |
| Insert size | 3679 | 281.2 | 455.5 |
| Average sequencing depth (rmdup) | 30.12 | 30.62 | 30.52 |
| Uniformity(>0.2f) | 98.07% | 98.02% | 98.17% |




Figure 1-4 SNP' s and Indel' s detection results

Whole Exome Sequencing (WES)

Case 2

Background

With the development of precision medicine, Whole Exome Sequencing (WES) offers advantages for analysing the data of the patient's cases where the diagnoses or treatment options are unclear. For example, for many cancer types, Whole Exome Sequencing (WES) has been able to place specific types on treatments that in some cases have led to a complete cure.

MGI's Genetic Sequencer DNBSEQ-G400RS* utilizes an innovative system that can support WES sequencing

Samples: Cell lines YH and NA12878

Library kit: MGIEasy Exome Capture V4 Probe Set, MGIEasy Exome Universal Library Prep Set

Sequencing Parameter: Genetic Sequencer DNBSEQ-G400RS* with HotMPS High-throughput Sequencing Set(G400 HM FCL PE100)**

Analysis Pipeline: MGI MegaBOLT (V2.2.3.0)

Result

- Raw depth of more than 96% of the samples is above 20X;
- Average uniformity (>0.2f) is over than 92%;
- High data quality of WES ensures accurate SNP/InDel calling(>90%).

Conclusion:

High capture rate, uniformity, and good SNP/InDel calling results show DNBSEQ-G400' can perform well in WES application.

Table 2 WES 200X data quality and analysis

| | YH-1 | YH-2 | NA12878 |
|-----------------------|--------|--------|---------|
| Capture_Rate_on_Reads | 68.00% | 69.54% | 65.89% |
| Capture_Rate_on_Bases | 57.47% | 58.98% | 55.63% |
| Average_depth(rmdup) | 205.52 | 206.04 | 208.61 |
| Coverage(>=1X) | 99.27% | 99.30% | 99.69% |
| Coverage(>=20X) | 96.36% | 96.30% | 97.32% |
| Coverage(>=100X) | 76.01% | 74.58% | 74.72% |
| Uniformity(>0.2f) | 93.03% | 92.47% | 92.85% |

Table 2a WES 200X SNP data quality and analysis

| | NA12878 |
|-------------------|---------|
| Precision_SNP | 99.64% |
| Sensitivity_SNP | 98.62% |
| F-measure_SNP | 99.12% |
| Precision_InDel | 91.91% |
| Sensitivity_InDel | 90.42% |
| F-measure_InDel | 91.16% |

RNA sequencing

Case 3

Background

RNA sequencing is a powerful method for comprehensive and rapid analysis of gene expression changes, examining rare and novel transcripts, discovery of alternate splicing events, gene fusions, SNPs and allele-specific expression in tissues or cells. Researchers aim to understand known or novel factors that alter gene expression during particular biological processes or development of diseases. Therefore, RNA-Seq has been widely applied to many fields such as biological research, disease research, drug development, agriculture and environmental research.


Samples: Universal Human Reference RNA (UHRR)

Library kit: MGIEasy RNA Library Prep Set

Sequencing Parameter: Genetic Sequencer DNBSEQ-G400RS* with HotMPS High-throughput Sequencing Set (G400 HM FCL PE100)**

Result

- Q30 ratios of Clean Reads in all samples are over 96%.
- Average mapping genome ratios is over 93%;
- RNA-seq and qPCR correlation coefficient indicate good correlation (pearson $r > 0.87$).



Conclusion

The high accuracy of DNBSEQ-G400* allows researchers to confirm and examine the transcriptome on an unprecedented scale.

Table 3 RNA-Seq data quality

| | RNA-UHRR_1 | RNA-UHRR_2 | RNA-UHRR_3 | RNA-UHRR_4 | RNA-UHRR_5 |
|-------------------------|------------|------------|------------|------------|------------|
| Read_length | 100:100 | 100:100 | 100:100 | 100:100 | 100:100 |
| Total Raw Reads(M) | 98.75 | 115.03 | 118.7 | 88.58 | 111.49 |
| Total Clean Reads(M) | 97.66 | 113.74 | 117.42 | 87.53 | 110.28 |
| Clean Reads Ratio(%) | 98.89% | 98.88% | 98.92% | 98.82% | 98.91% |
| Clean Reads Q30(%) | 96.86% | 96.76% | 96.65% | 96.64% | 96.52% |
| rRNA Ratio(%) | 3.25% | 3.00% | 3.06% | 2.93% | 3.28% |
| Mapping Genome Ratio(%) | 94.05% | 93.98% | 93.93% | 94.00% | 93.93% |
| Mapping Gene Ratio(%) | 52.98% | 52.79% | 52.46% | 52.96% | 53.42% |
| Gene Number | 21102 | 21184 | 21219 | 20929 | 21121 |
| Transcript Number | 40918 | 41342 | 41415 | 40501 | 41209 |
| Pearson_R | 87.20% | 87.00% | 86.90% | 87.00% | 87.50% |
| Spearman_R | 88.20% | 88.30% | 88.30% | 88.10% | 88.10% |

10X Genomics Single Cell 3' Gene Expression

Case 4

Background

The development of single-cell omics technology provides a powerful tool for the understanding of systems biology. Application of DNBSEQ-G400* with HotMPS** in 10x Genomics single-cell sequencing uses MGI library conversion technology to cooperate with 10x Genomics' Chromium system to achieve perfect compatibility between 10x Genomics single-cell sequencing technology and MGI's DNBSEQ platform*. MGI's sequencing products are compatible with the Chromium system while retaining its unique sequencing technology advantages.

Samples: Animal tissue

Library kit: Chromium Next GEM Single Cell 3' Reagent Kits v3.1 from 10X Genomics

Sequencing Parameter: Genetic Sequencer

DNBSEQ-G400RS* with HotMPS High-throughput Sequencing Set (G400 HM FCL PE100)**




Figure 4-1 % of Reads mapped to different regions

Result

- Higher data utilization (Number of reads on DNBSEQ-G400** is less than competitor but estimated number of cells detected are quite similar).
- Genome mapping rate is higher than competitor.

Conclusion

The single-cell RNA sequencing data on the DNBSEQ-G400* with the HotMPS** sequencing demonstrate the openness and flexibility of MGI's DNBSEQ sequencing platform, providing single-cell researchers with more choices.

Table 4 10X Genomics 3' RNA-Seq data quality and analysis

| Parameters | DNBSEQ-G400* | Competitor N |
|-----------------------------|--------------|--------------|
| Number of Reads | 494,078,254 | 719,160,550 |
| Valid Barcodes (%) | 97.70% | 98.30% |
| Valid UMIs (%) | 99.90% | 100% |
| Sequencing Saturation (%) | 90.00% | 92.80% |
| Q30 Bases in Barcode (%) | 92.80% | 97.50% |
| Q30 Bases in RNA Read (%) | 93.30% | 93.90% |
| Q30 Bases in UMI (%) | 92.80% | 97.40% |
| Estimated Number of Cells | 4,687 | 4,714 |
| Fraction Reads in Cells (%) | 91.20% | 91.20% |
| Mean Reads per Cell | 105,415 | 152,558 |
| Median Genes per Cell | 1,104 | 1,123 |
| Total Genes Detected | 16,722 | 16,309 |
| Median UMI Counts per Cell | 1,959 | 2,042 |

Pathogen Sequencing - COVID-19

Case 5

Background

The identification and tracking of viral variants and lineages have rapidly become a priority in the fight against SARS-CoV-2. MGI ATOplex RNA Library Prep Set is a 2-step multiplex PCR-based library preparation set, which provides a streamlined workflow for SARS-CoV-2 whole genome enrichment and amplification. Combined with DNBSEQ-based high-throughput sequencing platform, it can obtain the full-length genome sequences of SARS-CoV-2 and achieve relative quantification of SARS-CoV-2 for population-scale virus detection, surveillance and tracing.

Samples: SARS-CoV-2 RNA QC sample from GeneWell with 4000 copy/ μ L (4000cp) and 400copy / μ L (400cp), NA12878, water (25% Balanced library)

Library kit: ATOplex RNA Library Prep. Set (V1.1)

Sequencing Parameter: Genetic Sequencer DNBSEQ-G400RS* with HotMPS High-throughput Sequencing Set (G400 HM FCL PE100)**

Result

- High quality data (Q30 > 96%, estimated error rate is lower than 0.1%)
- In all libraries, mapping rate can reach 99.9%
- In the positive samples, uniformity can reach 98.5%
- There is no difference between repeats
- All libraries were detected accurately(100X coverage>98%; Relative virus percentage (Pct%) >95% for 4000cp and >70% for 400cp)
- 100X coverage and relative virus percentage (Pct%) were highly consistent between repeats




Figure 5-1 Clean Q30 curve of 2 repeats

Conclusion

Empowered by its multiplex PCR technology, MGI's proprietary ATOplex ultra-high multiplex platform together with the DNBSEQ-G400* provides effective whole-genome sequencing, surveillance, and monitoring of SARS-CoV-2 variants

| | 4000cp | 400cp | NA12878 | Water | Total |
|----------------------------|--------|-------|---------|-------|-------|
| Number of pooled libraries | 32 | 32 | 30 | 2 | 96 |

Table 5-1 Pooled Libraries

| Samples | Total Reads (M) | Q20% | Q30% | EstErr% | SplitRate (%) |
|---------|-----------------|-------|-------|---------|---------------|
| Pool 1 | 439.19 | 99.14 | 96.61 | 0.08 | 90.31 |
| Pool 2 | 458.74 | 99.14 | 96.57 | 0.08 | 91.58 |

Table 5-2 SARS-CoV-2 sequencing data quality




Figure 4 The mapping rate and uniformity (upper), 100X coverage and Pct% (below, Pct% refer to relative virus percentage of different samples)

Appendix

Hardware Parameters

| Model ¹ | Model ¹ | Intended Market |
|---|--------------------------------|-------------------------------|
| Model ¹ | DNBSEQ-G400* | IVD |
| | DNBSEQ-G400RS* | RUO |
| Dimensions | 1086 mm(L)×756 mm(W)×710 mm(H) | |
| Net Weight | 200 kg | |
| Power | Type | 100-240 V, 50/60 Hz |
| | Rated Power | 1200 VA |
| Operating Environment Requirements ² | Temperature | 19°C-25°C |
| | Relative Humidity | 20% RH-80% RH, non-condensing |
| | Atmospheric Pressure | 70 kPa-106 kPa |
| | Waterproof Rating | IPX0 |
| Control Computer Configurations ³ | CPU | Intel Xeon E5 10Core 2 2.2GHz |
| | Internal Storage | 256 GB RAM |
| | HDD | 16 Tb |
| | SSD | 480 Gb |
| | Operating System | Windows 10 Enterprise |

1. Only for model classification
2. For indoor use only; The Flow Cells can be stored and transported at room temperature. No liquid medium is needed
3. Supporting the computer's configurations and system updates

Ordering Information

| Cat. No. | Product Name |
|---------------|---|
| 900-000493-00 | Genetic Sequencer DNBSEQ-G400RS* |
| 900-000492-00 | Genetic Sequencer (DNBSEQ-G400)* |
| 1000028237 | HotMPS High-throughput Sequencing Set (G400 HM FCL SE50)** |
| 940-000089-00 | HotMPS High-throughput Sequencing Set (G400 HM FCL SE100)** |
| 940-000091-00 | HotMPS High-throughput Sequencing Set (G400 HM FCL PE100)** |
| 940-000092-00 | HotMPS High-throughput Sequencing Set (G400 HM FCL PE50)** |
| 1000028240 | HotMPS High-throughput Sequencing Set (G400 HM FCL sRNA SE50)** |
| 940-000098-00 | Universal Sequencing Reaction Kit (G400 HM FCL PE100(CE))** |
| 940-000099-00 | Universal Sequencing Reaction Kit (G400 HM FCL PE50(CE))** |
| 940-000094-00 | Universal Sequencing Reaction Kit (G400 HM FCL SE100(CE))** |
| 1000028239 | Universal Sequencing Reaction Kit (G400 HM FCL SE50(CE))** |

* This sequencer is only available in selected countries, and its software has been specially configured to be used in conjunction with MGI's HotMPS sequencing reagents exclusively.

**This sequencing reagent is only available in selected countries.

✓ Technical Support Globally

The technical support team has a complete global coverage including technical services centers and multiple locations in major international regions to maximize customer satisfaction.



Multiple local technical support centers around the world provide timely and effective technical support and training.



Spare part centers in Shenzhen, Wuhan, Qingdao, Tianjin, Hong Kong (China); Brisbane (Australia); and Riga (Latvia), to ensure sufficient supply of parts for machine maintenance.



Online technical support is accessible worldwide, with a fully functioning call center (Toll-Free Hotline: 4000-966-988) (9:00-12:00,13:00-18:00, Beijing time, workday) and multi-language online training courses coming soon.

✓ Comprehensive Instrument Service and Warranty Plans Globally



Warehouses in Shenzhen, Wuhan, Qingdao, Tianjin, Hong Kong, Taipei (China), Bangkok (Thailand, Asia-Pacific); Brisbane (Australia, Oceania); Riga (Latvia, Europe); and San Jose (the USA, America) are established to ensure sufficient supply of maintenance parts for major regions.



Free installation and system verification service (including the QC reagents and consumables) are provided to turn your investment into production quickly.



MGI is responsible for any manufacturing defects or faults on the system within the warranty. Warranty covers labor, parts and travel charges.



One Free instrument preventive maintenance provided with warranty, along with a variety of available extended warranty support plans.



Euroclone SpA Società a Socio Unico
Via Figino, 20/22 - 20016 Pero (MI) -
+39 02.381951 - +39 02.38101465 -
info@euroclone.it - www.euroclone.it
Quality Management System Certified
to ISO 9001 and ISO 13485
international standard