



## IMMUNOASSAY KITS



Bioprocess  
Cancer  
Cardiovascular  
Cell Death  
Cell Signaling  
Cyclic Nucleotides  
Cytokines  
Eicosanoids  
Endocrinology/Hormones

Epigenetics  
Immunology  
Immunity/Inflammation  
Metabolism  
Nephrology  
Neuroscience  
Oxidative Stress  
Proteostasis/Chaperones



# CONTENTS

|   |           |
|---|-----------|
| <b>Basic Immunoassay Formats</b>          | <b>5</b>  |
| <b>ELISAs By Target and Research Area</b> | <b>7</b>  |
| <b>ELISA Amplification Technology</b>     | <b>12</b> |
| <b>Research Areas</b>                     |           |
| Bioprocess Optimization                   | 14        |
| Cancer                                    | 16        |
| Cardiovascular                            | 18        |
| Cell Death                                | 20        |
| Cell Signaling                            | 21        |
| Cyclic Nucleotides                        | 22        |
| Endocrinology/Hormones                    | 24        |
| Epigenetics                               | 26        |
| Immunity/Inflammation Signaling           | 27        |
| Cytokines                                 | 29        |
| Eicosanoids                               | 31        |
| Metabolism                                | 33        |
| Nephrology                                | 35        |
| Neuroscience                              | 37        |
| Oxidative Stress                          | 40        |
| Proteostasis/Chaperones                   | 41        |
| <b>Good Practices</b>                     | <b>44</b> |
| <b>Troubleshooting</b>                    | <b>45</b> |



# ENTRUST YOUR PRECIOUS SAMPLES TO US

Assay kits can be convenient time-saving tools to help accelerate your research, but not all assays are created equal. Some that promise convenience fail to deliver biologically relevant sensitivity, or the reproducibility needed for long-term studies. This can cost you more time and money in the end.

Whether you use ELISA kits for biomarker detection, assessing cellular function with fluorescent probes, or screening enzyme modulators with activity assays, every Enzo assay kit includes something a lot of other companies out there today cannot match — Experience.


*Detect Small and Large Analytes Accurately and Efficiently*

It takes more than just an antibody pair to make an ELISA. Building an immunoassay requires screening multiple antibodies, selection of appropriate standards and conjugates, establishing proper sample preparation protocols, and validation of the assay in relevant matrices. Our expertise in developing ELISAs is further backed by four decades of manufacturing excellence. Strict validation criteria and state-of-the-art manufacturing facilities deliver reproducible assays that continue to be cited in peer-reviewed publications by scientists around the world.

Enzo Life Sciences offers hundreds of ELISA kits in both immunometric and competitive assay formats. As scientists and manufacturers of kits, we understand the critical nature of your research. Each kit is put through rigorous testing to ensure high precision, accuracy, sensitivity, and specificity. You can be confident that you will obtain reproducible results, day-after day and lot-after-lot.

*Strict QC Guidelines  
Ensure Consistent Results,  
Lot-After-Lot*

**Reliable Data Lot-After-Lot**




Graph demonstrates the robust and reproducible nature of the competitive cAMP ELISA kit (ADI-900-067A) showing standard curves from 10 lots manufactured over 10 years.

# THE BASICS OF IMMUNOASSAYS


Understanding the basic principles of immunoassays is easy. The essential components of antibody-based immunoassay systems are threefold: 1) an antigen to detect and perhaps quantitate, 2) a specific antibody to this antigen, and 3) a system to measure the amount of antigen in a given sample. Although it appears to be a very simple system, in many cases a number of other assay materials are necessary to allow for quick and convenient measurement.

## Immunometric/Sandwich ELISA




Immunometric assays, also known as sandwich ELISAs (enzyme-linked immunosorbent assays), use two antibodies specific to the antigen to capture or “sandwich” antigen in the well for detection. Immunometric assays exhibit a direct correlation between antigen concentration and substrate response. Immunometric assays typically employ a “capture” antibody coated on the plate to bind the antigen of interest. During a second incubation, the antigen is bound by a second “detection” antibody that is also specific to the antigen. The detection antibody can either be bound by a secondary antibody-enzyme conjugate, or the detection antibody itself is enzyme-conjugated. When chromogenic substrate is added to the assay to develop color, samples with high antigen concentration generate more signal than those with low antigen concentration, producing a signal directly proportional to the amount of antigen in the sample. This correlation can then be used to extrapolate the concentration of antigen in an unknown sample from a standard curve.

**Typical Standard Curve of a Sandwich ELISA**




Standard curve experiment using the HSP90α (human), ELISA Kit (ADI-EKS-895).

## Competitive ELISA



In competitive enzyme immunoassays, the antigen in a sample competes for limited antibody binding sites with antigen conjugated to a reporter enzyme. This produces an inverse relationship between antigen concentration and substrate turnover. Competitive ELISAs typically use a single antibody to a small molecular weight antigen, generally less than 10,000 Daltons. During incubation, samples with high antigen content result in unlabeled antigen being bound in greater amounts than conjugated antigen. When chromogenic substrate is added to the assay to develop color, samples with high antigen concentration generate a lower signal than those containing low antigen concentration, yielding the inverse correlation between antigen concentration in the sample and color development in the assay. This relationship can then be used to extrapolate antigen concentration in an unknown sample from a standard curve. This type of reaction is one of the few methods possible for small molecular weight antigens, such as steroids, drugs, lipids and peptides.

**Typical Standard Curve of a Competitive ELISA**



Standard curve experiment using Histamine ELISA Kit (ENZ-KIT140).

# ELISA KITS BY TARGET & RESEARCH AREA

| Target  | Alternative Names  | PAGE # | Bioprocess | Cancer | Cardiovascular | Cell Death | Cell Signaling | Cyclic Nucleotides | Cytokines | Eicosanoids | Endocrinology/ Hormones | Epigenetics | Immunology/ Inflammation | Metabolism | Nephrology | Neuroscience | Oxidative Stress | Proteostasis (HSP) |
|---|--|--------|------------|--------|----------------|------------|----------------|--------------------|-----------|-------------|-------------------------|-------------|--------------------------|------------|------------|--------------|------------------|--------------------|
| 11-dehydro-TXB <sub>2</sub>                   | 11-dehydro-Thromboxane B2  | 32     |            |        | •              |            |                |                    |           | •           | •                       |             | •                        |            |            |              |                  |                    |
| 15-deoxy-Δ <sup>12,14</sup> -PGJ <sub>2</sub> | 15-deoxy-delta-12,14-Prostaglandin J2  | 32     |            |        |                |            |                |                    |           | •           | •                       |             | •                        |            |            |              |                  |                    |
| 17β-Estradiol                                 | Estrogen, Oestradiol   | 25     |            |        |                |            |                |                    |           |             | •                       |             |                          |            |            |              | •                |                    |
| 24(S)-Hydroxycholesterol                      | 24(S)-OHC  | 38     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            | •            |                  |                    |
| 25(OH) Vitamin D                              | 25-Hydroxyvitamin D, Vitamin D   | 34     |            | •      | •              |            |                |                    |           |             | •                       |             |                          | •          |            |              |                  |                    |
| 5-Hydroxymethylcytosine                       | 5-hmC  | 26     |            | •      |                |            |                |                    |           |             |                         | •           |                          |            |            | •            |                  |                    |
| 5-Methylcytosine                              | 5-mC   | 26     |            | •      |                |            |                |                    |           |             |                         | •           |                          |            |            |              |                  |                    |
| 6-keto-PGF <sub>1α</sub>                      | Prostaglandin F <sub>1α</sub>  | 32     |            |        |                |            |                |                    |           | •           | •                       |             | •                        |            |            |              |                  |                    |
| ACTH  | Adrenocorticotrophic hormone, Corticotropin  | 25     |            |        |                |            |                |                    |           |             | •                       |             |                          |            |            | •            |                  |                    |
| Adiponectin                                   | ACRP30, ADIPOQ, GBP28, Adipocyte C1q and collagen domain-containing protein, Gelatin-binding protein, ACDC, APM1 | 34     |            |        |                |            |                |                    | •         |             | •                       |             |                          | •          |            |              |                  |                    |
| ADMA  |  | 40     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Akt   | PKB, Protein kinase B  | 20     |            |        |                | •          |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Aldosterone                                   |  | 25     |            |        |                |            | •              |                    |           |             | •                       |             |                          |            |            |              |                  |                    |
| Angiotensin                                   |  | 19     |            |        | •              |            | •              |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Annexin V                                     |  | 19     |            |        | •              | •          |                |                    |           |             |                         |             |                          | •          |            |              | •                |                    |
| Antileukoprotease                             | SLPI   | 17     | •          |        |                |            |                |                    |           |             |                         |             |                          |            |            | •            |                  |                    |
| ApoE  | Apolipoprotein E   | 34     |            |        |                |            |                |                    |           |             |                         |             |                          | •          |            |              |                  |                    |
| APP   | Amyloid precursor protein, APP neo, APP ΔC31   | 37     |            |        |                | •          |                |                    |           |             |                         |             |                          |            |            |              | •                |                    |
| Arg <sup>8</sup> -Vasopressin                 | AVP  | 19     |            |        | •              |            | •              |                    |           |             | •                       |             |                          |            |            |              |                  |                    |
| BAFF  | BLyS, TALL-1, THANK, zTNF4, TNFSF 13B/20, CD257  | 28     |            |        |                |            | •              |                    |           |             |                         |             | •                        |            |            |              |                  |                    |
| Bax   | BCL2L4, BCL2-associated X protein  | 20     |            |        |                | •          |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Bcl-2   |  | 20     |            |        |                | •          |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Big Endothelin-1                              |  | 19     |            |        | •              |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| BNP Fragment                                  |  | 19     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Bradykinin                                    |  | 19     |            |        | •              |            | •              |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Cadherin                                      |  | 21     |            |        |                |            | •              |                    |           |             |                         |             |                          |            |            | •            |                  |                    |
| CD14  |  | 28     |            |        |                |            | •              |                    |           |             |                         |             | •                        |            |            |              |                  |                    |
| CD40  | TNFSF5   | 28     |            |        |                |            | •              |                    |           |             |                         |             | •                        |            |            |              |                  |                    |
| CD40L   | CD154, TNFSF5  | 28     |            |        |                |            | •              |                    |           |             |                         |             | •                        |            |            |              |                  |                    |
| CD44  |  | 28     |            |        |                |            | •              |                    |           |             |                         |             | •                        |            |            |              |                  |                    |
| CHO HCP                                       | CHO Host Cell Protein  | 15     | •          |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| CINC-1/GRO                                    | Groα, CXCL1, rat IL-8  | 28     |            |        |                |            | •              |                    |           |             |                         |             | •                        |            |            |              |                  |                    |
| Complement C3a des Arg                        |  | 28     |            |        |                |            | •              |                    |           |             |                         |             | •                        |            |            | •            | •                | •                  |
| Complement C4a des Arg                        |  | 28     |            |        |                |            | •              |                    |           |             |                         |             | •                        |            |            | •            |                  |                    |
| Corticosterone                                |  | 25     |            |        |                |            |                |                    |           |             | •                       |             |                          |            |            | •            |                  |                    |
| Cortisol                                      |  | 25     |            |        |                |            |                |                    |           |             | •                       |             |                          |            |            |              | •                |                    |
| COX-2   | Cyclooxygenase-2   | 28     |            |        |                |            | •              |                    |           |             |                         |             | •                        |            |            |              |                  |                    |
| CRP   | C-reactive protein   | 19/28  |            |        | •              |            |                |                    |           |             |                         |             | •                        |            |            |              |                  |                    |

# ELISA KITS BY TARGET & RESEARCH AREA

| Target                    | Alternative Names  | PAGE # | Bioprocess | Cancer | Cardiovascular | Cell Death | Cell Signaling | Cyclic Nucleotides | Cytokines | Eicosanoids | Endocrinology/ Hormones | Epigenetics | Immunology/ Inflammation | Metabolism | Nephrology | Neuroscience | Oxidative Stress | Proteostasis (HSP) |
|---------------------------|--|--------|------------|--------|----------------|------------|----------------|--------------------|-----------|-------------|-------------------------|-------------|--------------------------|------------|------------|--------------|------------------|--------------------|
| Crystallin, $\alpha\beta$ | HSP25  | 42     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| cyclic AMP                | cAMP   | 23     |            | •      | •              | •          | •              | •                  |           |             |                         |             | •                        | •          |            |              |                  |                    |
| cyclic GMP                | cGMP   | 23     |            | •      | •              |            | •              | •                  |           |             |                         |             |                          |            |            |              |                  |                    |
| Cystatin                  |  | 35     |            |        | •              |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Cysteinyl leukotriene     | LTC <sub>4</sub> , LTD <sub>4</sub> , LTE <sub>4</sub>         | 32     |            |        |                |            |                |                    |           | •           | •                       |             | •                        |            |            |              |                  |                    |
| DHEA                      | Dehydroepiandrosterone   | 25     |            |        |                |            |                |                    |           |             | •                       |             |                          |            |            |              |                  |                    |
| DKK-1                     | Dickkopf-1   | 21     |            | •      |                |            | •              |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Dopamine                  | 3,4-dihydroxyphenethylamine                                    | 39     |            |        |                |            | •              |                    |           |             | •                       |             |                          |            |            | •            |                  |                    |
| <i>E. coli</i> HCP        | <i>E. coli</i> Host Cell Protein                               | 15     | •          |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Erk                       | Erk1, Erk2, MAPK   | 21     |            |        |                |            | •              | •                  |           |             |                         |             |                          |            |            |              |                  |                    |
| Estriol                   | Oestriol   | 25     |            |        |                |            |                |                    |           |             | •                       |             |                          |            |            |              |                  |                    |
| ET-1                      | Endothelin-1, Big ET-1   | 19     |            |        | •              |            | •              |                    |           |             | •                       |             |                          |            |            |              |                  |                    |
| FasL                      | AP0-1L, CD95L, CD178, TNFSF 6                                  | 28     |            |        |                |            | •              |                    |           |             |                         |             | •                        |            |            |              |                  |                    |
| Fibrinogen                |  | 19     |            |        | •              |            | •              |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Fibronectin               |  | 17     |            | •      |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| FSH                       | Follicle stimulating hormone                                   | 25     |            |        |                |            |                |                    |           |             | •                       |             |                          |            |            |              |                  | •                  |
| Gastrin                   |  | 25     |            |        |                |            |                |                    |           |             | •                       |             |                          | •          |            |              | •                | •                  |
| GLP-1                     | Glucagon-like peptide  | 13     |            |        |                |            |                |                    |           |             |                         |             |                          | •          |            |              |                  |                    |
| Grp75                     | Glucose-regulated protein                                      | 42     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            |              | •                | •                  |
| Grp78                     | BiP, Glucose-regulated protein, binding immunoglobulin protein | 42     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Grp94                     | Glucose-regulated protein                                      | 42     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            |              | •                | •                  |
| Haptoglobin               | Hp   | 19     |            |        | •              |            |                |                    |           |             |                         |             | •                        | •          |            |              |                  |                    |
| HEK239T HCP               | HEK239T Host Cell Protein                                      | 15     | •          |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Hepsin                    | TMPRSS1  | 16/17  |            | •      |                |            |                |                    |           |             |                         |             |                          |            |            |              | •                |                    |
| HETE                      | 12(S)-hydroperoxy tetraenoic eicosatetraenoic acid             | 32     |            |        |                |            | •              |                    |           | •           | •                       |             | •                        | •          |            |              |                  |                    |
| Histamine                 |  | 39     |            |        |                |            |                |                    |           |             |                         |             | •                        |            |            |              | •                |                    |
| H0-1                      | Heme Oxygenase-1   | 40     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  | •                  |
| HODE                      | 13-hydroxyoctadecadienoic acid                                 | 32     |            |        |                |            | •              |                    |           | •           | •                       |             | •                        | •          |            |              |                  | •                  |
| HSF1                      | Heat Shock Transcription Factor 1                              | 42     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  | •                  |
| HSP25                     | Heat Shock Protein 25  | 42     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  | •                  |
| HSP27                     | Heat Shock Protein 27  | 42     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  | •                  |
| HSP60                     | Heat Shock Protein 60  | 42     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  | •                  |
| HSP70                     | Heat Shock Protein 70, HSP72                                   | 42     |            | •      |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  | •                  |
| HSP70B'                   | Heat Shock Protein 70B'  | 42     |            | •      |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| HSP90 $\alpha$            | Heat Shock Protein 90  | 42     |            | •      |                |            |                |                    |           |             |                         |             |                          |            |            | •            |                  |                    |
| IFN- $\gamma$             | Interferon- $\gamma$   | 30     |            |        | •              |            |                |                    | •         |             | •                       |             | •                        |            |            |              |                  |                    |
| IGF-1                     | Insulin-like growth factor 1                                   | 28     |            |        | •              |            | •              |                    |           |             |                         |             | •                        | •          |            |              |                  |                    |
| IgG                       | Immunoglobulin G   | 28     |            |        |                |            |                |                    |           |             |                         |             | •                        |            |            |              |                  |                    |



# ELISA KITS BY TARGET & RESEARCH AREA

| Target           | Alternative Names   | PAGE # | Bioprocess | Cancer | Cardiovascular | Cell Death | Cell Signaling | Cyclic Nucleotides | Cytokines | Eicosanoids | Endocrinology/ Hormones | Epigenetics | Immunology/ Inflammation | Metabolism | Nephrology | Neuroscience | Oxidative Stress | Proteostasis (HSP) |
|------------------|---|--------|------------|--------|----------------|------------|----------------|--------------------|-----------|-------------|-------------------------|-------------|--------------------------|------------|------------|--------------|------------------|--------------------|
| IgM              | Immunoglobulin M  | 28     |            |        |                |            |                |                    |           |             |                         |             | •                        |            |            |              |                  |                    |
| IL-10            | Interleukin-10  | 30     |            |        | •              |            |                |                    | •         |             | •                       |             | •                        |            |            |              |                  |                    |
| IL-12p70         | Interleukin-12p70   | 30     |            |        | •              |            |                |                    | •         |             | •                       |             | •                        |            |            |              |                  |                    |
| IL-13            | Interleukin-13  | 30     |            |        | •              |            |                |                    | •         |             | •                       |             | •                        |            |            |              |                  |                    |
| IL-17A           | Interleukin-17A, CTLA-8, Cytotoxic T-lymphocyte-associated antigen 8  | 30     |            |        | •              |            |                |                    | •         |             | •                       |             | •                        |            |            |              |                  |                    |
| IL-1β            | Interleukin-1β  | 30     |            |        | •              |            |                |                    | •         |             | •                       |             | •                        |            |            |              |                  |                    |
| IL-2             | Interleukin-2   | 30     |            |        | •              |            |                |                    | •         |             | •                       |             | •                        |            |            |              |                  |                    |
| IL-33            | Interleukin-33, NF-HEV  | 30     |            |        | •              |            |                |                    | •         |             | •                       |             | •                        |            |            |              |                  |                    |
| IL-4             | Interleukin-4   | 30     |            |        | •              |            |                |                    | •         |             | •                       |             | •                        |            |            |              |                  |                    |
| IL-5             | Interleukin-5   | 30     |            |        | •              |            |                |                    | •         |             |                         |             |                          |            |            |              |                  |                    |
| IL-6             | Interleukin-6   | 30     |            |        | •              |            |                |                    | •         |             | •                       |             | •                        |            |            |              |                  |                    |
| IL-8             | Interleukin-8   | 30     |            |        | •              |            |                |                    | •         |             | •                       |             | •                        |            |            |              |                  |                    |
| Insulin          |   | 34     |            |        |                |            |                |                    |           |             |                         |             |                          | •          |            |              |                  |                    |
| Jnk              | SAPK, JNK 1, JNK 2  | 21     |            | •      |                | •          | •              |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Kallikrein       | KLK6, Neurosin, Protease M  | 17     |            | •      |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| KIM-1            | Kidney injury molecule-1, T-cell immunoglobulin and mucin domain-containing protein 1, TIMD-1, TIM-1, HAVCR1, Hepatitis A virus cellular receptor 1 | 34     |            |        |                |            | •              |                    |           |             |                         |             | •                        | •          | •          |              |                  |                    |
| LBP              | Lipopolysaccharide binding protein  | 28     |            |        |                |            |                |                    |           |             |                         |             | •                        |            |            |              |                  |                    |
| Leptin           | OB gene product   | 34     |            |        |                |            |                |                    |           |             | •                       |             |                          | •          |            |              |                  |                    |
| LH               | Lutenizing hormone  | 25     |            |        |                |            |                |                    |           |             | •                       |             |                          |            |            |              |                  |                    |
| LTB <sub>4</sub> | Leukotriene B <sub>4</sub>  | 32     |            |        |                |            |                |                    | •         |             | •                       |             | •                        |            |            |              |                  |                    |
| LVV Hemorphin 7  | LVVPWTQRF peptide   | 39     |            |        |                |            | •              |                    |           |             |                         |             |                          | •          |            |              |                  |                    |
| Matrptase        | TADG-15, ST14   | 17     |            | •      |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| MBL              | Mannan-binding lectin   | 28     |            | •      | •              |            |                |                    |           |             |                         |             | •                        | •          | •          |              |                  |                    |
| MCP-1            | Monocyte chemoattractant protein 1  | 30     |            |        |                |            | •              |                    |           |             |                         |             | •                        |            |            |              |                  |                    |
| MEK1             | MAP2K1, MKK1, MAPKK1, PRKMK   | 21     |            |        |                |            | •              | •                  |           |             |                         |             |                          |            |            |              |                  |                    |
| Melatonin        |   | 25     |            |        |                |            |                |                    |           |             | •                       |             |                          |            |            |              |                  |                    |
| Methotrexate     |   | 15/17  | •          | •      |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Microcystin      |   | 17     |            | •      |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Myeloperoxidase  | MPO   | 40     |            |        | •              |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Nampt            | Visfatin, Nicotinamide Phosphoribosyltransferase, PBEF1   | 34     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| NF-κB            | NF-κB p65 ELISA Kit   | 27     |            |        |                |            | •              |                    |           |             |                         |             | •                        |            |            |              |                  |                    |
| NGAL             | Lipocalin 2, Neutrophil gelatinase-associated lipocalin   | 34     |            | •      | •              |            | •              |                    |           |             |                         |             |                          | •          | •          |              |                  |                    |
| NT-proCNP        | C-type natriuretic peptide  | 19     |            |        | •              |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Osteopontin      | OPN, SPP1, Secreted phosphoprotein 1, Bone sialoprotein 1   | 21     |            | •      |                |            | •              |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Osteoprotegerin  | OPG, TNFRSF 11B, Osteoclastogenesis Inhibitory Factor, OCIF   | 30     |            |        |                |            | •              |                    | •         |             | •                       |             | •                        |            |            |              |                  |                    |
| Oxytocin         | OT  | 24/39  |            |        |                |            |                |                    |           |             | •                       |             |                          |            |            | •            |                  |                    |
| p27-Kip1         |   | 21     |            |        |                |            | •              |                    |           |             |                         |             |                          |            |            |              |                  |                    |

# ELISA KITS BY TARGET & RESEARCH AREA

| Target            | Alternative Names   | PAGE # | Bioprocess | Cancer | Cardiovascular | Cell Death | Cell Signaling | Cyclic Nucleotides | Cytokines | Eicosanoids | Endocrinology/ Hormones | Epigenetics | Immunology/ Inflammation | Metabolism | Nephrology | Neuroscience | Oxidative Stress | Proteostasis (HSP) |
|-------------------|---|--------|------------|--------|----------------|------------|----------------|--------------------|-----------|-------------|-------------------------|-------------|--------------------------|------------|------------|--------------|------------------|--------------------|
| p38               | SAPK2, SAP kinase 2   | 21     |            | •      |                |            | •              |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| p53               |   | 17     |            | •      |                | •          | •              |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| p53/MDM2          | Mouse double minute 2 homolog, E3 ubiquitin-protein ligase Mdm2 | 17     |            | •      |                | •          | •              |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| p62               | Sequestosome 1  | 20     |            | •      |                | •          |                |                    |           |             |                         |             |                          |            |            |              | •                | •                  |
| PDI               | Protein disulfide-isomerase                                     | 40     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| PEGylated protein | PEG, Polyethylene Glycol  | 15     | •          |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Peptide YY        |   | 34     |            |        |                |            |                |                    |           |             |                         |             |                          | •          |            | •            |                  |                    |
| PGE <sub>1</sub>  | Prostaglandin E <sub>1</sub>                                    | 32     |            |        |                |            |                |                    |           | •           | •                       |             | •                        |            |            |              |                  |                    |
| PGE <sub>2</sub>  | Prostaglandin E <sub>2</sub>                                    | 32     |            | •      | •              |            |                |                    |           | •           | •                       |             | •                        |            |            |              |                  |                    |
| PGF <sub>2α</sub> | Prostaglandin F <sub>2α</sub>                                   | 32     |            |        |                |            |                |                    |           | •           | •                       |             | •                        |            |            |              |                  |                    |
| Pin1              | Peptidyl-prolyl cis-trans isomerase, PPlase, Rotamase           | 17     |            | •      |                |            | •              |                    |           |             |                         |             |                          |            |            | •            |                  |                    |
| Plasminogen       |   | 19     |            |        | •              |            | •              |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| PMN-Elastase      | Polymorphonuclear elastase                                      | 28     |            |        |                |            |                |                    |           |             |                         |             | •                        |            |            |              |                  |                    |
| proANP            |   | 19     |            |        | •              |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Prolactin         | Luteotropic hormone, Luteotropin                                | 25     |            |        |                |            |                |                    |           |             | •                       |             |                          |            |            |              |                  |                    |
| Progesterone      |   | 25     |            |        |                |            |                |                    |           |             | •                       |             |                          |            |            |              |                  |                    |
| Progranulin       | Proepithelin, PEPI, PC Cell-derived Growth Factor               | 28     |            |        |                |            | •              |                    |           |             |                         |             | •                        |            |            |              |                  |                    |
| Proinsulin        |   | 34     |            |        |                |            |                |                    |           |             |                         |             |                          | •          |            |              |                  |                    |
| Prostacyclin      | PGI <sub>2</sub> , Prostaglandin I <sub>2</sub>                 | 32     |            |        |                |            |                |                    | •         |             | •                       |             | •                        |            |            |              |                  | •                  |
| Proteasome        |   | 42     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Protein A         |   | 14     | •          |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Protein carbonyl  |   | 40     | •          |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| PTX3              | Pentraxin 3, TSG14, TNF stimulated Gene-14                      | 28     |            |        |                |            | •              |                    |           | •           | •                       |             | •                        |            |            |              |                  |                    |
| RANKL             | ODF, OPL, TRANCE, TNFSF 11, CD254                               | 28     |            |        |                |            | •              |                    |           |             |                         |             | •                        |            |            |              |                  |                    |
| Sclerostin        |   | 21     |            |        |                |            |                |                    |           |             |                         |             |                          | •          |            |              |                  |                    |
| SDMA              | NG, N'-G-Dimethyl-L-arginine                                    | 40     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            |              | •                |                    |
| Serotonin         | 5HT, 5-hydroxytryptamine  | 25     |            |        |                |            | •              |                    |           |             | •                       |             |                          |            |            | •            |                  |                    |
| sHLA-G            | Soluble human leukocyte antigen-G                               | 28     |            |        |                |            |                |                    |           |             |                         |             | •                        |            |            |              | •                |                    |
| SLPI              | Secretory leukocyte proteinase inhibitor, ALP                   | 17     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            | •            |                  |                    |
| SMN               | Survival Motor Neuron   | 39     |            |        |                |            | •              |                    |           |             |                         |             |                          |            |            | •            |                  |                    |
| SOD               | Superoxide dismutase  | 40     |            |        |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Substance P       |   | 39     |            |        |                |            | •              |                    |           |             | •                       |             |                          |            |            |              |                  |                    |
| Survivin          | TIAP  | 17     |            | •      |                |            | •              |                    |           |             |                         |             |                          | •          |            |              |                  |                    |
| sVEGFR            | soluble Vascular Endothelial Growth Factor Receptor 1           | 21     |            | •      | •              |            | •              |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Testosterone      |   | 25     |            |        |                |            | •              |                    |           |             | •                       |             |                          |            |            |              |                  |                    |
| TGF-α             |   | 17     |            | •      |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| TGF-β1            | Transforming growth factor β1                                   | 30     |            |        |                |            | •              |                    | •         |             | •                       |             | •                        |            |            |              |                  |                    |
| TIMP-1            |   | 17     |            | •      |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |

## ELISA KITS BY TARGET & RESEARCH AREA

| Target           | Alternative Names                           | PAGE # | Bioprocess | Cancer | Cardiovascular | Cell Death | Cell Signaling | Cyclic Nucleotides | Cytokines | Eicosanoids | Endocrinology/ Hormones | Epigenetics | Immunology/ Inflammation | Metabolism | Nephrology | Neuroscience | Oxidative Stress | Proteostasis (HSP) |
|------------------|---|--------|------------|--------|----------------|------------|----------------|--------------------|-----------|-------------|-------------------------|-------------|--------------------------|------------|------------|--------------|------------------|--------------------|
| TL1A             | TNFSF 15, VEGI                              | 30     |            |        |                |            | •              |                    | •         |             | •                       |             | •                        |            |            |              |                  |                    |
| TNF-R1           | Tumor necrosis factor Receptor 1, TNFRSF 1A | 30     |            |        |                |            | •              |                    |           |             |                         |             | •                        |            |            |              |                  |                    |
| TNF- $\alpha$    | Tumor necrosis factor- $\alpha$ , TNFSF 2   | 30     |            |        | •              |            | •              |                    | •         |             | •                       |             | •                        |            |            |              |                  |                    |
| Total PSA        |   | 17     |            | •      |                |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Transferrin      |   | 15     | •          |        | •              |            |                |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| Troponin I       |   | 19     |            |        | •              |            | •              |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| TSH              | Thyroid stimulating hormone                 | 25     |            |        |                |            |                |                    |           |             | •                       |             |                          |            |            |              |                  |                    |
| VEGF             | Vascular endothelial growth factor          | 19     |            | •      | •              |            | •              |                    |           |             |                         |             |                          |            |            |              |                  |                    |
| XIAP             | X-linked inhibitor of apoptosis             | 20     |            |        |                | •          |                |                    |           |             |                         |             |                          | •          |            |              |                  |                    |
| $\beta$ -Catenin | CTNNB, Cadherin-associated protein          | 21     |            |        |                |            | •              |                    |           |             |                         |             |                          |            |            |              |                  |                    |

# ELISA AMPLIFICATION TECHNOLOGY


## AMP'D® SIGNAL AMPLIFICATION TECHNOLOGY

### Increase ELISA Sensitivity with the AMP'D® Signal Amplification Technology

The AMP'D® ELISA Signal Amplification system is designed to replace traditional alkaline phosphatase (AP) substrates, such as pNpp (p-Nitrophenyl phosphate), with a combination substrate and amplifier system that results in greater sensitivity when compared to a traditional substrate ELISA.

In a conventional detection system, enzyme bound to the microtiter plate interacts directly with the substrate producing a color change where the resulting absorbance is directly proportional to the amount of captured analyte. In the AMP'D® ELISA system, bound AP converts a substrate that is utilized in a second enzyme reaction system which is initiated by addition of the amplifier reagent. It is this amplification step that allows for greater (amplified) color production at lower analyte concentrations resulting in an increase in assay sensitivity.

In the AMP'D® ELISA system, the critical substrate component, NADPH, is added to wells containing AP and the AP reduces to NADH via release of a phosphate group. This reaction is allowed to proceed for an amount of time with the accumulated NADH being proportional to the amount of analyte/bound AP-conjugate. Upon the addition of the reconstituted amplifier reagent, this first reaction is quenched and the NADH feeds a second redox enzyme system. Here diaphorase utilizes NADH to reduce the idonitrotetrazolium salt into formazan (purple color) producing NAD+. A counter enzymatic reaction then occurs where the NAD+ is reduced to NADH while ethanol is oxidized to acetaldehyde via alcohol dehydrogenase. This set of enzymatic reactions is also allowed to proceed for a period of time, recycling the NADH thus amplifying the original AP/substrate reaction. The resulting color intensity is ultimately proportional to the amount of bound analyte.




## AMP'D® ELISA SIGNAL AMPLIFICATION KIT (ENZ-KIT-100)

### Increase ELISA Sensitivity with the Signal Amplification AMP'D® Technology

The AMP'D® ELISA Signal Amplification kit provides up to 50-fold increase in sensitivity over traditional ELISAs while detecting lower concentrations of target in samples.

- Quantify difficult-to-detect analytes
- Easy-to-use, simple procedure for sandwich format ELISAs containing HRP/TMB detection systems
- Convenient one or five 96-well plate formats for high-throughput analysis

### Detect up to 50-fold Increase in Sensitivity



A resulting purple color intensity is measured in the AMP'D® ELISA Signal Amplification kit (ENZ-KIT-100) that is proportional to the amount of bound analyte in your sample.

#### RELATED PRODUCTS

| Product Name   | Product #  | Size |
|--|------------|------|
| Goat anti-rabbit IgG, polyclonal antibody (AP conjugate) | ENZ-ABS257 | 5 mL |
| Streptavidin (AP Conjugate)                              | ENZ-PRT119 | 1 mL |

# ELISA AMPLIFICATION TECHNOLOGY

## AMP'D® HSP70 HIGH SENSITIVITY ELISA KIT (ENZ-KIT-101)

HSP70 functions in folding of newly synthesized proteins, re-folding of misfolded or denatured proteins, trafficking of proteins across cellular membranes, inhibiting protein aggregation, and coordinating proteins for degradation via the proteasomal pathway. Enzo's AMP'D® HSP70 high sensitivity ELISA kit provides the ability to use less sample and detect both baseline and upregulated levels of human, mouse, and rat HSP70 (HSP72).

Exposure of cells to oxidative and environmental stresses frequently results in the breakdown or oxidation of genomic DNA. Assays to evaluate the integrity of genomic DNA, or to assess the presence of oxidized DNA are frequently used as a means of verifying the onset of apoptosis or DNA damage.

- Detect as little as 7 pg/mL of HSP70 (HSP72) in < 4.5 hours
- Quantify both baseline and upregulated levels of HSP70 in serum and plasma
- Negligible reactivity from similar HSP70 family members

### Reliable Quantification, Even for Baseline Detection

| Reliable Quantification, Even for Baseline Detection |                     |               |          |             |
|--|---------------------|---------------|----------|-------------|
| Product  | Sensitivity (pg/mL) | Range (pg/mL) | Baseline | Upregulated |
| AMP'D® HSP70 Kit                                     | 7                   | 39-500        | ✓        | ✓           |
| Competitor A Kit                                     | 200                 | 780-50,000    |          | ✓           |
| Competitor R Kit                                     | 156                 | 156-10,000    |          | ✓           |
| Competitor S Kit                                     | 20                  | 550-35,000    |          | ✓           |
| Competitor U Kit                                     | 90                  | 200-12,500    |          | ✓           |


The sensitivity and range of HSP70 (HSP72) levels of the AMP'D® HSP70 high sensitivity ELISA kit (ENZ-KIT-101) was compared to 4 other similar ELISAs. Results indicate that only the AMP'D® HSP70 ELISA kit was able to detect baseline HSP70 (HSP72) levels.

## AMP'D® GLP-1 ELISA KIT (ENZ-KIT104)

GLP-1 functions in the regulation of blood glucose through insulin, and plays a role in obesity and diabetes. It is the most potent promoter of insulin secretion and an important gastrointestinal peptide hormone, or incretin. In response to eating, increased levels of GLP-1 stimulate insulin secretion and inhibit glucagon release leading to lower blood glucose levels. These positive effects of GLP-1 on glucose regulation make it an important target for the management of type 2 diabetes and the study of obesity.

- Ultra-sensitive measurement of GLP-1, detecting as little as 5.5 pg/mL
- Negligible reactivity from similar GLP-1 forms
- Compatible with serum and plasma human samples
- High-throughput format with results in 3 hours for up to 40 samples in duplicate
- Fully quantitative results that surpass semi-quantitative Western blot analysis

### Ultra-sensitive Detection of GLP-1



Standard curve of GLP-1 ELISA Kit (ENZ-KIT104).

### PRODUCT LISTING

| Product Name                            | Product #   | Size                     | Sensitivity (Range)                        | Species             | Sample Types        | Assay Time                     |
|---|-------------|--------------------------|--|---------------------|---------------------|--------------------------------|
| AMP'D® ELISA Signal Amplification Kit   | ENZ-KIT-100 | 1x96 wells<br>5x96 wells | Improves ELISA sensitivity ~10 to 50-fold. | ELISA Kit dependent | ELISA Kit dependent | ~30 Minutes (replacement time) |
| AMP'D® HSP70 High Sensitivity ELISA Kit | ENZ-KIT-101 | 1x96 wells               | 7 pg/mL (39-500 pg/mL)                     | H, M, R             | P, S                | 4.5 Hours                      |
| AMP'D® GLP-1 ELISA Kit                  | ENZ-KIT104  | 1x96 wells               | 5.5 pg/mL (7.8-250 pg/mL)                  | H                   | P, S                | 3 Hours                        |

#### LEGEND

Species: B= bovine, C = canine (dog), H = human, M = mouse, MK = monkey, P = porcine (pig), R= rat, RB = rabbit

Sample Types: AM = amniotic fluid, CL = cell lysates, CSF = cerebral spinal fluid, CS = culture supernatants, F = feces, M = milk, NL = nasal lavage, P = plasma, S = serum, SA = saliva, SS= sputum supernatant, T = tissue, TC = tissue culture, U = urine, WhB = white blood cells. Product tables includes both certified and cited species and sample types.



# BIOPROCESS OPTIMIZATION


Development of a single drug, whether it is a new chemical entity, biotherapeutic, or genetic/cellular therapy, requires significant investment of resources. Each step of the process from early discovery through production and delivery must be fully explored, characterized, and understood to prevent issues resulting from cell stress, cell death, protein aggregation, and factors affecting reliable manufacturing of the drug. Whether you work in drug discovery, upstream, or downstream bioprocessing, Enzo offers a range of products to help you maintain cell line viability, optimize and monitor product integrity, and maximize yield. Our ELISAs for protein optimization and contamination monitoring are easy to run, quantitative assays, compatible for both manual and automated workflows.

## PROTEIN A ELISA KIT (ADI-900-057)

The sensitive PEGylated Protein ELISA kit is ideal for drug development and pharmaceutical manufacturing applications including drug formulations, pharmacokinetics analysis, drug comparison, lead candidate identification, lot-release criteria and in-process QC studies.

- Validated for a wide range of molecular weight linear and branched PEGs, both in free form and when conjugated to proteins
- Sensitive assay measures < 1 ng/mL of PEGylated molecules
- Quantifies PEGylated target molecules in complex matrices to allow monitoring of drug levels or its accumulation in tissue

### Recognize All Commonly Used Protein A Constructs




Assay recognition of different Protein A constructs using the Protein A ELISA kit (ADI-900-057). A: Natural Protein A from *S. aureus* (Millipore), B: Recombinant Protein A from *E. coli* (Repligen), C: Recombinant Cys-Protein A from *E. coli* (GE), and D: Recombinant alkaline-resistant Protein A variant from *E. coli* (MabSelet SuRe from GE).

## PEGYLATED PROTEIN ELISA KIT (ADI-900-213)

The Protein A ELISA kit is a sensitive and reproducible sandwich assay to quantify Protein A residuals in monoclonal antibody preparations. This extensively validated ELISA kit enables efficient detection of natural and recombinant Protein A constructs with up to 100% recovery.

- Detects 1 ppm of Protein A residuals in human IgG
- Useful for contamination analysis and measurement of Protein A variants in monoclonal antibody preparations
- Produces results in < 3 hours with low cost per test

### Quantify a Variety of PEGylated Molecule Constructs



Standard curves using the PEGylated Protein ELISA kit (ADI-900-213) were generated to detect a variety of linear and branched PEGs.


## HOST CELL PROTEIN ELISA KITS

### Quantify HCP contamination with sensitive ELISAs for CHO- and *E. coli*-derived biologics

The production and purification of biopharmaceuticals by cell culture expression systems can produce contaminants in the form of endogenous cell line proteins, typically referred to as host cell proteins (HCPs). HCP impurities can elicit immune responses, therefore, contamination must be minimized and monitored in biologics. Enzo provides sensitive ELISAs for the detection of HCPs in bulk products expressed in CHO and *E. coli* expression systems.

- Quantitatively measure host cell protein contamination in CHO- and *E. coli*-derived biologics
- Sensitive measurement of HCPs, detecting as little as 10 ng/mL (CHO) and 30 ng/mL (*E. coli*)
- High-throughput format with results in 2.5 hours

#### Sensitive Measurement of HCPs



Typical standard curve for CHO HCP ELISA Kit (ENZ-KIT128).

## BIOPROCESS ELISA KITS

| Product Name                               | Product #   | Size       | Sensitivity (Range)            | Species                          | Sample Types                        | Assay Time |
|--|-------------|------------|--------------------------------|----------------------------------|-------------------------------------|------------|
| CHO Host Cell Protein ELISA Kit            | ENZ-KIT128  | 1x96 wells | 30 ng/mL (3-810 ng/mL)         | CHO Host Cell Protein            | CS                                  | 3 Hours    |
| <i>E. coli</i> Host Cell Protein ELISA Kit | ENZ-KIT127  | 1x96 wells | 10 ng/mL (3-810 ng/mL)         | <i>E. coli</i> Host Cell Protein | CS                                  | 3 Hours    |
| HEK293T host cell protein ELISA Kit        | ENZ-KIT162  | 1x96 wells | 37 ng/mL (37-27,000 ng/mL)     | HEK293T Host Cell Protein        | CS                                  | 6 Hours    |
| Methotrexate ELISA Kit                     | ENZ-KIT142  | 1x96 wells | 0.087 ng/mL (0.13-1,000 ng/mL) | H, M, R                          | P, S, U                             | 1.5 Hours  |
| PEGylated Protein ELISA Kit                | ADI-900-213 | 1x96 wells | < 1 ng/mL (1.75-225 ng/mL)     | Species Independent              | P, S, T, other biological samples   | 2 Hours    |
| Protein A ELISA Kit                        | ADI-900-057 | 1x96 wells | 9.01 pg/mL (15.63-1,000 pg/mL) | Species Independent              | Protein A purified IgG preparations | < 3 Hours  |
| Transferrin ELISA Kit                      | ENZ-KIT143  | 1x96 wells | 4.6 ng/mL (4.9-5,000 ng/mL)    | H, C                             | P, S                                | 2 Hours    |

#### LEGEND

Species: B= bovine, C = canine (dog), H = human, M = mouse, MK = monkey, P = porcine (pig), R= rat, RB = rabbit

Sample Types: AM = amniotic fluid, CL = cell lysates, CSF = cerebral spinal fluid, CS = culture supernatants, F = feces, M = milk, NL = nasal lavage, P = plasma, S = serum, SA = saliva, SS= sputum supernatant, T = tissue, TC = tissue culture, U = urine, WhB = white blood cells. Product tables includes both certified and cited species and sample types.

# CANCER RESEARCH

Your aim may be to investigate cancer's genetic origins, establish the mechanism of action or clinical implications of tumor biomarkers, or develop a targeted pharmaceutical or cellular therapy. At Enzo, our aim is to support your research. As cancer hallmarks have evolved over the last decade, so has our portfolio of tools to support cancer research. Over this time, we have called upon our diverse scientific and technical expertise to create innovative assays and reagents in the fields of epigenetics, autophagy, and proteostasis research. These novel tools, and a broad portfolio of products established for decades in peer-reviewed literature, will help support cancer discovery into the next decade, and beyond. We are continually focused on enabling a future with more hope, more collaboration, more discovery, and less disease. Choose from our sensitive ELISA kits to detect a variety of relevant biomarkers.

## Analyze Novel Serine Protease Cancer Biomarkers

|                 | Serine Protease Expression in Cancer |          |       |            |             |         |                  |          |      |       |      |         |                |            |          |                |      |
|-----------------|--------------------------------------|----------|-------|------------|-------------|---------|------------------|----------|------|-------|------|---------|----------------|------------|----------|----------------|------|
| Serine Protease | Breast                               | Cervical | Colon | Endothelia | Endometrial | Gastric | Gastrointestinal | Leukemia | Lung | Renal | Oral | Ovarian | Nasopharyngeal | Pancreatic | Prostate | Salivary Gland | Skin |
| Hepsin          | ↑                                    | ✓        |       |            | ↑           |         |                  |          | ✓    | ↑     |      | ↑       |                |            | ↑        |                |      |
| KLK6            |                                      |          |       |            |             |         |                  |          |      |       |      | ✓       |                |            | ✓        |                |      |
| KLK7            | ↑↓                                   | ↑        | ↑     |            |             |         |                  |          |      |       | ↑    | ↑       |                | ✓          | ↑        |                |      |
| KLK8            | ✓                                    | ↑        |       |            | ↑           |         |                  |          | ↑↓   |       | ↑    | ↑       |                |            |          | ↑              |      |
| Matriptase      | ↑↓                                   | ↑        | ↑↓    | ↑          | ↑           | ↑       | ↓                | ↑        | ↓    | ↑     |      | ↑↓      |                | ↑          | ↑        |                | ↑    |
| SLP1            | ✓                                    | ↓        | ✓     |            |             | ↑       |                  |          | ↑    |       |      | ↑↓      | ↓              |            | ↓        |                |      |
| TMPR553         | ✓                                    |          |       |            |             |         |                  |          | ✓    |       |      | ↑       |                | ↑          | ✓        |                |      |


Serine Protease ELISAs for early stage cancer detection. Checks (✓) indicate expression, while arrows (↑↓) indicate increase or decrease in the level of expression in the different cancer types.

## SERINE PROTEASE CANCER BIOMARKERS

Serine proteases are enzymes that mediate a variety of events relevant to fundamental processes of tumor invasion and metastasis. In ovarian cancer, they allow growth and spread of the cancer. Serine proteases are produced early in tumor formation and can be detected early, and when inhibited, they may possibly stop ovarian tumor formation. They do this by (1) breaking down extracellular components surrounding ovaries and (2) allowing tumor cells to rapidly divide. Then tumor cells can break off and travel to different locations, resulting in the spread of cancer. Enzo Life Sciences offers a complete portfolio of over 40 unique cancer biomarker products including ELISAs, monoclonal antibodies, and recombinant proteins. The serine protein ELISAs are rapid assays enabling the measurement of pivotal biomarkers of ovarian, breast, cervical, prostate and other cancers without the need for expensive equipment or non-quantitative procedures.

- Easy-to-use, rapid ELISA kits with results in ~3 hours
- Nanogram level detection of unique serine protease cancer biomarkers
- Fully quantitative results surpass semi-quantitative Western blot analysis

### Rapid, Sensitive ELISA for Matriptase Detection




Standard curve using the Matriptase ELISA Kit (ADI-900-221).

## TOTAL PSA ELISA KIT (ENZ-KIT146)

Prostate-specific antigen (PSA), also known as gamma-seminoprotein or kallikrein-3 (KLK3), is a kallikrein-like protease produced by the epithelial cells in the prostate gland, the lining of the urethra, and the bulbourethral gland. It is a single-chain glycoprotein with a molecular weight of about 30-33 kDa in human seminal plasma. PSA exists in serum in multiple forms and most PSA is bound to serum proteins. Increases in glandular size and tissue damage caused by benign prostatic hypertrophy, prostatitis, or prostate cancer may increase circulating total PSA levels. Monitoring the total PSA level in male blood has been used as an indicator for abnormal prostate function in the past decades. The Total PSA (human) ELISA kit is highly specific, with no cross reactivity with KLK2 and ACT.

- Highly-sensitive measurement of Total PSA, detecting as little as 0.287 pM (or 0.024 ng/mL) with negligible cross reactivity with similar proteins (KLK2, ACT)
- Compatible with various sample types (plasma, serum, urine, tissue culture media)
- High-throughput format with results in 2.5 hours for up to 40 samples in duplicate
- Fully quantitative results that surpass semi-quantitative Western blot analysis

### Sensitive Detection of Total PSA



Standard curve experiment using Total PSA (human) ELISA Kit (ENZ-KIT146).

## SELECT ELISA KITS FOR CANCER RESEARCH

| Product Name                                      | Product #   | Size       | Sensitivity (Range)                                    | Species | Sample Types                      | Assay Time  |
|---|-------------|------------|--|---------|-----------------------------------|---|
| Fibronectin ELISA Kit                             | ENZ-KIT135  | 1x96 wells | < 10 pg/mL (156-10,000 pg/mL)                          | H       | CS, S, P (Heparin, EDTA, Citrate) | 3.5 Hours   |
| Hepsin (human) ELISA Kit                          | ADI-900-220 | 1x96 wells | 0.2 U/mL (16-256 U/mL)                                 | H       | P, S                              | ~3 Hours  |
| IMMUNOSET® p53/MDM2 Complex ELISA Development Set | ADI-960-070 | 5x96 wells | 0.35 ng/mL p53 (0.78-50 ng/mL p53)                     | H, M, R | CL                                | Plate coating - Overnight + 1 Hour; Assay - 3 Hours |
| Kallikrein-6 (human) ELISA Kit                    | ADI-900-217 | 1x96 wells | 1 ng/mL (5-280 ng/mL)                                  | H       | S                                 | ~3 Hours  |
| Kallikrein-7 (human) ELISA Kit                    | ADI-900-218 | 1x96 wells | 0.5 ng/mL (35-560 ng/mL)                               | H       | S                                 | ~3 Hours  |
| Kallikrein-8 (human) ELISA Kit                    | ADI-900-219 | 1x96 wells | 0.1 ng/mL (35-560 ng/mL)                               | H       | S                                 | ~3 Hours  |
| Matriptase (human) ELISA Kit                      | ADI-900-221 | 1x96 wells | 1 ng/mL (5-280 ng/mL)                                  | H       | P, S                              | ~3 Hours  |
| Methotrexate ELISA Kit                            | ENZ-KIT142  | 1x96 wells | 0.087 ng/mL (0.13-1,000 ng/mL)                         | H, M, R | P, S, U                           | 1.5 Hours   |
| Microcystins (Adda specific) ELISA Kit            | ALX-850-319 | 1x96 wells | 0.1 ng/mL (0.15-5 ng/mL)                               | NA      | Not applicable                    | 2.5 Hours   |
| p53 (human) ELISA Kit                             | ALX-850-057 | 1x96 wells | 0.5 U/mL   | H       | CS, P, S                          | 3 Hours   |
| Pin1 ELISA Kit                                    | ADI-900-146 | 1x96 wells | 15.5 pg/mL (62.5-2,000 pg/mL)                          | H, M    | CL                                | 3 Hours   |
| SLPI (human) ELISA Kit                            | ADI-900-222 | 1x96 wells | 1 ng/mL (10-160 ng/mL)                                 | H       | S                                 | 2.5 Hours   |
| Survivin (human) ELISA Kit                        | ADI-900-111 | 1x96 wells | 4 pg/mL (31.25-1,000 pg/mL)                            | H       | CL, CS, P, S, U                   | 3 Hours   |
| TGF-α ELISA Kit                                   | ENZ-KIT132  | 1x96 wells | 9.375 pg/mL  | H       | P, S, CS, T                       | 3.5 Hours   |
| TIMP-1 (human) ELISA Kit                          | ENZ-KIT147  | 1x96 wells | ≤ 30 pg/mL (0.049-12 ng/mL)                            | H       | P, S, SA, TC, U                   | 2 Hours   |
| Total PSA (human) ELISA Kit                       | ENZ-KIT146  | 1x96 wells | 0.287 pM or 0.024 ng/mL (333-1.37 pM or 28-0.12 ng/mL) | H       | P, S, TC, U                       | 2.5 Hours   |

#### LEGEND

Species: B= bovine, C = canine (dog), H = human, M = mouse, MK = monkey, P = porcine (pig), R= rat, RB = rabbit

Sample Types: AM = amniotic fluid, CL = cell lysates, CSF = cerebral spinal fluid, CS = culture supernatants, F = feces, M = milk, NL = nasal lavage, P = plasma, S = serum, SA = saliva, SS= sputum supernatant, T = tissue, TC = tissue culture, U = urine, WhB = white blood cells. Product tables includes both certified and cited species and sample types.

# CARDIOVASCULAR RESEARCH


The cardiovascular or circulatory system provides the body with essential nutrients and oxygen. Besides its nourishing function, it also plays a crucial role in inflammation and fighting infections, maintaining homeostasis, transporting hormones, and stabilizing body temperature. It also assists the body in recovering from damage and trauma. Although the system is extremely resilient, repeated or continuous stress can cause long-term damage and lead to cardiovascular diseases such as heart failure and stroke. Cardiovascular disease is consistently ranked as the single largest cause of death in the world, and it may be caused by heart damage or vascular problems. Challenges in cardiovascular clinical development have caused researchers to take a closer look at cardiac function and metabolism. Enzo Life Sciences offers a comprehensive portfolio of products to enable discovery of cardiac risk factors as well as analysis of the cellular response to novel therapeutics for cardiovascular medicine.

## TROPONIN I (HUMAN) ELISA KIT (ADI-900-228)

Troponin I is an inhibitory troponin subunit that acts in muscle relaxation and which, in its cardiac isoform, is a reliable marker of cardiac injury when identified in elevated levels in the bloodstream.

- Most sensitive ELISA, detecting as little as 0.38 ng/mL of human troponin
- Flexible sample types, including cell lysates, plasma, and serum
- Simple protocol with results in just 3.5 hours

### Most Sensitive ELISA for Detection of Troponin I




Standard curve experiment using the Troponin I ELISA Kit (ADI-900-228).

## ARG<sup>8</sup>-VASOPRESSIN ELISA KIT (ADI-900-017A)

Highly sensitive ELISA kit for the detection of Arg<sup>8</sup>-Vasopressin in serum, plasma and tissue culture media.

- High sensitivity, detecting as little as 2.84 pg/mL
- Highly specific with negligible reactivity with oxytocin
- Simple and easy-to-use liquid color-coded reagents reduces error
- Fully quantitative results that surpass semi-quantitative Western blot analysis

### Typical Standard Curves for Arg<sup>8</sup>-Vasopressin



Standard curve experiment using Arg<sup>8</sup>-Vasopressin ELISA Kit (ADI-900-017A).



## SELECT ELISA KITS FOR CARDIOVASCULAR RESEARCH

| Product Name                            | Product #    | Size       | Sensitivity (Range)             | Species              | Sample Types             | Assay Time         |
|---|--------------|------------|---------------------------------|----------------------|--------------------------|--------------------|
| Angiotensin A ELISA Kit                 | ADI-900-207  | 1x96 wells | 8.3 pg/mL (9.8-10,000 pg/mL)    | Species Independent  | P, S                     | 3 Hours            |
| Angiotensin I ELISA Kit                 | ADI-900-203  | 1x96 wells | 4.3 pg/mL (3.9-10,000 pg/mL)    | Species Independent  | P, S                     | 3.5 Hours          |
| Angiotensin II ELISA Kit                | ADI-900-204  | 1x96 wells | 4.6 pg/mL (3.9-10,000 pg/mL)    | Species Independent  | P, S                     | 3.5 Hours          |
| Annexin V (human) ELISA Kit             | ALX-850-049  | 1x96 wells | 0.33 ng/mL (0.8-50 ng/mL)       | H                    | CS, S                    | 3.5 Hours          |
| Anti-Annexin V (human) ELISA Kit        | ALX-850-040  | 1x96 wells | 1.18 ng/mL (6.25-400 ng/mL)     | H                    | CS, S                    | 4 Hours            |
| Arg <sup>8</sup> -Vasopressin ELISA Kit | ADI-900-017A | 1x96 wells | 2.84 pg/mL (4.10-1,000 pg/mL)   | H                    | P, S, TC                 | Overnight + 1 Hour |
| Big Endothelin-1 (human) ELISA Kit      | ADI-900-022  | 1x96 wells | 0.23 pg/mL (0.78-100 pg/mL)     | H, RB                | CS, Lung Lavage Fluid, S | Overnight + 1 Hour |
| Big Endothelin-1 (rat) ELISA Kit        | ADI-900-073  | 1x96 wells | 0.72 pg/mL (0.78-100 pg/mL)     | R                    | CS, P, S                 | Overnight + 1 Hour |
| BNP Fragment ELISA Kit                  | ENZ-KIT152   | 1x96 wells | 171 pmol/L (0-6400 pmol/L)      | H                    | P, S                     |                    |
| Bradykinin ELISA Kit                    | ADI-900-206  | 1x96 wells | 24.8 pg/mL (11.7-30,000 pg/mL)  | Species Independent  | P, S, U                  | 3 Hours            |
| CRP (human) ELISA Kit                   | ENZ-KIT102   | 1x96 wells | 12.685 ng/mL (10-810 ng/mL)     | H                    | P, S                     | < 2 Hours          |
| Endothelin-1 ELISA Kit                  | ADI-900-020A | 1x96 wells | 0.40 pg/mL (0.78-100 pg/mL)     | H, M, R, B, C, P, RB | CL, CS, P, S, T          | 2 Hours            |
| Fibrinogen (human) ELISA Kit            | ADI-900-230  | 1x96 wells | < 7.63 ng/mL (15.6-1,000 ng/mL) | H                    | CL, P, S                 | ~3.5 Hours         |
| Haptoglobin (human) ELISA Kit           | ADI-900-229  | 1x96 wells | < 0.78 ng/mL (0.78-50 ng/mL)    | H                    | CL, P, S                 | ~3.5 Hours         |
| NT-proCNP ELISA Kit                     | ENZ-KIT154   | 1x96 wells | 0.2 pmol/L (0-40 pmol/L)        | H                    | P, S                     | 4.5 Hours          |
| Plasminogen (human) ELISA Kit           | ADI-900-231  | 1x96 wells | < 2.01 ng/mL (1-64 ng/mL)       | H                    | CL, P, S                 | ~3.5 Hours         |
| proANP ELISA Kit                        | ENZ-KIT153   | 1x96 wells | 0.05 nmol/L (0-10 nmol/L)       | H                    | CL, P, S, U              | 3.5 Hours          |
| Transferrin ELISA kit                   | ENZ-KIT143   | 1x96 wells | 4.6 ng/mL (4.9-5,000 ng/mL)     | H, C                 | P, S                     | 2 Hours            |
| Troponin I (human) ELISA Kit            | ADI-900-228  | 1x96 wells | < 0.38 ng/mL (0.38-25 ng/mL)    | H                    | CL, P, S                 | ~3.5 Hours         |
| VEGF (human) ELISA Kit                  | ENZ-KIT156   | 1x96 wells | 4 pg/mL (8-2,000 pg/mL)         | H                    | CS, P, S                 | 2.5 Hours          |
| VEGF-C (human) ELISA Kit                | ALX-850-306  | 1x96 wells | 0.057 ng/mL (0.23-15 ng/mL)     | H                    | CS, S                    | 4.5 Hours          |

### LEGEND

Species: B= bovine, C = canine (dog), H = human, M = mouse, MK = monkey, P = porcine (pig), R= rat, RB = rabbit

Sample Types: AM = amniotic fluid, CL = cell lysates, CSF = cerebral spinal fluid, CS = culture supernatants, F = feces, M = milk, NL = nasal lavage, P = plasma, S = serum, SA = saliva, SS= sputum supernatant, T = tissue, TC = tissue culture, U = urine, WhB = white blood cells. Product tables includes both certified and cited species and sample types.

# CELL DEATH RESEARCH

Different types of cell death are often defined by morphological criteria, without a clear reference to precise biochemical mechanisms. Apoptosis (or programmed cell death) is the most well-characterized type of cell death, being recognized as a critical regulator of development and immunity, as well as organ and tissue homeostasis. Apoptotic cells die in a controlled fashion in response to a variety of extrinsic or intrinsic signals (e.g., activation of TNF receptors, DNA damage, mitochondrial pathways).

Cells can also die due to necrosis, which does not follow the apoptotic signal transduction pathway. Rather, various receptors are activated that result in the loss of cell membrane integrity and an uncontrolled release of products of cell death into the intracellular space.


Autophagy (or autophagocytosis) is another common cell death pathway which involves cell degradation of unnecessary or dysfunctional cellular components through the actions of lysosomes. Autophagic activity is critical to the maintenance of cellular homeostasis and energy balance. Although typically low under basal conditions, autophagy can be markedly upregulated by a variety of physiological stimuli such as nutrient starvation, hypoxia, endoplasmic reticulum stress, as well as immune and hormonal stimulation. With our innovative assays and reagents, Enzo enables scientists to gain a deeper understanding of the benefits, and potential consequences, of altering autophagic activity.

## p62 ELISA KITS (ADI-900-212)

Enzo's sensitive p62 ELISA kit allows for quantitative, immunometric detection of the autophagy biomarker p62 (Sequestosome 1) in human, rat and mouse cell lysates. p62 functions as a scaffold protein, aiding in autophagy protein trafficking and degradation. This ELISA kit enables quantitative measurement of autophagy without the need for expensive equipment or long procedures.

- Sensitive assay measures as little as 100 pg/mL of p62
- High-throughput format allows analysis of up to 40 samples in duplicate in < 3 hours
- Easy-to-follow protocols and liquid color-coded reagents save time and reduce errors

## p62 ELISA Measures Induction of Macroautophagy



Human breast cancer cells were treated with either autophagy-inducing drug or vehicle. Cells were harvested 6 and 12 hours post-treatment, lysed, and analyzed with the p62 ELISA Kit (ADI-900-212), and for p62 and LC3-II by Western blot. Drug treatment correlated with the induction of autophagy as indicated by the decrease in p62 levels (A and B) and by elevation of LC3-II levels (C).

## SELECT ELISA KITS FOR CELL DEATH RESEARCH

| Product Name                               | Product #   | Size       | Sensitivity (Range)             | Species | Sample Types     | Assay Time |
|--|-------------|------------|---------------------------------|---------|------------------|------------|
| [pSer <sup>473/474</sup> ]Akt1/2 ELISA Kit | ADI-900-162 | 1x96 wells | 5.5 pg/mL (17.5-560 pg/mL)      | H, M, R | CL               | 3 Hours    |
| APOSTRAND™ ELISA Apoptosis Detection Kit   | BML-AK120   | 1x96 wells | Not applicable                  | N/A     | N/A              | < 4 Hours  |
| Bax (human) ELISA Kit                      | ADI-900-138 | 1x96 wells | 25.6 pg/mL (62.5 - 2,000 pg/mL) | H       | CL               | 3 Hours    |
| Bcl-2 (human) ELISA Kit                    | ADI-900-133 | 1x96 wells | 3.8 pg/mL (18.8-1,200 pg/mL)    | H       | CL               | 3 Hours    |
| p62 ELISA Kit                              | ADI-900-212 | 1x96 wells | 100 pg/mL (625-40,000 pg/mL)    | H, M, R | CL, PBMC lysates | 3 Hours    |
| XIAP (human) ELISA Kit                     | ADI-900-124 | 1x96 wells | 90.6 pg/mL (312.5-10,000 pg/mL) | H       | CL               | 3 Hours    |

### LEGEND

Species: B= bovine, C = canine (dog), H = human, M = mouse, MK = monkey, P = porcine (pig), R= rat, RB = rabbit

Sample Types: AM = amniotic fluid, CL = cell lysates, CSF = cerebral spinal fluid, CS = culture supernatants, F = feces, M = milk, NL = nasal lavage, P = plasma, S = serum, SA = saliva, SS= sputum supernatant, T = tissue, TC = tissue culture, U = urine, WhB = white blood cells. Product tables includes both certified and cited species and sample types.

# CELL SIGNALING RESEARCH


Cell signaling encompasses a complex network of highly coordinated communication pathways from ligands down to the effector molecule. From receptor- or ion channel-mediated transmission of signals across the cell membrane, to amplification through second messenger systems, kinase or proteolytic cascades or other post-translational modifications, Enzo Life Sciences offers an unrivaled catalog of cutting-edge tools for signal transduction research. From receptor to effector, a portfolio of assays, recombinant enzymes and substrates, small molecule reagents, and antibodies is available to facilitate thorough characterization of your pathway of interest.

## SENSITIVE WNT PATHWAY ELISA KITS

Obtain fully quantitative data with immunometric detection of  $\beta$ -Catenin, total and phosphorylated GSK-3 $\beta$ , and Dkk-1 in human, mouse and rat samples.

- Sensitive assays measure picogram levels compared to microgram levels in Western blot
- Suitable for cell lysates, culture supernatants, plasma, or serum sample types
- High-throughput format allows analysis of 40 samples in duplicate in < 3 hours

### Quantify Critical Wnt Pathway Proteins



## SELECT ELISA KITS FOR CELL SIGNALING RESEARCH

| Product Name   | Product #    | Size       | Sensitivity (Range)              | Species | Sample Types   | Assay Time  |
|--|--------------|------------|----------------------------------|---------|----------------|---|
| VE-cadherin, Soluble ELISA Kit                                   | ALX-850-059  | 1x96 wells | 0.15 ng/mL (0.16-10 ng/mL)       | H       | CS, S          | 3.25 Hours  |
| $\beta$ -Catenin ELISA Kit                                       | ADI-900-135  | 1x96 wells | 26.8 pg/mL (125-8,000 pg/mL)     | H, M, R | CL             | 3 Hours   |
| Dkk-1 (human) ELISA Kit  | ADI-900-151  | 1x96 wells | 0.98 pg/mL (7.81-500 pg/mL)      | H       | CS, P, S       | 3 Hours   |
| Dkk-1 (mouse) ELISA Kit  | ADI-900-172  | 1x96 wells | 116.7 pg/mL (125-4,000 pg/mL)    | M       | CS, P, S       | 3 Hours   |
| Dkk-1 (rat), ELISA Kit   | ADI-900-171  | 1x96 wells | 26.1 pg/mL (39.1-1,250 pg/mL)    | R       | CS, P, S       | 3 Hours   |
| [pThr <sup>202</sup> /Tyr <sup>204</sup> ]Erk1/2 ELISA Kit       | ADI-900-098A | 1x96 wells | 2.67 pg/mL (62.5-2,000 pg/mL)    | H, M, R | CL             | 3 Hours   |
| [pSer] <sup>9</sup> GSK-3 $\beta$ ELISA Kit                      | ADI-900-123A | 1x96 wells | 9.0 pg/mL (62.5-2,000 pg/mL)     | H, M, R | CL             | 3 Hours   |
| IGF-1 (human) ELISA Kit  | ADI-900-150  | 1x96 wells | 48.5 pg/mL (187-6,000 pg/mL)     | H       | P, S           | 4 Hours   |
| [pThr <sup>183</sup> /Tyr <sup>185</sup> ]Jnk1/2 ELISA Kit       | ADI-900-106  | 1x96 wells | 75.8 pg/mL (125-4,000 pg/mL)     | H, M, R | CL             | 3 Hours   |
| MEK1 ELISA Kit   | ADI-900-122A | 1x96 wells | 139.0 pg/mL (312.5-10,000 pg/mL) | H, R    | CL             | 3 Hours   |
| [pSer <sup>218</sup> /Ser <sup>222</sup> ]MEK1 ELISA Kit         | ADI-900-119  | 1x96 wells | 85.2 pg/mL (187.5-6,000 pg/mL)   | H, M, R | CL             | 3 Hours   |
| Osteopontin (human) ELISA Kit                                    | ADI-900-142  | 1x96 wells | 0.110 ng/mL (2-32 ng/mL)         | H       | CS, M, P, S, U | 3 Hours   |
| Osteopontin (rodent) ELISA Kit                                   | ADI-900-090A | 1x96 wells | 3.03 ng/mL (3.13-100 ng/mL)      | M, R    | CS, P, U       | 3 Hours   |
| IMMUNOSET <sup>®</sup> Osteopontin (human) ELISA Development Set | ADI-960-142  | 1x96 wells | 0.016 ng/mL (0.03-1 ng/mL)       | H       | CS, M, P       | Plate coating - Overnight + 1 Hour; Assay - 3 Hours |
| p27-Kip1 (human) ELISA Kit                                       | ADI-900-139  | 1x96 wells | 10.1 pg/mL (25-1,600 pg/mL)      | H       | CL             | 3 Hours   |
| [pThr <sup>180</sup> /Tyr <sup>182</sup> ]p38 ELISA Kit          | ADI-900-101  | 1x96 wells | 52.1 pg/mL (156-5,000 pg/mL)     | H, M    | CL             | 3 Hours   |
| Sclerostin ELISA Kit   | ENZ-KIT155   | 1x96 wells | 3.2 pmol/L (0-240 pmol/L)        | H       | P, S, U        | Overnight + 1.5 Hours                               |
| sVEGFR-1 (human) ELISA Kit                                       | ALX-850-264  | 1x96 wells | 0.06 ng/mL (0.16-10 ng/mL)       | H       | CS, S          | 3.5 Hours   |
| VE-cadherin, Soluble ELISA Kit                                   | ALX-850-059  | 1x96 wells | 0.15 ng/mL (0.16-10 ng/mL)       | H       | CS, S          | 3.25 Hours  |

# CYCLIC NUCLEOTIDE ELISA KITS

## MOST SENSITIVE cAMP/cGMP ELISA KITS

The second messengers adenosine 3', 5'-cyclic monophosphate (cyclic AMP; cAMP) and guanosine 3', 5'-cyclic monophosphate (cyclic GMP; cGMP) are important intracellular regulators downstream from GPCRs.


**cAMP** is also involved in regulating neuronal, glandular, cardiovascular, and immune functions amongst others. A number of hormones are known to activate cAMP through the action of the enzyme adenylate cyclase which converts ATP to cAMP. The diversity of receptors known to impact cAMP-mediated signaling includes those known to affect cardiovascular and nervous systems, immune mechanisms, cell growth, differentiation, and general metabolism. Thus, there remains considerable interest in the measurement of intracellular cAMP concentrations in tissues and cell culture.

**cGMP** is typically present at levels 10- to 100-fold lower than cAMP in most tissues and is formed by the action of the enzyme guanylate cyclase on GTP. It is most well recognized as a key signaling trigger in the relaxation of smooth muscle cells (resulting in vasodilation) in response to nitric oxide, and as a sodium ion channel regulator (via its degradation by phosphodiesterases) critical to phototransduction. Vasodilators and hormones, such as acetylcholine, insulin and oxytocin, as well as serotonin and histamine, cause an increase in cGMP levels.

In addition to measuring cAMP and cGMP, detection of downstream GPCR signaling events such as activation of Erk1/2, Akt, Src and Stat3 are routinely assayed using phospho-specific antibodies to distinguish phosphorylation of signaling proteins at specific activating or inactivating serine, threonine, and tyrosine residues. We offer the most sensitive and complete colorimetric ELISA kits for quantification of intracellular or extracellular cAMP or cGMP in a variety of sample types.

- Enhance sensitivity by 10-fold with optional acetylation protocol
- Cited in peer-reviewed publications
- Simple, efficient and well-established sample handling protocols

### Reliable Results with Consistent Lot-to-Lot Performance



Lot-to-lot consistency: Graph demonstrates the robust and reproducible nature of the competitive cAMP ELISA Kit (ADI-900-067A) showing standard curves from 10 lots manufactured over 10 years.

|      |          | VALIDATED SAMPLE TYPES |                |                     |       |        |       |        |
|------|----------|------------------------|----------------|---------------------|-------|--------|-------|--------|
|      |          | Intracellular          |                | Extracellular       |       |        |       |        |
|      |          | Cell Lysates           | Tissue Extract | Culture Supernatant | Serum | Saliva | Urine | Plasma |
| cAMP | Complete | ✓                      | ✓              | ✓                   | ✓     | ✓      |       |        |
|      | Direct   | ✓                      | ✓              |                     |       |        |       |        |
|      | Standard |                        |                | ✓                   | ✓     | ✓      |       |        |
| cGMP | Complete | ✓                      | ✓              | ✓                   | ✓     | ✓      |       |        |
|      | Direct   | ✓                      | ✓              |                     |       |        |       |        |
|      | Standard |                        |                | ✓                   | ✓     | ✓      | ✓     | ✓      |

## CYCLIC NUCLEOTIDE ELISA KITS

| Product Name            | Product #                    | Size                     | Sensitivity (Range)  | Species             | Sample Types               | Assay Time |
|-------------------------|------------------------------|--------------------------|--|---------------------|----------------------------|------------|
| cAMP Complete ELISA Kit | ADI-900-163A<br>ADI-901-163A | 1x96 wells<br>5x96 wells | Non-acetylated: 0.49 pmol/ml (assay buffer) / 1.18 pmol/ml (HCl) (range (non-acetylated) 0.78 - 200 pmol/ml)<br>Acetylated: 0.027 pmol/ml (assay buffer) / 0.006 pmol/ml (HCl) (range (acetylated) 0.078 - 20 pmol/ml) | Species Independent | CL, CS, S, SA, T           | 3 Hours    |
| Direct cAMP ELISA Kit   | ADI-900-066A<br>ADI-901-066A | 1x96 wells<br>5x96 wells | Non-acetylated: 1.18 pmol/ml (range 0.78 - 200 pmol/ml)<br>Acetylated: 0.006 pmol/ml (range 0.078 - 20 pmol/ml)  | Species Independent | CL, CS, T, TC              | 3 Hours    |
| cAMP ELISA Kit          | ADI-900-067A<br>ADI-901-067A | 1x96 wells<br>5x96 wells | Non-acetylated: 0.49 pmol/ml (range 0.78 - 200 pmol/ml)<br>Acetylated: 0.027 pmol/ml (range 0.078 - 20 pmol/ml)  | Species Independent | CL, CS, S, SA, T, U        | 3 Hours    |
| cGMP Complete ELISA Kit | ADI-900-164<br>ADI-901-164   | 1x96 wells<br>5x96 wells | Assay Buffer: 0.42 pmol/mL (non-acetylated) 0.043 pmol/mL (acetylated) HCl: 0.604 pmol/mL (non-acetylated) 0.059 pmol/mL (acetylated) (0.8-500 pmol/mL (non-acetylated); 0.08-50 pmol/mL (acetylated)                  | Species Independent | CL, CS, P, S, SA, T, U     | 3 Hours    |
| Direct cGMP ELISA Kit   | ADI-900-014<br>ADI-901-014   | 1x96 wells<br>5x96 wells | Non-acetylated: 0.604 pmol/mL; (2 hour acetylated) 0.059 pmol/mL; (overnight acetylated) 0.025 pmol/mL (non-acetylated) 0.8-500 pmol/mL; (acetylated) 0.08-50 pmol/mL)   | Species Independent | CL, CSF, Microdialysate, T | 2 Hours    |
| cGMP ELISA Kit          | ADI-900-013<br>ADI-901-013   | 1x96 wells<br>5x96 wells | Non-acetylated: 0.37 pmol/mL; (acetylated) 0.088 pmol/mL (non-acetylated) 0.16-500 pmol/mL; (acetylated) 0.16-100 pmol/mL)   | Species Independent | CS, P, S, SA, U            | 3 Hours    |

### LEGEND

Species: B= bovine, C = canine (dog), H = human, M = mouse, MK = monkey, P = porcine (pig), R= rat, RB = rabbit

Sample Types: AM = amniotic fluid, CL = cell lysates, CSF = cerebral spinal fluid, CS = culture supernatants, F = feces, M = milk, NL = nasal lavage, P = plasma, S = serum, SA = saliva, SS= sputum supernatant, T = tissue, TC = tissue culture, U = urine, WhB = white blood cells. Product tables includes both certified and cited species and sample types.



# ENDOCRINOLOGY/HORMONE RESEARCH


All multicellular organisms produce hormones, which are regulatory biochemicals part of the endocrine system and serve as a major form of communication between different organs and tissues. Hormones regulate a variety of physiological and behavioral activities, including digestion, metabolism, respiration, tissue function, sensory perception, sleep, excretion, lactation, stress, growth and development, movement, reproduction, and mood. Animal hormones are classified by their chemical types (e.g. peptide- or lipid-based). The peptide hormone family includes vasopressin, insulin, luteinizing hormone, follicle-stimulating hormones, and several others, while the lipid hormones family includes phospholipid-derived hormones from the arachidonic pathway (e.g. eicosanoids) and steroid hormones like testosterone and cortisol. Enzo Life Sciences provides a variety of high sensitivity ELISA kits for researchers looking at hormone-related research.

## OXYTOCIN ELISA KIT (ADI-900-153A)

Commonly known as the “love” hormone, Oxytocin is a key neuromodulator in the brain, with defined roles in social behavior including parental nurturing, social pair-bonding, trust, and management of stress experiences. It's also a key hormone during mammalian birthing and lactation. Join the many scientists worldwide who entrust their research to Enzo's Oxytocin ELISA kit.

- Detect as low as 15 pg/mL of oxytocin in a variety of sample types
- Negligible detection of vasopressin, providing confidence in assay results
- Faster and less costly than LC/MS methods
- Widely published in peer-reviewed literature

### Most Sensitive Oxytocin ELISA Kit




Standard curve experiment using the Oxytocin ELISA Kit (ADI-900-153A).

## CORTICOSTERONE ELISA KIT (ADI-900-097)

Highly sensitive and species-independent Corticosterone ELISA kit used for animal stress research. The ELISA kit has low cross-reactivity with related steroids. The broad dynamic range makes this kit ideal for a wide variety of sample matrices from any species.

- Sensitive measurement of corticosterone, detecting as little as 27 pg/mL
- High-throughput format with results in < 3 hours for up to 39 samples in duplicate
- Broad dynamic range
- Species-independent
- Routinely used to monitor stress levels in zoo animals

### Sensitive Measurement of Corticosterone



Standard curve experiment using the Corticosterone ELISA Kit (ADI-900-097).

## SELECT ELISA KITS FOR ENDOCRINOLOGY RESEARCH

| Product Name                             | Product #                    | Size                     | Sensitivity (Range)            | Species             | Sample Types                | Assay Time         |
|--|------------------------------|--------------------------|--------------------------------|---------------------|-----------------------------|--------------------|
| 17β-Estradiol ELISA Kit                  | ADI-900-008<br>ADI-901-008   | 1x96 wells<br>5x96 wells | 28.5 pg/mL (29.3-30,000 pg/mL) | Species Independent | CS, SA                      | < 3 Hours          |
| 17β-Estradiol High Sensitivity ELISA Kit | ADI-900-174<br>ADI-901-174   | 1x96 wells<br>5x96 wells | 14.0 pg/mL (15.6-1,000 pg/mL)  | Species Independent | P, S                        | 3 Hours            |
| Arg <sup>8</sup> -Vasopressin ELISA Kit  | ADI-900-017<br>ADI-901-017   | 1x96 wells<br>5x96 wells | 3.39 pg/mL (4.10-1,000 pg/mL)  | Species Independent | CS, P, S, U                 | Overnight + 1 hour |
| ACTH ELISA Kit                           | ENZ-KIT138                   | 1x96 wells               | 0.46 pg/mL (0.5-165 pg/mL)     | H                   | P                           | 4.5 Hours          |
| Aldosterone ELISA Kit                    | ADI-900-173<br>ADI-901-173   | 1x96 wells<br>5x96 wells | 4.7 pg/mL (3.9-250 pg/mL)      | Species Independent | P, S, U                     | Overnight + 1 hour |
| Corticosterone ELISA Kit                 | ADI-900-097<br>ADI-901-097   | 1x96 wells<br>5x96 wells | 27.0 pg/mL (32-20,000 pg/mL)   | Species Independent | CL, CS, F, P, S, SA, U, WhB | 3 Hours            |
| Cortisol ELISA Kit                       | ADI-900-071<br>ADI-901-071   | 1x96 wells<br>5x96 wells | 56.72 pg/mL (156-10,000 pg/mL) | Species Independent | CS, F, P, S, SA, U          | 3 Hours            |
| DHEA ELISA Kit                           | ADI-900-093<br>ADI-901-093   | 1x96 wells<br>5x96 wells | 2.9 pg/mL (12.21-50,000 pg/mL) | Species Independent | CS, P, S, SA, U             | 5 Hours            |
| Estriol ELISA Kit                        | ADI-900-100                  | 1x96 wells               | 59.6 pg/mL (122-500,000 pg/mL) | Species Independent | CS, P, S, SA, U             | 3 Hours            |
| FSH ELISA Kit                            | ENZ-KIT108                   | 1x96 wells               | 0.5 mIU/mL (0.78-100 mIU/mL)   | H                   | P, S, TC                    | 2 Hours            |
| Histamine ELISA Kit                      | ENZ-KIT140                   | 1x96 wells               | 0.03 ng/mL (0.098-25 ng/mL)    | H, M, R, C          | P, TC, U                    | 2 Hours            |
| LH ELISA Kit                             | ENZ-KIT107                   | 1x96 wells               | 5.2 mIU/mL (1.2-280 mIU/mL)    | H, R                | P, S, TC                    | 2 Hours            |
| Melatonin ELISA Kit                      | ENZ-KIT150                   | 1x96 wells               | 0.162 ng/mL                    | Species Independent | S, SA, P, Fruit             | 2 Hours            |
| Oxytocin ELISA Kit                       | ADI-900-153A<br>ADI-901-153A | 1x96 wells<br>5x96 wells | 15 pg/mL (15.6-1,000 pg/mL)    | Species Independent | CS, CSF, M, P, S, SA, T, U  | Overnight + 1 hour |
| Progesterone ELISA Kit                   | ADI-900-011<br>ADI-901-011   | 1x96 wells<br>5x96 wells | 8.57 pg/mL (15.62-500 pg/mL)   | Species Independent | CS, P, S, SA                | < 3 Hours          |
| Prolactin ELISA Kit                      | ENZ-KIT161                   | 1x96 wells               | 19.2 pg/mL (31.3-2000 pg/mL)   | H                   | P, S                        | 3 Hours            |
| Serotonin ELISA Kit                      | ADI-900-175                  | 1x96 wells               | 0.293 ng/mL (0.49-500 ng/mL)   | Species Independent | CS, Platelets, P, S, U      | 3 Hours            |
| Testosterone ELISA Kit                   | ADI-900-065<br>ADI-901-065   | 1x96 wells<br>5x96 wells | 5.67 pg/mL (7.81-2,000 pg/mL)  | Species Independent | CS, F, P, S, SA, T          | 3 Hours            |
| Testosterone High Sensitivity ELISA Kit  | ADI-900-176                  | 1x96 wells               | 2.6 pg/mL (3.9-1,000 pg/mL)    | H                   | P, S, U                     | 3 Hours            |
| TSH ELISA Kit                            | ENZ-KIT131                   | 1x96 wells               | 0.5 µIU/mL (0.5-40 µIU/mL)     | H                   | S                           | 1.25 Hours         |

### LEGEND

Species: B= bovine, C = canine (dog), H = human, M = mouse, MK = monkey, P = porcine (pig), R= rat, RB = rabbit

Sample Types: AM = amniotic fluid, CL = cell lysates, CSF = cerebral spinal fluid, CS = culture supernatants, F = feces, M = milk, NL = nasal lavage, P = plasma, S = serum, SA = saliva, SS= sputum supernatant, T = tissue, TC = tissue culture, U = urine, WhB = white blood cells. Product tables includes both certified and cited species and sample types.

# EPIGENETICS RESEARCH

Epigenetics can be defined as the structural adaptation of chromosomal regions so as to register, signal or perpetuate altered activity states. A variety of chromatin-altering post-translational modifications (PTMs) to histone proteins including acetylation, methylation, and ubiquitinylation, as well as direct modifications to DNA, are known to turn gene transcription off or on. To date, acetylation/deacetylation of histone lysine residues, and methylation of DNA have proven to be of greatest clinical significance amongst epigenetic changes. Enzo's epigenetics kits highlighted below focus on DNA methylation.

## Convenient Kits for Sample Conversion and Detection


DNA methylation can alter gene expression, cell differentiation and often results in unidirectional changes to genomic DNA. Its involvement in so many cellular processes may lead to errors resulting in abnormal DNA methylation patterns which can lead to disease. Expanding our expertise in epigenetic modification analysis, Enzo Life Sciences now offers a portfolio of products that enable detection of DNA methylation which includes 5-Methylcytosine and 5-Hydroxymethylcytosine DNA ELISA kits.

## 5-METHYLCYTOSINE DNA ELISA KIT (ADI-900-224)

5-methylcytosine is a key component in the regulation of gene transcription, formed by the action of DNA methyltransferases. In epigenetics, the expression of methylated genes is altered due to the presence of 5-methylcytosine in methylated DNA. Quantification of 5-methylcytosine is critically important for the determination of epigenetic effects. The 5-Methylcytosine DNA ELISA kit is a complete, colorimetric, indirect immunoassay kit for the quantitative determination of 5-methylcytosine in DNA with results in less than 3 hours.

- Accurately quantitate 5-mC in any DNA sample in < 3 hours
- Ideal for high-throughput analysis
- High specificity comparable to LC-MS/MS-MRM analysis

## ELISA Results Closely Correlate to Mass Spectrometry Analysis



Genomic DNA Mass Spectrometry versus ELISA analysis: The 5-methylcytosine DNA ELISA kit (ADI-900-224) quantifies 5-mC in numerous DNA samples with close correlation to LC-MS/MS-MRM analysis. Genomic DNA samples include: human brain (HB), human kidney (HK), human embryonic stem cell (HESC), mouse brain (MB), mouse kidney (MK), and mouse testes (MT).

## SELECT ELISA KITS FOR EPIGENETIC RESEARCH

| Product Name                          | Product #    | Size       | Sensitivity (Range)   | Species             | Sample Types | Assay Time |
|---------------------------------------|--------------|------------|---|---------------------|--------------|------------|
| 5-Methylcytosine DNA ELISA Kit        | ADI-900-224A | 1x96 wells | ~0.5% 5-methylcytosine per 100 ng single-stranded DNA (5-100% (100 ng/μL))        | Species Independent | DNA          | < 3 Hours  |
| 5-Hydroxymethylcytosine DNA ELISA Kit | ADI-900-225  | 1x96 wells | < 0.02% 5-hydroxymethylcytosine DNA per 100 ng input DNA (0.03-0.55% (100 ng/μL)) | Species Independent | DNA          | ~3 Hours   |

### LEGEND

Species: B= bovine, C = canine (dog), H = human, M = mouse, MK = monkey, P = porcine (pig), R= rat, RB = rabbit

Sample Types: AM = amniotic fluid, CL = cell lysates, CSF = cerebral spinal fluid, CS = culture supernatants, F = feces, M = milk, NL = nasal lavage, P = plasma, S = serum, SA = saliva, SS= sputum supernatant, T = tissue, TC = tissue culture, U = urine, WhB = white blood cells. Product tables includes both certified and cited species and sample types.

# IMMUNITY/INFLAMMATION SIGNALING RESEARCH

Inflammation is one of the first responses of the immune system to infection or harmful stimuli, such as pathogens, damaged cells or irritants. It is a beneficial host response that protects tissues and involves a complex biological cascade of molecular and cellular signals that alter physiological responses, ultimately resulting in the familiar clinical symptoms. In the end, the body uses its vascular and immune system to try and remove the injurious stimuli and initiate the healing process. When inflammation is left unchecked, chronic inflammation can lead to various diseases, such as Alzheimer's disease, atherosclerosis, autoimmune disorders, cancer and others. There are many intracellular regulators of inflammation, one of which is the transcription factor NF- $\kappa$ B which upregulates the expression of cytokines, eicosanoids and adhesion molecules. Other inflammation regulators include the complement system of proteins which are part of the innate immunity system, blood pressure regulators (e.g. vasoconstrictors like ET-1, vasodilators like nitric oxide), coagulation regulators (e.g. haptoglobin, fibrinogen), and general inflammation regulators (e.g. cytokines like IL-6, IL-8, and TNF- $\alpha$ ; eicosanoids like prostaglandins). Overall, Enzo provides a plethora of products for the area of immunity and inflammation signaling research.

## MBL OLIGOMER ELISA KIT (BPD-KIT-029)


MBL (Mannose-Binding Lectin) is a carbohydrate-binding protein produced in the liver and secreted into the blood. The protein plays a key role in innate immunity, functioning as a player in the complement activation cascade and defending the body against invading microorganisms (bacteria, viruses, protozoa and fungi). MBL deficiency is common and associated with increased susceptibility to infections. Due to the high genetic variation in the MBL gene, at least 12% of the average Caucasian population has insufficient levels of functional MBL. Diagnosing MBL-deficiency is possible, but unfortunately no direct treatment exists, and it is difficult to manage. MBL evaluations are relevant for many patient groups, including children with recurrent infections and adults with recurrent and severe infections. Other relevant patient groups include women who have experienced recurrent spontaneous abortions; or immunocompromised patients undergoing cancer chemotherapy or immunosuppression (transplant), or those with an autoimmune disease. The main supportive therapeutics for these patients include regular antibiotics, antivirals and antifungals. For patients with combined immunodeficiencies or recurrent spontaneous abortions, some doctors prefer IV immunoglobulin treatment. To determine low levels of MBL in MBL-deficient people, Enzo Life Sciences offers a sensitive MBL oligomer ELISA Kit to quantify the oligomerized MBL in human serum and plasma.

- Detect low levels of oligomerized MBL in deficient persons
- Ready-to-use, pre-coated plate and individual calibrators save time and reduce errors
- Analyze 40 samples in duplicate in just 4 hours

### Research Areas

- Cancer
- Cardiovascular
- Dengue
- Diabetes
- Hepatitis B
- HIV
- Immunology
- Pulmonary Fibrosis
- Renal Disease
- Rheumatoid Arthritis

### Detect Low Levels of MBL Functional Form in MBL-deficient People



MBL exists at different levels of polymerization — Non-functional and functional forms. A: a single chain monomer; B: 3 monomers joined in a collagen helix forming a subunit “flower”; or C: up to 6 “flowers” joined together to form a hexameric structure comprising a total of 18 MBL monomers. The MBL oligomer ELISA Kit (BPD-KIT-029) allows you to measure the functional, oligomeric form of MBL.

## SELECT ELISA KITS FOR IMMUNITY/INFLAMMATION RESEARCH

| Product Name  | Product #   | Size       | Sensitivity (Range)          | Species            | Sample Types | Assay Time          |
|---|-------------|------------|------------------------------|--------------------|--------------|---------------------|
| BAFF Soluble (mouse) Detection Set                  | APO-54N-013 | 4x96 wells | 0.3 ng/mL (0-20 ng/mL)       | M                  | CS, S        | 2 Hours             |
| CD14, Soluble (human) ELISA Kit                     | ALX-850-302 | 1x96 wells | 5-50 ng/mL                   | H                  | P, S         | < 2.5 Hours         |
| CD14, Soluble (mouse) ELISA Kit                     | ALX-850-303 | 1x96 wells | 5-50 ng/mL                   | M                  | P, S         | < 2 Hours           |
| CD40, Soluble (human) ELISA Kit                     | ALX-850-262 | 1x96 wells | 6.92 pg/mL (7.8-500 pg/mL)   | H                  | CS, S        | < 3.25 Hours        |
| CD40L, Soluble (human) ELISA Kit (High Sensitivity) | ALX-850-311 | 1x96 wells | 0.005 ng/mL (0.08-5 ng/mL)   | H                  | CS, S        | < 3.5 Hours         |
| sCD44std, Soluble (human) ELISA Kit                 | ALX-850-053 | 1x96 wells | 0.015 ng/mL (0.12-4 ng/mL)   | H                  | CS, S        | < 3.5 Hours         |
| Complement C3a des Arg (human) ELISA Kit            | ADI-900-058 | 1x96 wells | 0.120 ng/mL (0.313-20 ng/mL) | H                  | P            | 3 Hours             |
| Complement C4a des Arg (human) ELISA Kit            | ADI-900-059 | 1x96 wells | 0.76 ng/mL (0.78-200 ng/mL)  | H                  | P            | 3 Hours             |
| COX-2 (human) ELISA Kit                             | ADI-900-094 | 1x96 wells | 0.25 ng/mL (1.09-70 ng/mL)   | H                  | CL           | 2 Hours             |
| CRP (human) ELISA Kit                               | ENZ-KIT102  | 1x96 wells | 12.685 ng/mL (10-810 ng/mL)  | H                  | P, S         | < 2 Hours           |
| FasL, Soluble (human) ELISA Kit                     | ALX-850-246 | 1x96 wells | 0.07 ng/mL (0.16-10 ng/mL)   | H                  | CS, S        | < 3.5 Hours         |
| GRO $\alpha$ /CINC-1 (rat) ELISA Kit                | ADI-900-074 | 1x96 wells | 1.99 pg/mL (4.7-300 pg/mL)   | R                  | CS, P, S     | 2 Hours             |
| Histamine ELISA Kit                                 | ENZ-KIT140  | 1x96 wells | 0.03 ng/mL (0.098-25 ng/mL)  | H, M, R, C         | P, S, TC, U  | 2 Hours             |
| sHLA-G ELISA Kit                                    | ALX-850-309 | 1x96 wells | 1.0 U/mL (1.95-125 U/mL)     | H                  | CS, P, S     | ~20 Hours           |
| IgG1 (mouse) ELISA Kit                              | ADI-900-109 | 1x96 wells | 0.064 ng/mL (7.81-250 ng/mL) | M                  | CS, S        | 1.5 Hours           |
| IgG2a (mouse) ELISA Kit                             | ADI-900-113 | 1x96 wells | 318.8 pg/mL (7.81-250 ng/mL) | M                  | CS, S        | 1.5 Hours           |
| IgG2b (mouse) ELISA Kit                             | ADI-900-110 | 1x96 wells | 1.05 ng/mL (7.81-250 ng/mL)  | M                  | CS, S        | 1.5 Hours           |
| IgM (mouse) ELISA Kit                               | ADI-900-120 | 1x96 wells | 0.60 ng/mL (3.91-250 ng/mL)  | M                  | AS, CS, S    | 1.5 Hours           |
| LBP, Soluble (mouse) ELISA Kit                      | ALX-850-305 | 1x96 wells | 5-50 ng/mL                   | M, R               | P, S         | < 2.5 Hours         |
| LBP, Soluble ELISA Kit                              | ALX-850-304 | 1x96 wells | 5-50 ng/mL                   | H, B, C, RB, P, HR | P, S         | < 2.5 Hours         |
| MBL Oligomer ELISA Kit                              | BPD-KIT-029 | 1x96 wells | 2 ng/mL (0.50-40 ng/mL)      | H                  | P, S         | < 4 Hours           |
| NF $\kappa$ B p65 ELISA Kit (chemiluminescent)      | ADI-EKS-446 | 2x96 wells | Not available                | H, M, R            | CL           | 3 Hours             |
| PMN-Elastase (human) ELISA Kit                      | ALX-850-265 | 1x96 wells | 1.98 ng/mL (0.16-10 ng/mL)   | H                  | CS, P        | < 2.5 Hours         |
| Progranulin (human) ELISA Kit                       | ALX-850-376 | 1x96 wells | 18 pg/mL (75-2,500 pg/mL)    | H                  | P, S         | 2.5 Hours           |
| PTX3 (human) Detection Set                          | ALX-850-299 | 1 Set      | 75 pg/mL                     | H, B               | P, S         | ~5 Hours            |
| RANKL (total) Soluble (human) ELISA Kit             | ALX-850-019 | 1x96 wells | ~1.56 pg/mL (2.2-60 pmol/L)  | H                  | CS, P, S     | Overnight + 3 Hours |

### LEGEND

Species: B= bovine, C = canine (dog), H = human, M = mouse, MK = monkey, P = porcine (pig), R= rat, RB = rabbit

Sample Types: AM = amniotic fluid, CL = cell lysates, CSF = cerebral spinal fluid, CS = culture supernatants, F = feces, M = milk, NL = nasal lavage, P = plasma, S = serum, SA = saliva, SS= sputum supernatant, T = tissue, TC = tissue culture, U = urine, WhB = white blood cells. Product tables includes both certified and cited species and sample types.



# CYTOKINE RESEARCH

Cytokines, including interleukins, interferons, and chemokines, are secreted signaling molecules that mediate and regulate immunity, inflammation and some developmental processes during embryogenesis. They also play an important role in angiogenesis and cancer. Cytokines can have pleiotropic, overlapping, and sometimes contradictory functions. This is dependent upon their concentration, the cell type they are acting on, and the presence of other cytokines and mediators. Enzo Life Sciences offers a wide selection of cytokine antibodies, proteins, and kits for innate and adaptive immune research including our wide selection of high-specificity human and mouse cytokine and interleukin ELISA kits.

| CYTOKINE/INTERLEUKIN ELISA KITS |          |                    |               |                    |
|---------------------------------|----------|--------------------|---------------|--------------------|
| Human                           |          |                    | Mouse         |                    |
| Adiponectin                     | IL-4     | Leptin             | IFN- $\gamma$ | Leptin             |
| CD40                            | IL-6     | Lipocalin 2 (NGAL) | IL-1 $\beta$  | Lipocalin 2 (NGAL) |
| CD40L                           | IL-8     | Nampt              | IL-2          | Nampt              |
| FasL                            | IL-10    | Osteoprotegerin    | IL-4          | TGF- $\beta$ 1     |
| IFN- $\gamma$                   | IL-12p70 | RANKL              | IL-6          | TNF- $\alpha$      |
| IGF-1                           | IL-13    | TNF- $\alpha$      |               |                    |
| IL-1 $\beta$                    | IL-17A   | TNF-R1             |               |                    |
| IL-2                            | IL-33    |                    |               |                    |


## IL-6 (HUMAN) HIGH SENSITIVITY ELISA KIT (ENZ-KIT178)

Interleukin-6 (IL-6) is a cytokine critical to the regulation of the immune and hematopoietic systems. Human IL-6 is expressed by T-cells, mast cells, monocytes, macrophages, fibroblasts, endothelial cells, keratinocytes, and many tumor cell lines. It plays a role in acute phase reactions and response to injury and inflammation. Human IL-6 also stimulates the differentiation of B-cells for antibody production, promotes expansion of activated T-cells, expands hematopoietic cell production, and induces the expression of acute phase proteins.

We are proud to provide a comprehensive portfolio of products to detect and modulate the complex processes of immunological disease. Our IL-6 High Sensitivity ELISA Kit can detect concentrations as low as 0.057 pg/mL of IL-6 in culture supernatants, plasma, serum, and urine.

- High-sensitivity ELISA enables quantification of low amounts of IL-6
- Low cross-reactivity with other interleukins and similar proteins
- Easy-to-use liquid color-coded reagents and pre-coated plates reduce errors and save time
- Rapid quantitative results in 3 hours
- High-throughput immunoassay measures up to 40 samples in duplicate per 96-well kit

### Highly Sensitive Measurement of IL-6



Standard curve for IL-6 High-Sensitivity ELISA Kit (ENZ-KIT178).

## SELECT CYTOKINE ELISA KITS

| Product Name                            | Product #                                 | Size                                    | Sensitivity (Range)   | Species       | Sample Types | Assay Time |
|---|---|---|---|---------------|--------------|------------|
| IFN- $\gamma$ (human) ELISA Kit         | ADI-900-136                               | 1x96 wells                              | < 2 pg/mL (25.6-1,000 pg/mL)  | H             | CS, P, S, U  | 3 Hours    |
| IFN- $\gamma$ (mouse) ELISA Kit         | ADI-900-137                               | 1x96 wells                              | < 10 pg/mL (37-3,000 pg/mL)   | M             | CS, P, S, U  | 4 Hours    |
| IL-10 (human) ELISA Kit                 | ADI-900-036<br>ADI-901-036                | 1x96 wells<br>5x96 wells                | 3.75 pg/mL (7.81-500 pg/mL)   | H             | CS, P, S     | < 3 Hours  |
| IL-10 (mouse) ELISA Kit                 | ADI-900-148                               | 1x96 wells                              | < 12 pg/mL (37-3,000 pg/mL)   | M             | CS, P, S     | 5 Hours    |
| IL-12p70 (human) ELISA Kit              | ADI-900-202                               | 1x96 wells                              | 0.9 pg/mL (7.81-500 pg/mL)  | H             | CS, P, S     | 3 Hours    |
| IL-13 (human) ELISA Kit                 | ADI-900-208                               | 1x96 wells                              | 1.71 pg/mL (1.56-100 pg/mL)   | H             | CS, P, S     | 2 Hours    |
| IL-17A (human) ELISA Kit                | ADI-900-177                               | 1x96 wells                              | 0.201 pg/mL (2.34-75 pg/mL)   | H             | CS, P, S     | 3 Hours    |
| IL-1 $\beta$ (human) ELISA Kit          | ADI-900-130                               | 1x96 wells                              | < 1 pg/mL (10.24-400 pg/mL)   | H             | CS, P, S, U  | 4 Hours    |
| IL-1 $\beta$ (mouse) ELISA Kit          | ADI-900-132A                              | 1x96 wells                              | < 3 pg/mL (31.25-1,000 pg/mL)   | M             | CS, P, S     | 4 Hours    |
| IL-1 $\beta$ (rat) ELISA Kit            | ADI-900-131                               | 1x96 wells                              | < 12 pg/mL (25.6-2,500 pg/mL)   | R             | CS, P, S, T  | < 4 Hours  |
| IL-2 (human) ELISA Kit                  | ADI-900-118A                              | 1x96 wells                              | 6.6 pg/mL (7.81-500 pg/mL)  | H             | CS, P, S     | 3 Hours    |
| IL-2 (mouse) ELISA Kit                  | ADI-900-042                               | 1x96 wells                              | 3.12 pg/mL (7.81-1,000 pg/mL)   | M             | CS, S        | < 4 Hours  |
| IL-33 (human) ELISA Kit                 | ADI-900-201                               | 1x96 wells                              | 0.60 ng/mL (3.91-250 ng/mL)   | M             | AS, CS, S    | 90 Minutes |
| IL-33 soluble (human) Detection Set     | APO-54N-025                               | 5x96 wells                              | 0.60 ng/mL (3.91-250 ng/mL)   | H, M          | CS, S        | 2.5 Hours  |
| IL-4 (human) ELISA Kit                  | ADI-900-145A                              | 1x96 wells                              | < 2 pg/mL (10.24-400 pg/mL)   | H             | CS, P, S     | 3.5 Hours  |
| IL-4 (mouse) ELISA Kit                  | ADI-900-043                               | 1x96 wells                              | 4.34 pg/mL (7.81-1,000 pg/mL)   | M             | CS, P, S     | 3 Hours    |
| IL-5 (human) ELISA Kit                  | ENZ-KIT139                                | 1x96 wells                              | < 2 pg/mL (7.8-500 pg/mL)   | H             | CS, P, S, U  | 3 Hours    |
| IL-6 (human) ELISA Kit                  | ADI-900-033<br>ADI-901-033                | 1x96 wells<br>5x96 wells                | 6.01 pg/mL (7.81-500 pg/mL)   | H             | CS, P, S, U  | < 3 Hours  |
| IL-6 (human) High Sensitivity ELISA Kit | ENZ-KIT178                                | 1x96 wells                              | 0.057 pg/mL   | H             | CS, P, S, U  | 3 Hours    |
| IL-6 (mouse) ELISA Kit                  | ADI-900-045                               | 1x96 wells                              | 1.01 pg/mL (7.81-1,000 pg/mL)   | M             | CS, S, T     | < 3 Hours  |
| IL-8 (human) ELISA Kit                  | ADI-900-156<br>ADI-901-156<br>ADI-902-156 | 1x96 wells<br>5x96 wells<br>10x96 wells | 0.64 pg/mL (7.8-1,000 pg/mL)  | H             | CS, P, S     | < 3 Hours  |
| MCP-1 (rat) ELISA Kit                   | ADI-900-077                               | 1x96 wells                              | 20.45 pg/mL (50-3,200 pg/mL)  | R             | CS, S        | 2 Hours    |
| Osteoprotegerin (human) ELISA Kit       | ALX-850-280A                              | 1x96 wells                              | 1.4 pg/mL (0-400 pg/mL)   | H             | P, S         | ~5 Hours   |
| TGF- $\beta$ 1 ELISA Kit                | ADI-900-155                               | 1x96 wells                              | 3.3 pg/mL in Assay Buffer 13;<br>10.8 pg/mL in Assay Buffer 30<br>(31.25-1000 pg/mL)                          | H, M,<br>R, B | CS, P, S     | 4 Hours    |
| TL1A soluble (human) Detection Set      | APO-54N-024                               | 5x96 wells                              | 20 pg/mL (0-25 ng/mL)   | H             | CS, S        | < 3 Hours  |
| TL1A soluble (human) ELISA Kit          | APO-54N-027                               | 1x96 wells                              | 15 pg/mL (39-2,500 pg/mL)   | H             | CS, S        | < 3 Hours  |
| TNF-R1 (Soluble) (human) ELISA Kit      | ALX-850-047                               | 1x96 wells                              | 53 pg/mL (0.08-5 ng/mL)   | H             | CS, S        | 2.5 Hours  |
| TNF- $\alpha$ (human) ELISA Kit         | ADI-900-099<br>ADI-901-099                | 1x96 wells<br>5x96 wells                | 8.43 pg/mL (15.63-1,000 pg/<br>mL)  | H             | CS, P, S     | 4 Hours    |
| TNF- $\alpha$ (mouse) ELISA kit         | ADI-900-047                               | 1x96 wells                              | 3.9 pg/mL (15.63-2,000 pg/mL)   | M             | CS, P, S, T  | 5 Hours    |
| TNF- $\alpha$ (rat) ELISA kit           | ADI-900-086A                              | 1x96 wells                              | 12.0 pg/mL in Assay Buffer;<br>26.7 pg/mL in culture media<br>supplemented with 10% FBS<br>(31.3-2,000 pg/mL) | R             | CS           | 3 Hours    |

### LEGEND

Species: B= bovine, C = canine (dog), H = human, M = mouse, MK = monkey, P = porcine (pig), R= rat, RB = rabbit

Sample Types: AM = amniotic fluid, CL = cell lysates, CSF = cerebral spinal fluid, CS = culture supernatants, F = feces, M = milk, NL = nasal lavage, P = plasma, S = serum, SA = saliva, SS= sputum supernatant, T = tissue, TC = tissue culture, U = urine, WhB = white blood cells. Product tables includes both certified and cited species and sample types.

# EICOSANOIDS RESEARCH


When inflammation is left unchecked, it can contribute to the pathogenesis of numerous diseases which include rheumatoid arthritis, chronic bronchitis, emphysema, asthma, psoriasis, cancer, and colitis. Molecules in the arachidonic acid enzyme cascade play a number of important biological roles, both normal and pathological. The derivatives of arachidonic acid are known as eicosanoids, a group of biologically active oxygenated unsaturated fatty acids. They include prostaglandins, thromboxanes, leukotrienes, hydroxyeicosatetraenoic acids (HETEs), and lipoxins. Eicosanoids have been shown to enhance, as well as attenuate inflammation and have also been linked to carcinogenesis. Enzo Life Sciences provides a wide selection of high sensitivity ELISA kits for detection of eicosanoids, as well as small molecules, antibodies and a bioactive lipid library for studying the role of eicosanoids in inflammation and immunity.

## PROSTAGLANDIN E<sub>2</sub> ELISA KITS (ADI-900-001)

PGE<sub>2</sub> is an hormone-like lipophilic eicosanoid important in inducing fever, causing uterine contractions during labor, and stimulating osteoblasts to release factors that cause bone reabsorption. Enzo Life Sciences offers the most sensitive and complete colorimetric ELISA kits for quantification of prostaglandins in a wide variety of sample types.

- Ultra-sensitive colorimetric ELISAs to measure as little as 8.26 pg/mL PGE<sub>2</sub>
- Widely cited in peer-reviewed literature
- Available to use for cell lysates, culture supernatants, serum, saliva, urine, and many more sample types
- High-throughput capabilities with chemiluminescent and fluorescent format options

### Consistent Lot-to-Lot Performance for Over 5 Years




Lot-to-lot consistency: Graph demonstrates the robust and reproducible nature of the PGE<sub>2</sub> ELISA Kit (ADI-900-001) showing standard curves from 10 lots manufactured over 5 years.

## LTB<sub>4</sub> ELISA KIT (ADI-900-068)

Leukotrienes are major products of 5-lipoxygenase metabolism of arachidonic acid. Leukotriene B<sub>4</sub> (LTB<sub>4</sub>) stimulates leukocyte functions including lysosomal enzyme release, adhesion, and aggregation of polymorphonuclear leukocytes. This eicosanoid has also been implicated as a potent mediator of inflammatory diseases and immunoregulation.

- Detect as low as 5.63 pg/mL of LTB<sub>4</sub>
- Broad dynamic range suitable for a large variety of samples
- Rapidly assay up to 39 samples in duplicate in just 4 hours
- Reproducible results day-after-day and lot-after-lot

### Ultra-sensitive Detection of Leukotriene B<sub>4</sub>



Standard curve experiment using the LTB<sub>4</sub> ELISA Kit (ADI-900-068).

## SELECT EICOSANOID ELISA KITS

| Product Name  | Product #                  | Size                     | Sensitivity (Range)  | Species             | Sample Types  | Assay Time                          |
|---|----------------------------|--------------------------|--|---------------------|---|-------------------------------------|
| 12(S)-HETE ELISA Kit                                    | ADI-900-050<br>ADI-901-050 | 1x96 wells<br>5x96 wells | 146 pg/mL (195-50,000 pg/mL)   | Species Independent | CS, renal interstitial fluid, P, S, T   | 5 Hours                             |
| 15(S)-HETE ELISA Kit                                    | ADI-900-051<br>ADI-901-051 | 1x96 wells<br>5x96 wells | 69.21 pg/mL (78.1-20,000 pg/mL)  | Species Independent | CS, P, S, U   | Overnight + 1 Hour                  |
| 13(S)-HODE ELISA Kit                                    | ADI-900-108                | 1x96 wells               | 1.6 ng/mL (3.9-1,000 ng/mL)  | Species Independent | CS, CL, P, S, SA, T, U  | 4 Hours                             |
| LTB <sub>4</sub> ELISA Kit                              | ADI-900-068<br>ADI-901-068 | 1x96 wells<br>5x96 wells | 5.63 pg/mL (11.7-3,000 pg/mL)  | Species Independent | Bronchioalveolar Lavage Fluid, CL, CS, Crevicular Fluid, P, SA, T, U, WhB                 | 4 Hours                             |
| Cysteinyl leukotriene ELISA Kit                         | ADI-900-070<br>ADI-901-070 | 1x96 wells<br>5x96 wells | 26.6 pg/mL (78.1-2,500 pg/mL)  | Species Independent | CS, U   | 4 Hours                             |
| PGE <sub>1</sub> ELISA Kit                              | ADI-900-005<br>ADI-901-005 | 1x96 wells<br>5x96 wells | 5.58 pg/mL (4.88-5,000 pg/mL)  | Species Independent | CS, Simulated Interstitial Lung Fluid, P, S, SA, U  | < 3 Hours                           |
| PGE <sub>2</sub> CLIA Kit                               | ADI-910-001                | 1x96 wells               | 6.03 pg/mL (7.81-1,000 pg/mL)  | Species Independent | CS, S, SA, U, WhB   | 3 Hours                             |
| PGE <sub>2</sub> ELISA Kit                              | ADI-900-001<br>ADI-901-001 | 1x96 wells<br>5x96 wells | 13.4 pg/mL (39.1-2,500 pg/mL)  | Species Independent | CL, CS, CSF, Dialysate, Gingival Crevicular Fluid, P, S, SA, T, U, WhB                    | < 3 Hours                           |
| PGE <sub>2</sub> FPIA Kit                               | ADI-920-001                | 100 Tests                | 684 pg/mL (1,562-100,000 pg/mL)  | Species Independent | CS  | 30 Minutes                          |
| PGE <sub>2</sub> High-Sensitivity ELISA Kit             | ADI-930-001<br>ADI-931-001 | 1x96 wells<br>5x96 wells | 8.26 pg/mL (7.8-1,000 pg/mL)   | Species Independent | CS, P, S, SA, U, WhB  | Overnight + 1 hour                  |
| 6-keto-PGF <sub>1α</sub> ELISA Kit                      | ADI-900-004<br>ADI-901-004 | 1x96 wells<br>5x96 wells | 1.40 pg/mL (3.2-50,000 pg/mL)  | Species Independent | CS, S, SA, U  | < 3 Hours                           |
| PGF <sub>2α</sub> ELISA Kit                             | ADI-900-069<br>ADI-901-069 | 1x96 wells<br>5x96 wells | 6.71 pg/mL (3.05-50,000 pg/mL)   | Species Independent | CL, CS, P, S, SA, U, WhB  | < 3 Hours                           |
| PGF <sub>2α</sub> High-Sensitivity ELISA Kit            | ADI-930-069<br>ADI-931-069 | 1x96 wells<br>5x96 wells | 0.98 pg/mL (1.95-2,000 pg/mL)  | Species Independent | CS, M, P, S, SA, U  | Overnight + 3 Hours                 |
| 8-iso-PGF <sub>2α</sub> ELISA Kit                       | ADI-900-010<br>ADI-901-010 | 1x96 wells<br>5x96 wells | 16.3 pg/mL (6.1-100,000 pg/mL)   | Species Independent | CS, P, T, U   | < 3 Hours                           |
| Direct 8-iso-PGF <sub>2α</sub> ELISA Kit                | ADI-900-091<br>ADI-901-091 | 1x96 wells<br>5x96 wells | (2 hour) 103.2 pg/mL;<br>(overnight) 40.0 pg/mL<br>(160-100,000 pg/mL) | Species Independent | CL, P, S, T, U  | < 3 Hours or Overnight + 45 Minutes |
| 15-deoxy-Δ <sup>12,14</sup> -PGJ <sub>2</sub> ELISA Kit | ADI-900-023<br>ADI-901-023 | 1x96 wells<br>5x96 wells | 36.8 pg/mL (195-200,000 pg/mL)   | Species Independent | CS, P, SA, U  | 5 Hours                             |
| TXB <sub>2</sub> ELISA Kit                              | ADI-900-002<br>ADI-901-002 | 1x96 wells<br>5x96 wells | 10.54 pg/mL (13.7-10,000 pg/mL)  | Species Independent | CL, Coronary Effluent, CS, Liver Perfusate, Platelets, P, Rectal Dialysate, S, SA, U, WhB | < 3 Hours                           |
| 11-dehydro-TXB <sub>2</sub> ELISA Kit                   | ADI-900-092<br>ADI-901-092 | 1x96 wells<br>5x96 wells | 4.31 pg/mL (9.8-10,000 pg/mL)  | Species Independent | CS, U, S  | 3 Hours                             |
| Urinary Prostacyclin ELISA Kit                          | ADI-900-025<br>ADI-901-025 | 1x96 wells<br>5x96 wells | 6.58 pg/mL (7.81-2,000 pg/mL)  | Species Independent | CS, U   | < 3 Hours                           |

### LEGEND

Species: B= bovine, C = canine (dog), H = human, M = mouse, MK = monkey, P = porcine (pig), R= rat, RB = rabbit

Sample Types: AM = amniotic fluid, CL = cell lysates, CSF = cerebral spinal fluid, CS = culture supernatants, F = feces, M = milk, NL = nasal lavage, P = plasma, S = serum, SA = saliva, SS= sputum supernatant, T = tissue, TC = tissue culture, U = urine, WhB = white blood cells. Product tables includes both certified and cited species and sample types.

# METABOLISM RESEARCH


Metabolism includes enzyme-catalyzed reactions within cells that allow the organism to maintain their structure, grow, reproduce, and respond to their environment. Metabolism also involves the breakdown of organic substances (e.g. proteins, lipids, carbohydrates), harvesting of energy by way of cellular respiration, and using that energy to construct cellular components (e.g. proteins, nucleic acids). When the metabolic system is unbalanced, this can lead to metabolic diseases such as obesity, diabetes, insulin resistance, hypertension, atherosclerosis, and cancer. Enzo Life Sciences' ELISAs for metabolism research include kits for the analysis of popular metabolic biomarkers such as kidney failure (e.g. NGAL, KIM-1), bone mineralization (e.g. vitamin D), obesity (e.g. leptin) and several others.

## 25(OH) VITAMIN D ELISA KIT (ADI-900-215)

Known as the “sunshine vitamin”, Vitamin D is responsible for metabolic intestinal absorption of calcium and phosphate, and is critical to the process of bone mineralization. Vitamin D deficiency is associated with a variety of diseases, including osteoporosis, rheumatoid arthritis, diabetes, and cancer. With a rapid assay time and easy-to-use protocol, Enzo's 25(OH) Vitamin D ELISA kit is the fastest, user-friendly assay on the market. Reduce your sample preparation and assay time without sacrificing sensitivity or reproducibility. Discover the benefits of this colorimetric, competitive immunoassay kit for quantifying 25(OH) Vitamin D<sub>2</sub> and D<sub>3</sub> in human plasma and serum samples.

- Measures as little as 1.98 ng/mL of 25(OH) Vitamin D<sub>2</sub> and D<sub>3</sub> in just 1.5 hours
- Easy-to-use protocol with rapid dissociation step reduces errors and bench time
- Convenient alternative to labor-intensive LC-MS method

### Detect Low Concentrations of Common Forms of Vitamin D (D<sub>2</sub> & D<sub>3</sub>)




Standard curve experiment using the 25(OH) Vitamin D ELISA Kit (ADI-900-215).

## INSULIN ELISA KIT (ENZ-KIT141)

Insulin is a peptide hormone made in the pancreas by beta cells and promotes the uptake of glucose from the blood into tissues, where it is stored in the form of glycogen and fat. By controlling glucose levels, insulin serves as the central regulator of fat and carbohydrate metabolism. The disease diabetes mellitus results when control of glucose levels is lost due to lack of insulin production (type 1 diabetes) or when sufferers develop an inability to respond to insulin (type 2 diabetes).

- Highly sensitive measurement of insulin, detecting as little as 21.6 pg/mL
- Negligible cross reactivity with proinsulin
- High-throughput format with results in 3 hours for up to 39 samples in duplicate
- Fully quantitative results that surpass semi-quantitative Western blot analysis

### Highly Sensitive Measurement of Insulin



Standard curve experiment using Insulin ELISA Kit (ENZ-KIT141).

## SELECT ELISA KITS FOR METABOLISM RESEARCH

| Product Name                                     | Product #     | Size       | Sensitivity (Range)  | Species | Sample Types   | Assay Time |
|--|---------------|------------|--|---------|----------------|------------|
| 25(OH) Vitamin D ELISA Kit                       | ADI-900-215   | 1x96 wells | 1.98 ng/mL (0.5-1,010 ng/mL)                                       | H       | P, S           | ~1.5 Hours |
| Adiponectin (human) ELISA Kit                    | ALX-850-377   | 1x96 wells | 0.156 ng/mL (CSF, urine); 0.47 ng/mL (plasma, serum) (1-150 ng/mL) | H       | CSF, P, S, U   | < 3 Hours  |
| AMP'D® GLP-1 ELISA Kit                           | ENZ-KIT104    | 1x96 wells | 5.5 pg/mL (range 7.8-250 pg/mL)                                    | H       | P, S           | 3 Hours    |
| ApoE ELISA Kit                                   | ENZ-KIT134    | 1x96 wells | 14.063 ng/mL (23.438-1,500 ng/mL)                                  | H       | CS, P, S, T    | 3.5 Hours  |
| Gastrin I (human) ELISA Kit                      | ADI-900-026   | 1x96 wells | 7.27 pg/mL (39.1-10,000 pg/mL)                                     | H       | CS, P, S       | 5 Hours    |
| Gastrin I (rat) ELISA Kit                        | ADI-900-149   | 1x96 wells | 78.1 pg/mL (78.1-5,000 pg/mL)                                      | R       | CS, P, S       | 5 Hours    |
| Insulin ELISA Kit                                | ENZ-KIT141    | 1x96 wells | 21.6 pg/mL (15.6-500 pg/mL)  | H, P    | P, S           | 3 Hours    |
| KIM-1 (human) ELISA Kit                          | ADI-900-226   | 1x96 wells | 1.279 pg/mL (7.813-500 pg/mL)                                      | H       | U              | < 2 Hours  |
| Leptin (human) ELISA Kit                         | ADI-900-028A  | 1x96 wells | 23.4 pg/mL (31.3-2,000 pg/mL)                                      | H       | CS, P, S       | 3 Hours    |
| Leptin (mouse) ELISA Kit                         | ADI-900-019A  | 1x96 wells | 25.4 pg/mL (50-3,200 pg/mL)  | M       | CS, P, S, T    | 3 Hours    |
| Leptin (rat) ELISA Kit                           | ADI-900-015A  | 1x96 wells | 67.2 pg/mL (100-6,400 pg/mL)                                       | R       | CS, P, S       | 3 Hours    |
| Nampt (Visfatin/PBEF) (human) ELISA Kit          | AG-45A-0006EK | 1x96 wells | 30 pg/mL (0.125-8 ng/mL)   | H       | S              | 2.5 Hours  |
| Nampt (Visfatin/PBEF) (mouse/rat) Dual ELISA Kit | AG-45A-0007EK | 1x96 wells | 50 pg/mL (0.5-32 ng/mL)  | M, R    | S              | < 4 Hours  |
| NGAL (dog) ELISA Kit                             | BPD-KIT-043   | 1x96 wells | 0.56 pg/mL (4-400 pg/mL)   | C       | CS, P, S, U, T | 4 Hours    |
| NGAL (human) ELISA Kit                           | BPD-KIT-036   | 1x96 wells | 4 pg/mL (10-1,000 pg/mL)   | H       | CS, P, S, U, T | < 4 Hours  |
| NGAL (monkey) ELISA Kit                          | BPD-KIT-045   | 1x96 wells | 1.5 pg/mL (10-200 pg/mL)   | MK      | CS, P, S, U, T | < 4 Hours  |
| NGAL (mouse) ELISA Kit                           | BPD-KIT-042   | 1x96 wells | 0.75 pg/mL (10-1,000 pg/mL)  | M       | CS, P, S, U, T | 4 Hours    |
| NGAL (pig) ELISA Kit                             | BPD-KIT-044   | 1x96 wells | 1 pg/mL (10-400 pg/mL)   | P       | CS, P, S, U, T | 4 Hours    |
| NGAL (rat) ELISA Kit                             | BPD-KIT-046   | 1x96 wells | 0.5 pg/mL (4-400 pg/mL)  | R       | CS, P, S, U, T | 4 Hours    |
| NGAL rapid (human) ELISA Kit                     | BPD-KIT-037   | 1x96 wells | < 0.1 ng/mL (0.2-20 ng/mL)   | H       | P, U           | 1 Hour     |
| Peptide YY ELISA Kit                             | ENZ-KIT133    | 1x96 wells | 18.75 pg/mL  | H       | CS, P, S, T    | 3.5 Hours  |
| Proinsulin ELISA Kit                             | ENZ-KIT149    | 1x96 wells | 0.17 pM (1.56-50 pM)   | H       | P, S           | 3 Hours    |

### LEGEND

Species: B= bovine, C = canine (dog), H = human, M = mouse, MK = monkey, P = porcine (pig), R= rat, RB = rabbit

Sample Types: AM = amniotic fluid, CL = cell lysates, CSF = cerebral spinal fluid, CS = culture supernatants, F = feces, M = milk, NL = nasal lavage, P = plasma, S = serum, SA = saliva, SS= sputum supernatant, T = tissue, TC = tissue culture, U = urine, WhB = white blood cells. Product tables includes both certified and cited species and sample types.



# NEPHROLOGY RESEARCH

Chronic kidney disease (CKD) is a major public health problem, as 10% of the worldwide population is affected by CKD. There have been numerous studies focused on this disease over the years but there remains a need for improved therapeutics and identifying and/or predicting patient outcomes. One protein that has been identified as playing an important role in kidney disease is Kidney Injury Molecule-1 (KIM-1). It is the most highly upregulated protein in the proximal tubule of the injured kidney. It exists in very low levels in normal kidneys but when upregulated during injury, it is detectable in urine in a wide variety of human diseases. KIM-1 is a potential biomarker for renal injury, which would suggest it has great importance in various kidney diseases and disorders, such as chronic kidney disease (as mentioned above), as well as acute tubular necrosis and acute kidney failure.


## KIM-1 (HUMAN) ELISA KIT (ADI-900-226)

Acute kidney injury (AKI), the precursor to acute renal failure (ARF) is the rapid loss of the kidney's ability to remove metabolic waste and help balance fluids and electrolytes in your body. This sudden loss of kidney function, which is recoverable, can be the result of illness, physical trauma, or nephrotoxic drugs. Numerous studies have identified a number of proteins and metabolites present in blood or urine with varied utility as biomarkers of these various stages of kidney injury, disease progression, or nephrotoxicity. Enzo provides a variety of tools to detect the early stages of acute kidney injury or renal failure, including sensitive ELISA kits for detection of KIM-1 and NGAL (lipocalin-2). These rapid, high-sensitivity kits are validated for a variety of species and samples types.

The KIM-1 (human) ELISA kit provides ultra-sensitive quantification enabling reduced input sample and matrix interference.

- Sensitive ELISA measures as little as 1.2 pg/mL of KIM-1 in urine
- Low cross-reactivity to TIM-3 and TIM-4
- Rapid assay with results in < 2 hours

### Rapidly Quantify KIM-1, an Early Biomarker of AKI




Parallelism analysis using the KIM-1 (human) ELISA Kit (ADI-900-226) indicates antigen binding characteristics are similar to native KIM-1 in urine samples with no matrix interference at the dilutions tested.

## BIOPORTO® NGAL ELISA KITS

Neutrophil gelatinase-associated lipocalin (NGAL, lipocalin-2, siderocalin) is a small, robust protein detectable in urine and blood within 2 hours upon kidney tubular damage.


- Sensitive (pg/mL) ELISAs validated for human samples and a broad range of model systems
- Rapid human NGAL assay with results in less than 1 hour
- Large range of sample types including: culture supernatants, plasma, serum, tissue, and urine

### Detect NGAL Across Diverse Model Systems



Calibration curve for NGAL (monkey) ELISA Kit (BPD-KIT-045).

## Versatile NGAL ELISA Kits for Detection of NGAL in a Variety of Sample Types



| SELECT ELISA KITS FOR NEPHROLOGY RESEARCH |             |            |                                   |         |                |            |
|---|-------------|------------|-----------------------------------|---------|----------------|------------|
| Product Name                              | Product #   | Size       | Sensitivity (Range)               | Species | Sample Types   | Assay Time |
| KIM-1 (human) ELISA Kit                   | ADI-900-226 | 1x96 wells | 1.279 pg/mL (7.813-500 pg/mL)     | H       | U              | < 2 Hours  |
| NGAL (human) ELISA Kit                    | BPD-KIT-036 | 1x96 wells | 4 mL pg/mL (10-1,000 mL pg/mL)    | H       | CS, P, S, T, U | < 4 Hours  |
| NGAL (dog) ELISA Kit                      | BPD-KIT-043 | 1x96 wells | 0.56 pg/mL (4-400 pg/mL)          | C       | CS, P, S, T, U | 4 Hours    |
| NGAL (monkey) ELISA Kit                   | BPD-KIT-045 | 1x96 wells | 1.5 pg/mL (range 10-200 pg/mL)    | MK      | CS, P, S, T, U | < 4 Hours  |
| NGAL (mouse) ELISA Kit                    | BPD-KIT-042 | 1x96 wells | 0.75 pg/mL (range 10-1,000 pg/mL) | M       | CS, P, S, T, U | 4 Hours    |
| NGAL (pig) ELISA Kit                      | BPD-KIT-044 | 1x96 wells | 1 pg/mL (range 10-400 pg/mL)      | P       | CS, P, S, T, U | 4 Hours    |
| NGAL (rat) ELISA Kit                      | BPD-KIT-046 | 1x96 wells | 0.5 pg/mL (range 4-400 pg/mL)     | R       | CS, P, S, T, U | 4 Hours    |

### LEGEND

Species: B= bovine, C = canine (dog), H = human, M = mouse, MK = monkey, P = porcine (pig), R= rat, RB = rabbit

Sample Types: AM = amniotic fluid, CL = cell lysates, CSF = cerebral spinal fluid, CS = culture supernatants, F = feces, M = milk, NL = nasal lavage, P = plasma, S = serum, SA = saliva, SS= sputum supernatant, T = tissue, TC = tissue culture, U = urine, WhB = white blood cells. Product tables includes both certified and cited species and sample types.

# NEUROSCIENCE RESEARCH


Neuroscience encompasses the study of the nervous system and includes neurodegeneration, which involves the progressive loss of structure or function of neurons. It is estimated the aging worldwide population will result in over 100 million sufferers of dementia by the year 2050, a postulate that continues to drive major research efforts in neurodegenerative diseases and the associated loss of cognitive function. Diseases such as Alzheimer's (AD), Parkinson's, Huntington's, and Amyotrophic Lateral Sclerosis are considered diseases of protein homeostasis (proteostasis), characterized by loss of specific neuronal populations, and the presence of inclusion bodies consisting of insoluble, unfolded proteins. Drug development programs for treatment of these diseases include modulators of key proteins and enzymes that regulate proper protein folding, modification, and clearance, seeking to reverse or prevent the accumulation of protein aggregates and toxic intermediates. Enzo offers several neuroscience ELISAs which include unique targets such as APP  $\Delta$ C31, LVV Hemorphin 7, and SMN to understand specific neurodegenerative diseases.

## APP $\Delta$ C31 ELISA KIT (ADI-900-227)

The progressive neurodegenerative disease, Alzheimer's disease, is characterized by senile plaques, neurofibrillary tangles and loss of synapses and neurons. AD has been largely viewed as a disease of toxicity being mediated by the accumulation of the amyloid beta ( $A\beta$ ) peptide as plaques within the brain resulting in damage to brain cells from the binding of damaging metals, reactive oxygen species production and direct damage to cellular membranes. Recent research has suggested that the  $A\beta$  peptide is a multifunctional peptide with non-pathological effects. Its association with AD, in conjunction with other proteins such as the amyloid precursor protein (APP), results in the imbalance between the processes of memory formation and normal forgetting. Enzo Life Sciences offers a novel ELISA kit to measure APP  $\Delta$ C31, an important amyloid precursor protein fragment with a unique pro-apoptotic mechanism leading to AD.

- Sensitive detection of an important APP fragment with a unique pro-apoptotic mechanism
- Validated for human cell lysates and cerebral spinal fluid sample types
- High specificity with low cross-reactivity to similar APP isoforms

### Reliable Screening of APP $\Delta$ C31 Inhibition




Inhibition experiment by ELISA (ADI-900-227) and Western blot analysis indicates both methods are in agreement using treatment of increasing concentrations of caspase inhibitor which reduces the production of APP  $\Delta$ C31.

## SEROTONIN ELISA KIT (ADI-900-175)

Serotonin (5-hydroxytryptamine, 5-HT) is a monoamine found in the central nervous system, gastrointestinal tract, and blood with broad physiological functions in neurotransmission, gastric motility, hemostasis, and cardiovascular integrity.

- Highly sensitive measurement of serotonin, detecting as little as 0.293 ng/mL
- No acylation required for standards and samples
- Thoroughly validated in multiple complex sample matrices
- High-throughput format with results in just 3 hours for up to 39 samples in duplicate

### Highly Sensitive Measurement of Serotonin




The Serotonin ELISA Kit (ADI-900-175) is a colorimetric competitive immunoassay kit with results in 3 hours.

## 24(S)-HYDROXYCHOLESTEROL ELISA KIT (ADI-900-210)

24(S)-Hydroxycholesterol (24-OHC), an enzymatically-generated side chain-hydroxylated derivative of cholesterol, is a pivotal marker in the study of cerebral cholesterol homeostasis. Cholesterol is unable to cross the blood-brain barrier, however, Cyp46 enzyme converts cholesterol to the more soluble 24-OHC, and this hydroxylated form of cholesterol is able to cross the blood-brain barrier. This conversion allows for the reduction of cholesterol in the brain and the efflux of 24-OHC from the brain into cerebral spinal fluid and blood. The flux of 24-OHC has been seen in patients with a variety of neurodegenerative diseases. In the instance of Alzheimer's disease, changes in 24-OHC concentrations may be indicative of different pathogenetic mechanisms and/or the progression of the disease. As in the case of multiple sclerosis, concentrations of 24-OHC have been shown to decrease, likely due to the loss of neuronal cells responsible for the synthesis.

- Measure as little as 0.78 ng/mL of 24(S)-Hydroxycholesterol in just 2 hours
- Convenient, user-friendly alternative to mass spectrometry
- Low cross-reactivity with structurally related molecules

### Detect Normal and Diseased Levels of 24(S)-Hydroxycholesterol




Normal human cerebral spinal fluid samples were diluted 1:2 in assay buffer and analyzed in the assay (ADI-900-210) for 24(S)-Hydroxycholesterol levels.

## SMN ELISA KIT (ADI-900-209)

Survival Motor Neuron (SMN) is a ~38 kDa protein produced primarily by the SMN1 gene. Deletion or mutation of the SMN1 gene results in a reduced level of full-length SMN protein and manifests as a range of neuromuscular phenotypes in humans as the disease spinal muscular atrophy (SMA). The kit provides reproducible, fully quantitative results for the detection of mouse and human SMN in cell lysates, surpassing semi-quantitative Western blot analysis.

- Detect as little as 50 pg/mL of SMN
- Results in just 3 hours from up to 39 samples in duplicate
- ELISA kit developed in collaboration with the SMA Foundation

### First Ready-to-use ELISA for Quantifying Human or Mouse SMN



Parallelism analysis using the SMN ELISA Kit (ADI-900-209) indicates antigen binding characteristics are similar to native SMN in human and mouse cell lysate samples with no matrix interference at the dilutions tested.

## SELECT ELISA KITS FOR NEUROSCIENCE RESEARCH

| Product Name                       | Product #                    | Size                     | Sensitivity (Range)           | Species             | Sample Types               | Assay Time         |
|------------------------------------|------------------------------|--------------------------|-------------------------------|---------------------|----------------------------|--------------------|
| 24(S)-Hydroxycholesterol ELISA Kit | ADI-900-210                  | 1x96 wells               | 0.78 ng/mL (0.39-100 ng/mL)   | Species Independent | CSF, CS, T                 | 2 Hours            |
| APP ΔC31 ELISA Kit                 | ADI-900-227                  | 1x96 wells               | 0.92 pM (11.72-1,500 pM)      | H                   | CL, CSF                    | ~2 Hours           |
| Dopamine ELISA Kit                 | ENZ-KIT188                   | 1x96 wells               | 0.938 ng/mL (1.56-100 ng/mL)  | Species Independent | P, S, T, TC                | 2 Hours            |
| Histamine ELISA Kit                | ENZ-KIT140                   | 1x96 wells               | 0.03 ng/mL (0.098-25 ng/mL)   | H, M, R, C          | P, TC, U                   | 2 Hours            |
| LVV Hemorphin 7 ELISA Kit          | ADI-900-205                  | 1x96 wells               | 6.1 pg/mL (9.8-10,000 pg/mL)  | Species Independent | S, T                       | 3 Hours            |
| Melatonin ELISA Kit                | ENZ-KIT150                   | 1x96 wells               | 0.162 ng/mL                   | Species Independent | S, SA, P, Fruit            | 2 Hours            |
| Oxytocin ELISA Kit                 | ADI-900-153A<br>ADI-901-153A | 1x96 wells<br>5x96 wells | 15 pg/mL (15.6-1,000 pg/mL)   | Species Independent | CS, CSF, M, P, S, SA, T, U | Overnight + 1 Hour |
| Serotonin ELISA Kit                | ADI-900-175                  | 1x96 wells               | 0.293 ng/mL (0.49-500 ng/mL)  | Species Independent | P, S, U                    | 3 Hours            |
| SMN ELISA Kit                      | ADI-900-209                  | 1x96 wells               | 50 pg/mL (50-3,200 pg/mL)     | H, M                | CL                         | 3 Hours            |
| Substance P ELISA Kit              | ADI-900-018A<br>ADI-901-018A | 1x96 wells<br>5x96 wells | 5.3 pg/mL (9.76-10,000 pg/mL) | Species Independent | CS, P, S, SA, T, U         | 3 Hours            |

### LEGEND

Species: B= bovine, C = canine (dog), H = human, M = mouse, MK = monkey, P = porcine (pig), R= rat, RB = rabbit

Sample Types: AM = amniotic fluid, CL = cell lysates, CSF = cerebral spinal fluid, CS = culture supernatants, F = feces, M = milk, NL = nasal lavage, P = plasma, S = serum, SA = saliva, SS= sputum supernatant, T = tissue, TC = tissue culture, U = urine, WhB = white blood cells. Product tables includes both certified and cited species and sample types.

# OXIDATIVE STRESS RESEARCH


The term oxidative stress reflects an imbalance in free radical formation within a cell or organism, most commonly in the form of reactive oxygen or nitrogen species (ROS/RNS). ROS/RNS such as superoxide anions, hydroxyl radicals, hydrogen peroxide, nitric oxide, and peroxynitrite originate from a variety of sources including changes in aerobic metabolism, immune activation, UV radiation, heme accumulation, and hypoxia. Failure of the cell's defense mechanisms to compensate for accumulating insults such as mitochondrial dysfunction, DNA damage, misfolded proteins, and lipid peroxidation can trigger programmed cell death pathways, and has been linked to clinically relevant diseases including cancer, cardiovascular disease, asthma, ischemia, diabetes, and neurodegenerative disease. Choose one of our sensitive ELISAs to detect a variety of oxidative stress biomarkers.

## DNA DAMAGE (8-OHDG) ELISA KIT (ADI-EKS-350)

Exposure of cells to oxidative and environmental stresses frequently results in the breakdown or oxidation of genomic DNA. Assays to evaluate the integrity of genomic DNA, or to assess the presence of oxidized DNA are frequently used as a means of verifying the onset of apoptosis or DNA damage.

- Quantify levels of < 1ng/mL in less than 2.5 hours
- Tested in a variety of biofluids (urine, serum, saliva)
- Convenient, colorimetric 96-well plate formats

### Rapidly Detect < 1 ng/mL 8-OHDG



Standard curve experiment using the DNA Damage ELISA Kit (ADI-EKS-350).

### SELECT ELISA KITS FOR OXIDATIVE STRESS RESEARCH

| Product Name                                  | Product #    | Size       | Sensitivity (Range)   | Species             | Sample Types                                      | Assay Time  |
|---|--------------|------------|---|---------------------|---|---|
| ADMA Direct (mouse/rat) ELISA Kit             | ALX-850-327  | 1x96 wells | 0.05 $\mu$ mol/L  | M, R                | P, S  | 15 to 20 Hours + 1 Hour                             |
| ADMA Direct ELISA Kit                         | ALX-850-323  | 1x96 wells | 0.05 $\mu$ mol/L (expected value: 0.45 $\mu$ mol/L $\pm$ 0.19 $\mu$ mol/L)    | H                   | P, S  | 15 to 20 Hours + 1 Hour                             |
| Cytochrome C (human) ELISA Kit                | ADI-900-141  | 1x96 wells | 6.03 pg/mL (28.13-900 pg/mL)  | H                   | CL  | 3 Hours 15 Minutes                                  |
| DNA Damage ELISA Kit                          | ADI-EKS-350  | 1x96 wells | 0.59 ng/mL (1.875-60 ng/mL)   | Species Independent | CL, CS, P, S, SA, U, seminal fluids, DNA extracts | < 2.5 Hours   |
| HO-1 (human) ELISA Kit                        | ADI-EKS-800  | 1x96 wells | 0.78 ng/mL (0.78-25 ng/mL)  | H                   | CL, P, S, T                                       | < 2.5 Hours   |
| HO-1 (rat) ELISA Kit                          | ADI-EKS-810A | 1x96 wells | 0.036 ng/mL (0.195-12.5 ng/mL)  | R                   | CL, P, S, T, bronchoalveolar lavage fluid         | < 3 Hours   |
| IMMUNOSET® HO-1 (human) ELISA Development Set | ADI-960-800  | 5x96 wells | 49 pg/mL (0.195-12.5 ng/mL)   | H                   | CL, T   | Plate coating - Overnight + 1 Hour; Assay - 3 Hours |
| IMMUNOSET® HO-1 (mouse) ELISA Development Set | ADI-960-071  | 5x96 wells | 96 pg/mL (0.195-12.5 ng/mL)   | M                   | CL, P, S, T                                       | Plate coating - Overnight + 1 Hour; Assay - 3 Hours |
| IMMUNOSET® HO-1 (rat) ELISA Development Set   | ADI-960-810  | 5x96 wells | 39 pg/mL (0.195-12.5 ng/mL)   | R                   | CL, P, S, T, bronchoalveolar lavage fluid         | Plate coating - Overnight + 1 Hour; Assay - 3 Hours |
| Myeloperoxidase (human) ELISA Kit             | ADI-900-115  | 1x96 wells | Assay Buffer 13: 0.028 ng/mL; Assay Buffer 31: 0.019 ng/mL (0.195-12.5 ng/mL) | H                   | CS, Li-heP, NL, P, SS, U                          | 3 Hours   |
| IMMUNOSET® PDI ELISA Development Set          | ADI-960-072  | 5x96 wells | 3.93 ng/mL (7.8-250 ng/mL)  | H, M, R             | CL, CS, P, T                                      | Plate coating - Overnight + 1 Hour; Assay - 3 Hours |
| Protein Carbonyl ELISA Kit                    | ALX-850-312  | 1x96 wells | Not applicable  | N/A                 | P, T  | 4.5 Hours   |
| SDMA (human) ELISA Kit                        | ALX-850-331  | 1x96 wells | 0.05 $\mu$ mol/L (0.1-2 $\mu$ mol/L)  | H                   | P, S  | 15 to 20 Hours + 1 Hour                             |
| Superoxide Dismutase (Cu/Zn) ELISA Kit        | ALX-850-033  | 1x96 wells | 0.04 ng/mL (0.08-5 ng/mL)   | H                   | CS, S   | < 2 Hours   |



# PROTEOSTASIS/CHAPERONES RESEARCH

Protein homeostasis or 'proteostasis' is the process that regulates proteins within the cell in order to maintain the health of both the cellular proteome and the organism itself. Proteostasis involves a highly complex interconnection of pathways that influence the fate of a protein from synthesis to degradation. As individual components are affected, the others adjust accordingly to maintain normal function. Disruption of one or more of these proteostasis influencers can manifest in pathologies such as Alzheimer's disease, cancer, and diabetes.

To fuel your proteostasis research, Enzo Life Sciences has an unrivaled portfolio of ubiquitin/proteasome and heat shock protein research reagents through our acquisition of Biomol International, Stressgen Bioreagents and Assay Designs. Our development efforts within the areas of protein synthesis, modification, and degradation continue to produce cutting-edge tools for dissecting key biological processes at the forefront of academic discourse in the proteostasis field. This commitment is evidenced most recently by first-to-market immunoassays, cell-based assays, and compound libraries for autophagy, an emerging pathway of interest in cellular regulation of protein turnover and its relation to disease.

## CHAPERONES/HEAT SHOCK PROTEINS


Building proteins within the cell is a complex process involving post-translational modification and protein folding. Because of the critical role proteins play, incorrect folding can have devastating consequences. Enzo has a number of assay kits and reagents directed at monitoring the delicate balance of protein synthesis, folding, and degradation. These include highly conserved heat shock proteins (HSPs) that are constitutively expressed and function as molecular chaperones which facilitate the synthesis and folding of proteins. HSPs also participate in protein assembly, export, turnover and regulation. Under stressful conditions such as heat shock, pH shift or hypoxia, increased expression of HSPs protects the cell by stabilizing unfolded proteins, giving the cell time to repair or re-synthesize damaged proteins. Enzo offers a number of kits, antibodies, proteins, inhibitors and more reagents to study HSPs/chaperones, HSF's, GRP's, ERp's, clusterin and other stress-related molecules. Our ELISAs for proteostasis and chaperone research are sensitive, widely published and available for multiple species and sample types.

## HSP90 $\alpha$ (HUMAN) ELISA KIT (ADI-EKS-895)

The HSP90 family of heat shock proteins represents one of the most abundantly expressed and highly conserved families of cellular chaperones whose expression can be upregulated under conditions of cellular stress. This family includes cytoplasmic (HSP90 $\alpha$  and  $\beta$ ), ER (GRP94), and mitochondrial (TRAP1) localized members. Structurally, HSP90 is characterized by an N-terminal ATP-binding domain, a medial substrate binding domain, and a C-terminal dimerization motif. HSP90 dimers function in cooperation with co-chaperones (e.g. HSP40, HSP70, Hop, p23) to stabilize a multitude of client protein substrates, including steroid hormone receptors, protein kinases, and transcription factors. The essential binding and hydrolysis of ATP by HSP90 is inhibited by ansamycin drugs (e.g. geldanamycin, 17-AAG) which occupy the N-terminal HSP90 nucleotide binding pocket. Many HSP90 client proteins such as erbB2/Her-2, c-raf, BCR-Abl, p53, and hTERT, are members of well characterized oncogenic pathways, making HSP90 inhibitors useful anticancer agents.

- Detect as little as 50 pg/mL of human HSP90 $\alpha$  in just 3 hours
- Highly reproducible with less than 10% variation between assays
- No cross-reactivity with HSP90 $\beta$ , GRP94, HSP60, or HSP70 (HSP72)
- Obtain quantitative results for cell lysates, serum and tissue sample types

### Most Sensitive Detection of Human HSP90 $\alpha$



Standard curve experiment using the HSP90 $\alpha$  (human) ELISA Kit (ADI-EKS-895).

#### LEGEND

Species: B= bovine, C = canine (dog), H = human, M = mouse, MK = monkey, P = porcine (pig), R= rat, RB = rabbit

Sample Types: AM = amniotic fluid, CL = cell lysates, CSF = cerebral spinal fluid, CS = culture supernatants, F = feces, M = milk, NL = nasal lavage, P = plasma, S = serum, SA = saliva, SS= sputum supernatant, T = tissue, TC = tissue culture, U = urine, WhB = white blood cells. Product tables includes both certified and cited species and sample types.

# PROTEOSTASIS/CHAPERONES RESEARCH

## SELECT ELISA KITS FOR PROTEOSTASIS/CHAPERONES RESEARCH

| Product Name  | Product #    | Size       | Sensitivity (Range)             | Species    | Sample Types       | Assay Time   |
|---|--------------|------------|---------------------------------|------------|--------------------|--|
| IMMUNOSET® αB-Crystallin ELISA Development Set                  | ADI-960-074  | 5x96 wells | 0.59 ng/mL (1.25-40 ng/mL)      | H, M, R, B | CL                 | Plate coating - Overnight + 1 Hour; Assay - 3 Hours            |
| IMMUNOSET® Grp75 ELISA Development Set                          | ADI-960-143  | 5x96 wells | 0.762 ng/mL (3.125-100 ng/mL)   | H, M, R    | CL                 | Plate coating - Overnight + 1 Hour; Assay - 3 Hours            |
| Grp78/BiP ELISA Kit   | ADI-900-214  | 1x96 wells | 8.4 ng/mL (1.4-4,500 ng/mL)     | H, M, R    | CL, S              | 2.5 Hours  |
| IMMUNOSET® Grp94 ELISA Development Set                          | ADI-960-077  | 5x96 wells | 1.29 ng/mL (6.25-200 ng/mL)     | H, M, R, C | CL, T              | Plate coating - Overnight + 1 Hour; Assay - 3 Hours            |
| HSF1 ELISA kit  | ADI-900-198  | 1x96 wells | 61 pg/mL (0.39-12.5 ng/mL)      | H, M       | CL, T              | 3 Hours  |
| [pSer <sup>326</sup> ]HSF1 ELISA Kit                            | ADI-900-199  | 1x96 wells | 35 pg/mL (0.39-12.5 ng/mL)      | H, M       | CL, T              | 3 Hours  |
| IMMUNOSET® HSP25 (rodent) ELISA Development Set                 | ADI-960-075  | 5x96 wells | 0.38 ng/mL (0.78-25 ng/mL)      | M, R       | CL                 | Plate coating - Overnight + 1 Hour; Assay - 3 Hours            |
| HSP27 (human) ELISA kit   | ADI-EKS-500  | 1x96 wells | 0.39 ng/mL (0.39-25 ng/mL)      | H          | CL, P, S, T        | < 3 Hours  |
| [pSer <sup>15</sup> ]HSP27 (human) ELISA Kit                    | ADI-900-170  | 1x96 wells | 10.15 pg/mL (31.25-1,000 pg/mL) | H          | CL, P, S           | 3 Hours  |
| [pSer <sup>78</sup> ]HSP27 (human) ELISA Kit                    | ADI-900-165  | 1x96 wells | 4.30 pg/mL (31.25-1,000 pg/mL)  | H          | CL, P, S           | 3 Hours  |
| IMMUNOSET® HSP27 High Sensitivity (human) ELISA Development Set | ADI-960-076  | 5x96 wells | 97 pg/mL (0.1-3.2 ng/mL)        | H          | CL, P, S           | Plate coating - Overnight + 1 Hour; Assay - 2 Hours 45 Minutes |
| HSP60 (human) ELISA Kit   | ADI-EKS-600  | 1x96 wells | 3.125 ng/mL (3.125-100 ng/mL)   | H          | CL, S, T           | < 3 Hours  |
| Anti-HSP60 IgG/A/M (human) ELISA Kit                            | ADI-EKS-650  | 1x96 wells | 2.88 ng/mL (7.81-250 ng/mL)     | H          | S                  | < 4 Hours  |
| HSP70 ELISA Kit   | ADI-EKS-700B | 1x96 wells | 200 pg/mL (780-50,000 pg/mL)    | H, M, R, F | CL, S, T           | 4.5 Hours  |
| HSP70 High Sensitivity ELISA Kit                                | ADI-EKS-715  | 1x96 wells | 90 pg/mL (0.20-12.5 ng/mL)      | H, M, R    | CL, S, P, U        | 4.5 Hours  |
| AMP'D® HSP70 High Sensitivity ELISA Kit                         | ENZ-KIT-101  | 1x96 wells | 7 pg/mL (0.039-5 ng/mL)         | H, M, R    | S, P               | 4.5 Hours  |
| Anti-HSP70 IgG/A/M (human) ELISA Kit                            | ADI-EKS-750  | 1x96 wells | 6.79 ng/mL (31.25-1,000 ng/mL)  | H          | S                  | < 4 Hours  |
| HSP70B' ELISA Kit   | ADI-EKS-725A | 1x96 wells | 62 pg/mL (0.156-10 ng/mL)       | H          | CL, S, T           | < 4.5 Hours  |
| HSP90α (human) ELISA Kit  | ADI-EKS-895  | 1x96 wells | 50 pg/mL (62.5-4,000 pg/mL)     | H          | CL, S, T           | < 3 Hours  |
| Proteasome ELISA Kit  | BML-PW0575   | 1x96 wells | Not applicable                  | N/A        | Biological samples | < 3.5 Hours  |

### LEGEND

Species: B= bovine, C = canine (dog), H = human, M = mouse, MK = monkey, P = porcine (pig), R= rat, RB = rabbit

Sample Types: AM = amniotic fluid, CL = cell lysates, CSF = cerebral spinal fluid, CS = culture supernatants, F = feces, M = milk, NL = nasal lavage, P = plasma, S = serum, SA = saliva, SS= sputum supernatant, T = tissue, TC = tissue culture, U = urine, WhB = white blood cells. Product tables includes both certified and cited species and sample types.

# %CV IN ELISA: HOW TO REDUCE THEM AND WHY THEY ARE IMPORTANT

## INTRODUCTION

For ELISA users, having a low coefficient of variability (CV or %CV) between sample replicates is crucial in demonstrating an assay was well-run and the resultant data is precise and accurate. Reliable assay results are assessed by standardized measures such as coefficient of variability.

The coefficient of variability is a dimensionless numerical ratio used to describe the level of variability within a population independently of the absolute values of the observations. In statistical analysis of numerical data, if your absolute values are similar, sample populations can be assessed by using standard deviations; when absolute values vary, you must consider using an approach such as %CV, to assess the precision of a laboratory technique. CV is calculated by dividing the standard deviation ( $\sigma$ ) of a set of measurements by the mean ( $\mu$ ) of the set which is then expressed as a percentage of variation to the mean.

$$C_v = \frac{\sigma}{\mu}$$

General mathematical formula for Coefficient of Variability (%CV).

In ELISA data interpretation, %CV can highlight inconsistencies among sample replicates which is demonstrated in the data as variation among Optical Density (OD) readouts post-assay. These directly reflect the performance of the assay in the hands of the end-user. There are two types of %CVs that are used to express the precision of immunoassay results: intra-assay CV and inter-assay CV. Intra-assay CV is a measure of the variance between data points within an assay, meaning sample replicates ran within the same plate. Inter-assay CV is a measure of the variance between runs of sample replicates on different plates that can be used to assess plate-to-plate consistency. As a general guideline, to gauge the overall reliability of your immunoassay results, inter-assay %CV should be less than 15% while intra-assay %CV should be less than 10%.

It is important to identify the causes of high %CV in ELISAs. Human technical error can play a role such as inaccurate pipetting technique, splashing of reagents between wells, drying out of the wells, inconsistent sample handling e.g. variability due to freeze-thaw cycles among samples, and use of differing filter settings to start. Additionally, high %CV can be the result of machine error such as usage of uncalibrated automated machine pipettes, uncalibrated plate readers, and inappropriate plate reader software settings to analyze samples in wells. Lastly, plate, sample, and reagent contamination can lead to a high %CV. Cross-contamination between reagents can occur from handling errors leading to bacterial or fungal contamination of samples and reagents derived from compromised sterility.

## How To Reduce %CV

1. Use fresh pipette tips for each addition. Discarding old tips for new tips between each addition can prevent cross-contamination between wells which can keep the background and %CV low between sample replicates.
2. It is not recommended to pour excess reagent from a reservoir back into the original bottle as this can unintentionally contaminate your stock solution which will be reflected in your readouts and could produce a high %CV.
3. Pre-wetting pipette tips 2-3 times in the solution that is to be pipetted also can help to improve %CV.
4. Proper pipette technique can also be a major factor in lowering %CV. Some general practices include: holding the pipette vertically (and not at an angle), aspirating slowly and smoothly, making sure the tip touches the vessel when withdrawn to avoid extra liquid outside the tip, and ensuring the tip is just under the surface of the liquid in the reservoir when dispensing.
5. Regular re-calibration of autopipettes is a good practice that can extend to other instruments including the plate washer and reader.

| Group: |        |        |        |        |             |           |             |
|--------|--------|--------|--------|--------|-------------|-----------|-------------|
| Sample | Values | Net OD | %Bound | Result | Mean Result | Std. Dev. | CV%         |
| 100    | 0.541  | 0.543  | 58.68  | 18.449 | 17.172      | 1.806     | <b>10.5</b> |
|        | 0.565  | 0.567  | 61.276 | 15.894 |             |           |             |
| 102    | 0.531  | 0.533  | 57.599 | 19.626 | 21.196      | 2.22      | <b>10.5</b> |
|        | 0.507  | 0.509  | 55.003 | 22.766 |             |           |             |
| 103    | 0.399  | 0.401  | 43.321 | 45.333 | 41.39       | 5.576     | <b>13.5</b> |
|        | 0.428  | 0.43   | 46.458 | 37.447 |             |           |             |
| 141    | 0.451  | 0.453  | 48.945 | 32.312 | 31.709      | 0.852     | <b>2.7</b>  |
|        | 0.457  | 0.459  | 49.594 | 31.106 |             |           |             |
| 150    | 0.471  | 51.109 | 28.482 | 26.733 | 2.473       |           | <b>9.3</b>  |
|        | 0.492  | 0.494  | 53.38  | 24.985 |             |           |             |
| 151    | 0.45   | 0.452  | 48.837 | 32.518 | 34.155      | 2.315     | <b>6.8</b>  |
|        | 0.435  | 0.437  | 47.215 | 35.792 |             |           |             |
| 152    | 0.357  | 0.359  | 38.778 | 60.662 | 57.035      | 5.129     | <b>9</b>    |
|        | 0.375  | 0.377  | 40.725 | 53.408 |             |           |             |
| 157    | 0.362  | 0.364  | 39.319 | 58.529 | 60.036      | 2.132     | <b>3.6</b>  |
|        | 0.355  | 0.357  | 38.561 | 61.543 |             |           |             |
| 159    | 0.425  | 0.427  | 46.133 | 38.184 | 35.664      | 3.563     | <b>10</b>   |
|        | 0.447  | 0.449  | 48.513 | 33.144 |             |           |             |

An example of high and low intra-assay %CV values from data using Aldosterone ELISA Kit. Bolded in red are high% CVs (over 10) seen amongst sample replicates. Bolded in green are low %CV (underneath or at 10) seen between sample replicates

1. Always read and follow the product instructions. The information contained in the kit manual is your first source of knowledge for your particular assay. If there is any special handling information that you need to be aware of, it will be recorded in the kit insert. The insert has been written to provide you with the maximum amount of useful information in one convenient location. NOTE: Do not use the same kit insert repeatedly since the kits and inserts are periodically improved. The information that comes packaged with each kit contains information about any improvements and is the most current information available.
2. Check the calibration of the pipets that will be used for the assay. Each manufacturer should supply the recommended calibration protocol and performance specifications with each instrument. You need to verify that the volumes you are measuring are correct.
3. Do not mix components from different lots or kits. Each component is lot-specific and designed to give you a defined level of performance when used before the expiration date.
4. Always store the kit components at the recommended conditions. In some instances, special components may need to be frozen for maximum stability. Each component will have specific storage information printed on its label. Kits containing materials that require special storage should be clearly marked on the outside of the box.
5. Do not use the kit past the expiration date marked on the box. While many of the individual components may have expiration dates that extend beyond the kit date, the total kit expiration date has been determined based on the stability of all the materials stored at the suggested conditions.
6. Allow the plate to warm to room temperature before opening. The wells in the microtiter plate are coated with an antibody solution. While this coating is stable and robust, moisture will corrupt it causing a decrease in performance. The coated plate is shipped in a reusable mylar bag containing a mini-desiccant packet to ensure dry conditions. When not in use, the wells should be stored at 4°C, tightly sealed in the mylar bag with desiccant. Do not throw the plate frame away after your first assay. You will need it for the remaining wells stored in the foil bag.
7. Allow all reagents to warm to room temperature before opening unless instructed otherwise. Many of the reagents contain temperature-dependent components that may come out of solution when cold. Using the reagents at room temperature ensures that you have a consistent formulation from assay to assay.
8. Pre-rinse pipet tips with the reagent before transferring a volume. To pre-rinse, insert the end of the pipet tip just under the surface of the reagent and draw up a volume according to the pipet manufacturer's directions using a slow even motion. Expel the volume using the same steady motion. Typically you should pre-rinse the tip three times, using the fourth "draw" as your transfer volume. This is good lab practice regardless of the tip manufacturer's claims of low retention. Precision and accuracy will be maximized by taking the time to pre-rinse tips on a consistent basis. You need to pre-rinse a tip only once if you are using it to make multiple transfers. Always change the pipet tip if it has been accidentally contaminated.
9. Always use a fresh pipet tip for each standard, sample or reagent. Do not use the same tip for your standards even when pipetting from the lowest to highest concentrations as you will sacrifice precision.
10. Pipet the first reagent into the bottom of each well. Pipet subsequent volumes in different locations on the sides of each well, being careful to avoid cross-contaminating reagents. A simple way to avoid contamination is to pipet one reagent into all of the wells at the same location (e.g.: left side wall), then turn the plate 90° and pipet the next reagent at the same position (e.g.: left side wall). Because you rotated the plate, this position will put your pipet tip at a new location within the well.
11. Be careful handling the microtiter wells and reagents to prevent contamination. You have a variety of enzymes on your skin including proteases and endogenous alkaline phosphatase. If you accidentally touch a pipet tip, you can transfer these enzymes to your reagents resulting in unusable kit components. This is especially important for your substrate. The use of contaminated substrate will result in a less sensitive assay.
12. If culture media (CM) samples are being analyzed, then the standards should be diluted in non-conditioned media. Do not assume that all CM are the same. Not only does the medium differ in formulation, but serum supplements will also differ by type and lot. You should expect to find a basal level of endogenous analyte in media containing serum supplements. In addition, other medium supplements can contribute to the overall detection in an assay. If non-conditioned media is used as the standard diluent, then all samples will be relative to the standard curve. Using non-conditioned media usually results in a modest depression of the standard curve relative to standards diluted in kit assay buffer. This depression is acceptable since it is more important for the samples to be relative to the appropriate standard curve.
13. It is simply good lab practice to run duplicates for repeatability verification. In addition, you maximize your results when you run all wells in duplicate. It is easier to edit one outlying data point than to run the experiment over again because you were not able to get an appropriate measurement for some standard or sample. An outlying measurement can result from an improper amount of a reagent being added to a well or losing a portion of the volume during preparation or incubation. While one never thinks that this will happen to them, unforeseen things sometimes happen with unfortunate results.

### Weak Color Development

**Was substrate added at the correct point in the assay?**

See the assay procedure provided in the instruction manual.

**Were the antibody and conjugate added at the correct time?**

See the assay procedure provided in the instruction manual. The antibody and conjugate provided are often color-coded for convenience and to help reduce laboratory errors.

**Did all of the components belong to the specific kit being used?**

When customers order more than one type of kit, they can sometimes confuse the reagents between kits.

**How long was the substrate incubation?**

It is possible that Stop Solution was added to the plate without allowing the full substrate incubation.

**What were the conditions of the substrate incubation?**

If a plate is left to incubate on a cold lab bench or under a drafty area during ambient incubations, signal values (e.g. optical density) may be lower than expected.

**Were reagents brought to room temperature prior to use?**

It is important to ensure that all reagents are brought to room temperature prior to use, or as mentioned in the product specific instruction manual. Usually leaving the kit out on the bench top at ambient temperature for about half an hour prior to setting up the assay will be sufficient, when the reagents can be stored at 4°C. Frozen volumes take a little more time to come to room temperature. Do not thaw frozen reagents in a water bath. If a different standard/sample diluent is used (such as culture media) this must also be warmed.

**Was the plate read at the correct wavelength?**

See the assay procedure provided in the instruction manual to ensure you're reading the plate at the correct wavelength. It may be necessary to check the filters in your plate reader and the program using during reading. If others used the instrument, they may make changes to the settings for their experiment.

**Were the proper volumes of reagents added?**

See the assay procedure provided in the instruction manual.

**What were the conditions of the incubations?**

If the incubation times and temperatures are not observed, this can lead to lower than expected signal values (e.g. optical density). Pay attention that in air-conditioned rooms the temperature does not drop below 21°C.

**How was the plate shaken during incubations (if required)?**

If customers do not have a plate shaker, they will often use an orbital flask shaker or some other piece of equipment. This is not a problem as long as the liquid is vigorously displaced about 3/4 of the way up the sides of the wells without coming out. It is very important that the plate is secured into place. If the plate is not shaken and it is required in the procedure, a longer incubation may be necessary to bring the reagents to equilibrium.

**How long after the addition of Stop Solution was the plate read?**

The plate needs to be read at the correct wavelength as soon as possible after the addition of the Stop Solution. We generally recommend that the plate be read within 10 minutes.

### Poor Standard Curve

#### **What was used as the standard diluent?**

Diluents other than the supplied assay buffer may contain interfering substances that can affect the standard curve.

#### **How was the precision of the standard curve?**

If the %CV values for the standard curve signal values (e.g. optical density) are consistently above 5%, it may be a good idea to pay particular attention to pipetting technique. If the standard curve signal values were acceptable but the sample precision was not, the problem relates to the sample. Also, see recommendations under "Poor Precision".

#### **Were the Blank and NSB values subtracted out?**

If the net signal values (e.g. optical density) are not used, the signal values will appear higher than those presented in the sample data in the instruction manual.

#### **How were the standard dilutions prepared?**

It is important that test tubes of an appropriate size and material are used. Standard dilutions must be properly mixed (e.g. vortexed) while preparing the serial dilutions. It is also crucial that the standard dilutions be prepared and used within the time specified in the product specific instruction manual. Never store unused standard dilutions for a later use.

### Poor Precision

#### **Were the wells washed properly?**

All wells receive the same treatment during the wash step. If some are washed less than others, this can translate to poor precision. It is important that the plate is washed thoroughly. If plate washing is troublesome, a squirt bottle can be filled with diluted Wash Buffer and all of the wells completely filled from this. The plate contents can be dumped into the sink and shaken to remove excess buffer. This should be repeated for the number of times recommended in the instruction manual. It is important to remember that adding too little Wash Buffer can result in high background, while adding too much is not a problem. The contents of the wells should be aspirated and the plate tapped dry on lint-free paper towels. Plate tapping should consist of a few taps since excessive tapping may lead to plate drying and inconsistent assay.

#### **Were the wells aspirated sufficiently after the wash steps?**

It is very important that as little Wash Buffer as possible remains in the wells after aspiration. Residual buffer can cause dilution of subsequent reagents. After the last wash step, it is a good idea to hit the plate several times over a piece of paper towel to remove excess buffer.

#### **How were reagents pipetted into wells?**

In order to eliminate precision error, customers need to remember to pre-rinse all pipet tips used in the assay. We usually recommend that the customer draw up the liquid into the tip and aspirate it three times prior to addition into the well. Regular pipet calibration and maintenance is also essential to ensure that the tips fit properly and that the correct volumes are dispensed. Be sure reagents are not splashed between wells or outside of the wells during pipetting (especially if using repeater pipets).

### High Background

#### **How was the plate washed?**

It is important that the plate is washed thoroughly. If plate washing is troublesome, a squirt bottle can be filled with diluted Wash Buffer and all of the wells completely filled from this. The plate contents can be dumped into the sink and shaken to remove excess buffer. This should be repeated for the number of times recommended in the instruction manual. It is important to remember that adding too little Wash Buffer can result in high background, while adding too much is not a problem. The contents of the wells should be aspirated and the plate tapped dry on lint-free paper towels.

#### **What were the incubation times and temperatures?**

If the plate was incubated for too long or at a higher than recommended temperature, high background could result.



### Edge Effects

**Where was the plate incubated?**

Oftentimes the conditions for ambient incubations can be less than ideal. If there is a draft in the area or the plate is incubated on a cold lab bench, this can lead to uneven color development.

**If multiple plates were run, were they stacked on top of each other during incubation?**

Multiple plates should only be incubated in a single layer. This will assure that no area of the plate is at a different temperature than any other.

**If a non-ambient incubation was required, was the plate properly sealed?**

Making sure that the plate sealer is tightly covering all of the wells will help to discourage uneven evaporation of the well contents, or condensation for colder incubation conditions.

### Drift

**Were reagents brought to room temperature prior to use?**

If the reagents are not at a constant temperature prior to their addition into the wells, the results from one side of the plate to the other can differ depending on the temperatures at addition.

**Was the set-up of the assay interrupted?**

If the assay is interrupted at any point during the addition of reagents, it is possible that differing results will be seen before the interruption versus after. The wells that had reagents added before the interruption will have been incubating for longer than those after.

Have specific troubleshooting or technical questions?  
Call our tech support team! We are here to help!

|        |                 |  |
|--------|-----------------|--|
| U.S.   | 800.942.0430    | <a href="mailto:info-usa@enzolifesciences.com">info-usa@enzolifesciences.com</a> |
| Europe | +41 61 926 8989 | <a href="mailto:info-eu@enzolifesciences.com">info-eu@enzolifesciences.com</a>   |



**Global Headquarters**  
**ENZO LIFE SCIENCES, INC.**  
10 Executive Blvd.  
Farmingdale, NY 11735  
Ph: 800.942.0430  
Fax: 631.694.7501  
[info-usa@enzolifesciences.com](mailto:info-usa@enzolifesciences.com)

**European Sales Office**  
**ENZO LIFE SCIENCES (ELS) AG**  
Industriestrasse 17  
CH-4415 Lausen, Switzerland  
Ph: +41 61 926 8989  
Fax: +41 61 926 8979  
[info-eu@enzolifesciences.com](mailto:info-eu@enzolifesciences.com)

#### LOCAL EUROPEAN OFFICES

**Belgium, The Netherlands & Luxembourg**

Enzo Life Sciences BVBA  
Avenue Louise 65/Box 11  
1050 Bruxelles  
Belgium  
Ph: +32 3 466 0420  
Fax: +32 3 808 7033  
[info-be@enzolifesciences.com](mailto:info-be@enzolifesciences.com)

**France**

Enzo Life Sciences (ELS) AG  
Branch Office Lyon  
13, avenue Albert Einstein,  
F-69100 Villeurbanne, France  
Ph: +33 472 440 655  
Fax: +33 481 680 254  
[info-fr@enzolifesciences.com](mailto:info-fr@enzolifesciences.com)

**Germany**

Enzo Life Sciences GmbH  
Basler Strasse 57a  
DE-79540 Lörrach  
Germany  
Ph: +49 7621 5500 526  
Fax: +49 7621 5500 527  
[info-de@enzolifesciences.com](mailto:info-de@enzolifesciences.com)

**UK & Ireland**

Enzo Life Sciences (UK) Ltd.  
1 Colleton Crescent  
Exeter EX2 4DG  
Ph: 0845 601 1488 (UK customers)  
Ph: +44 1392 825900  
Fax: +44 1392 825910  
[info-uk@enzolifesciences.com](mailto:info-uk@enzolifesciences.com)

For local distributors and detailed product information visit us online:  
**[www.enzolifesciences.com](http://www.enzolifesciences.com)**