



EuroClone®

serving science through innovation

Composition of HBSS

Code: ECB4006

Inorganic salt :	
Calcium chloride Dihydrate	185 mg/l
Magnesium chloride Hexahydrate	100 mg/l
Magnesium sulfate anhydre	97,67 mg/l
Potassium chloride	400 mg/l
Potassium phosphate monobasic anhydre	60 mg/l
Sodium bicarbonate	350 mg/l
Sodium chloride	8000 mg/l
Sodium phosphate dibasic anhydrous	47.88 mg/l
Other components :	
D-Glucose	1000 mg/l
Phenol red	11 mg/l

EuroClone S.p.A.
Via Figino 20/22, 20016 Pero (MI) - Italy
Phone: +39.02.38195.1 - Fax +39.02.38101465
e-mail: info@euroclone.it



IOAIR®



www.euroclone.it