



EuroClone®

serving science through innovation

Composition of PBS Dulbecco's Phosphate Buffered Saline

Code: ECB4053

Inorganic salt :

Calcium chloride Dihydrate	133 mg/l
Magnesium chloride Hexahydrate	100 mg/l
Potassium chloride	200 mg/l
Potassium phosphate monobasic anhydrous	200 mg/l
Sodium chloride	8000 mg/l
Sodium phosphate dibasic anhydrous	1150 mg/l

ECB4053_Ed1_0218_691

