

PROMO FEDELTA'

EuroClone®
serving science through innovation



Piastre e strip per PCR e Real Time PCR: EuroClone premia la tua FEDELTA' con una doppia opportunità*

Acquista una selezione dei prodotti in offerta e approfitta dei vantaggi a te riservati.

FRAMESTAR® PLATES

Codice	Descrizione	Formato	Prezzo
ECPCR0384C	Primo® FrameStar® 384 well, PCR plates clear wells	50 plates	
ECPCR0960C	Primo® FrameStar® 96 PCR plates clear frame	50 plates	
ECPCR0770C	Primo® FrameStar® 96well PCR Plate, semi-skirted	50 plates	
ECPCR0730C	Primo® FrameStar® 96 PCRplates semi skirted with upstand clear	50 plates	
ECPCR0910C	Primo® 96 PCR plates semi skirted recommended for ABI Fast cyclers	50 plates	

STANDARD PLATES

Codice	Descrizione	Formato	Prezzo
ECPCR0750C	Primo® 96 well plate non-skirted, clear wells	50 plates	
ECPCR0760	Primo® 96 well semi skirted PCR plate, natural	50 plates	
ECPCR0740C	Primo® standard 96 PCR plates skirted low profile	50 plates	
ECPCR0955	Primo® 96 well semi skirted plate for Roche LC480 low profile, white well	50 plates	

PCR STRIPS & CAPS

Codice	Descrizione	Formato	Prezzo
ECPCR0208D	Primo® 0.2ml, 8 tubes/strip. Clear wells + domed caps	125 strips	
ECPCR0208F	Primo® 0.2ml, 8 tubes/strip. Clear wells + flat optical caps	125 strips	
ECPCR0754C	Primo® Low profile PCR strip of 8 tubes, clear+ optical caps	120 strips	

SEALING TAPES & STRIPS

Codice	Descrizione	Formato	Prezzo
ECPCR0500	Primo® Transparent Adhesive PCR Seal	100 films	
ECPCR0550	Primo® Aluminium Foil Adhesive Seal	100 films	
ECPCR0560	Primo® qPCR Adhesive Seal	100 films	
ECPCR0751	Primo® strips of 8 flat optical caps for PCR plates	300 strips	
ECPCR0752	Primo® strips of 8 domed caps for PCR plates	300 strips	

*

1 ORDINA almeno € 150,00 e NON paghi le spese di gestione ordine e di ghiaccio secco!

2 ORDINA almeno € 800,00 e richiedi, in sconto merce, una tra le due pipette Primo® qui indicate.

**Pipetta multicanale
8 canali 0.5 – 1 µl
Codice EPC80010**

**Pipetta multicanale
8 canali 5 – 50 µl
Codice EPC80050**



Riporta sull'ordine il codice della multicanale che hai scelto in sconto merce.

- Offerta valida fino al **31/12/2017**
- L'indicazione del codice **PROMOPLATES** sull'ordine è indispensabile per avere diritto allo sconto sulle spese e per ricevere lo strumento scelto in sconto merce!
- Non cumulabile con altri sconti e/o offerte in essere.

continua...

... qui sotto due tabelle utili per scegliere "NOSTRE" plastiche **Primo®**

TABELLA COMPATIBILITA' CON LA STRUMENTAZIONE qPCR PIU' DIFFUSA	Abi Life Technologies					Bio-Rad			BioRad MJ Research	Eppendorf	Roche	Stratagene/Agilent						
	7000, 7300, 7500, 7700, 7900	7900 HT	FAST 7500, 7900, 7900HT	StepOne Plus™	QuantStudio™	ViiA7™	CFX 384	CFX96	iCycler™ IQ, iQ™ 4, iQ™ 5	MiyiQ™	Chromo4™	Opticon™, Opticon2™	Mastercycler™ ep realplex	LC96	Nano	Mx3000P™	Mx3005P™	Mx4000™
FrameStar® Plates																		
ECPCR0384C		X			X	X	X											
ECPCR0730C	X	X			X	X						X						X
ECPCR0770C	X				X	X		X	X			X				X	X	X
ECPCR0910C			X	X	X	X												
ECPCR0960C								X		X	X	X						
Standard Plates																		
ECPCR0740C								X		X	X	X					X	
ECPCR0750C								X	X	X		X				X	X	X
ECPCR0760	X				X	X						X				X	X	X
ECPCR0955													X					
Strips																		
ECPCR0208F	X							X	X			X				X	X	X
ECPCR0754C/ECPCR0754														X				

TABELLA DI COMPATIBILITA' PER L'UTILIZZO DELLE PIASTRE SUI TERMOCICLATORI PIU' COMUNI																																							
	Abi Life Technologies						Bioer				BioRad				BioRad MJ Research	Eppendorf	EuroClone Bioer	EuroClone PeqLab	Stratagene	Thermo Fisher Scientific																			
	Veriti 384 well Block	Veriti 0.1 ml 96 well Block	Veriti 0.2 ml 96 well Block	GeneAmp® 2400	GeneAmp® 2700/2720/9600	GeneAmp®9700	GeneAmp®9800 FAST Block	XP Cycler	Gene Pro	Little Genius	Life Express	GeneQ	LifePro	C1000	Genecycler	iCycler™	MyCycler™	S1000	T100	Mini Gradient	PTC100™ (96-well block only)	DNA Engine™, DNA Dyad™, DNA Tetrad™	M384	MasterCycler®	MasterCycler® EP Gradient	MasterCycler® Gradient	EuroCycler Flex	EuroCycler Mini	EuroCycler 96 Gradient	EuroClone PeqStar 384	EuroClone PeqStar 96 Universal/Gradient	Gradient Cycler	MultiBlock System & MBS	PCR Sprint					
FrameStar® Plates																																							
ECPCR0384C	X					X		X	X				X					X			X	X					X			X					X		X		
ECPCR0960C								X				X	X							X	X					X	X		X				X		X	X	X		
ECPCR0770C			X		X	X		X					X		X	X	X	X			X	X			X	X					X								
ECPCR0730C			X		X	X							X		X	X	X	X			X	X			X	X													
ECPCR0910C		X					X																																
Standard Plates																																							
ECPCR0750C					X	X				X	X	X	X	X	X	X	X	X	X	X	X	X			X	X		X						X	X	X	X		
ECPCR0760			X		X	X		X					X		X	X	X	X			X	X			X	X					X								
ECPCR0740C								X				X	X					X			X	X			X	X		X	X	X									
Tubes and Strips																																							
ECPCR0208D				X	X	X		X	X	X	X	X	X	X	X	X	X	X	X	X	X	X			X	X	X	X	X	X	X			X		X	X	X	
ECPCR0208F				X	X	X		X	X	X	X	X	X	X	X	X	X	X	X	X	X	X			X	X	X	X	X	X			X		X		X	X	
ECPCR0754C		X						X	X									X			X	X			X	X	X	X	X			X			X	X	X	X	
ECPCR0754		X						X	X									X			X	X			X	X	X	X	X			X			X	X	X	X	

