



EuroClone®

serving science through innovation

Composition of HBSS w/o Ca⁺⁺ w/o Mg⁺⁺

Code: ECB4007L

Inorganic salt :

Potassium chloride	400 mg/l
Potassium phosphate monobasic anhydrous	60 mg/l
Sodium bicarbonate	350 mg/l
Sodium chloride	8000 mg/l
Sodium phosphate dibasic anhydrous	47.88 mg/l

Other components:

D-Glucose	1000 mg/l
Phenol red	11 mg/l

Euroclone S.p.A.

Via Figino 20/22, 20016 Pero (MI) - Italy
Phone: +39.02.38195.1 - Fax +39.02.38101465
e-mail: info@euroclone.it



www.euroclone.it