

Promotion

Corning® CoolCell® Container and Corning Cryogenic Vials



Corning® CoolCell® technology utilizes a thermo-conductive alloy core and highly insulated outer material to control the rate of heat removal and provide reproducible cell cryopreservation.

Unique features of controlled-rate freezing with Corning® CoolCell® include:

- ▶ Ease of use
- ▶ Alcohol and fluid-free freezing
- ▶ Lower total cost in use than alcohol-based devices
- ▶ High cell recovery and cell viability
- ▶ Reproducibility
- ▶ Simple, consistent way to standardize controlled-rate freezing
- ▶ Ready to use again in 5 minutes
- ▶ No lid stick due to frozen condensation

Corning Cat. No.	Description	Capacity (Vials)	Exposed Vial Tops	Qty/Cs	Promotional Price (€)	Compatible with 1.2 - 2 mL vials	Compatible with 4 - 5 mL vials
CC432001	CoolCell LX, purple	12	Yes	1	161,00	■	
CC432002	CoolCell LX, green	12	Yes	1	161,00	■	
CC432003	CoolCell LX, orange	12	Yes	1	161,00	■	
CC432004	CoolCell LX, pink	12	Yes	1	161,00	■	
CC432005	CoolCell 5 mL LX, purple	12	Yes	1	215,00		■
CC432006	CoolCell FTS30, purple	30	Yes	1	439,00	■	
CC432007	CoolCell FTS30, orange	30	Yes	1	439,00	■	
CC432008	CoolCell FTS30, green	30	Yes	1	439,00	■	
CC432009	CoolCell FTS30, pink	30	Yes	1	439,00	■	
CC432136	Cryogenic vial grippers			5	45,00		



Cryogenic vial grippers
Cat. Number CC432136



Corning Cryogenic Vials

- ▶ Larger marking spot
- ▶ Black graduations
- ▶ Certified RNase-/DNase-free and nonpyrogenic
- ▶ Vials have a silicone washer for a secure seal



External Thread Cryogenic Vials

Corning Cat. No.	Capacity (mL)	Style	Self-Standing	Qty/Pk	Qty/Cs	Promotional Price (€)
CC430658	1.2	Conical bottom	Yes	50	500	122,00
CC430659	2.0	Round bottom	Yes	50	500	122,00
CC430661	2.0	Round bottom	No	50	500	122,00
CC430662	4.0	Round bottom	Yes	50	500	192,00
CC430663	5.0	Round bottom	Yes	50	500	175,00



Internal Thread Cryogenic Vials

Corning Cat. No.	Capacity (mL)	Style	Self-Standing	Qty/Pk	Qty/Cs	Promotional Price (€)
CC430487	1.2	Conical bottom	Yes	50	500	122,00
CC430488	2.0	Round bottom	Yes	50	500	122,00
CC430489	2.0	Round bottom	No	50	500	122,00
CC430490	4.0	Round bottom	No	50	500	175,00
CC430492	5.0	Round bottom	No	50	500	220,00
CC430656	5.0	Round bottom	Yes	50	500	220,00
CC431386	2.0	Round bottom	Yes	50	250	96,00



Offer valid until July 31st,2018



Offer valid in Italy only and void where prohibited by government, law, or company policy. Offer valid for purchases till July 31,2018 or while supply lasts. Please allow 4-6 weeks for delivery.

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

Corning Incorporated
Life Sciences Europe
Corning BV
Fogstraat 12
1060 LJ Amsterdam
The Netherlands
Phone: +31 (0) 20 659 60 51
Fax: +31 (0) 20 659 76 73
CSEurope@corning.com
www.corning.com/lifesciences

CORNING | **FALCON**® **AXYGEN**® **GOSSSELIN**™

For additional product or technical information, visit www.corning.com/lifesciences, or contact our Scientific Support Team at ScientificSupportEMEA@corning.com.

For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners.

Distributed by:

EuroClone S.p.A.
Via Figino, 20/22
20016 Pero (MI)
Phone: +39 02 38.19.51
Fax: +39 02 38.10.14.65
Email: info@euroclone.it