



EuroClone®

serving science through innovation

Composition of MEM Vitamins 100X w/o L-Glutamine

Code: ECM0556

Inorganic salt :	
Sodium chloride	8500 mg/l
Vitamins :	
Choline chloride	100 mg/l
D-Ca pantothenate	100 mg/l
Folic acid	100 mg/l
Myo-inositol	200 mg/l
Nicotinamide	100 mg/l
Pyridoxal + HCl	100 mg/l
Riboflavin	10 mg/l
Thiamine + HCl	100 mg/l
Water	

EuroClone S.p.A.

Via Figino 20/22, 20016 Pero (MI) - Italy
Phone: +39.02.38195.1 - Fax +39.02.38101465
e-mail: info@euroclone.it



www.euroclone.it