

“I was looking
for efficiency,
quality and price.
I chose Euroclone”

Primo® Pipette
&
Primo® Tips

Primo® mechanical pipettes: high quality devices that guarantee maximum precision and reproducibility of measurement.

- Fully autoclavable
- UV resistant
- Ultra low pipetting forces
- Easy recalibration system
- 3 years warranty



Cat. N°	Volume	Inaccuracy %	Imprecision %
Single channel			
ECP10002	0.2 - 2 µl	± 12.0 - ± 1.5	± 6.0 - ± 0.7
ECP10010	0.5 - 10 µl	± 4.0 - ± 0.5	± 2.8 - ± 0.4
ECP10020	2 - 20 µl	± 3.0 - ± 0.8	± 1.5 - ± 0.3
ECP10100	10 - 100 µl	± 1.6 - ± 0.8	± 0.8 - ± 0.2
ECP10200	20 - 200 µl	± 1.2 - ± 0.6	± 0.6 - ± 0.2
ECP11000	100 - 1000 µl	± 1.6 - ± 0.6	± 0.4 - ± 0.15
8 channels			
ECP80010	0.5 - 10 µl	± 10.0 - ± 2.0	± 8.0 - ± 1.2
ECP80050	5 - 50 µl	± 4.0 - ± 1.6	± 2.5 - ± 0.6
ECP80200	20 - 200 µl	± 3.0 - ± 1.0	± 1.5 - ± 0.6
ECP80300	50 - 300 µl	± 1.6 - ± 1.0	± 1.5 - ± 0.6
12 channels			
ECP12010	0.5 - 10 µl	± 10.0 - ± 2.0	± 8.0 - ± 1.2
ECP12050	5 - 50 µl	± 4.0 - ± 1.6	± 2.5 - ± 0.6
ECP12200	20 - 200 µl	± 3.0 - ± 1.0	± 1.5 - ± 0.6
ECP12300	50 - 300 µl	± 1.6 - ± 1.0	± 1.5 - ± 0.6

Cat. N°	Description
ECP1KIT1	Kit Single Channel Primo® 1x (0.5 - 10 µl), 1x (10 - 100 µl), 1x (100 - 1000 µl)
ECP1KIT2	Kit Single Channel Primo® 1x (2 - 20 µl), 1x (20 - 200 µl), 1x (100 - 1000 µl)
ECP1KIT3	Kit Single Channel Primo® 2x (0.5 - 10 µl), 1x (10 - 100 µl)

Primo® filter tips:

- With high quality without cellulose additives, made to avoid cross- contamination
- Crystal clear quality
- Accurate graduation marks
- Low Retention properties
- Certified DNase, RNase, Human DNA, Pyrogen, PCR Inhibitors Free
- Sterile packaging
- Compatible with most pipettes available on the market



Code	Description	Format
ECTD00010	Primo® filter tips 0.1 - 10 µl, Sterile, Low Retention, racked	10x96 pcs
ECTD00011	Primo® filter tips 0.1 - 10 µl, Long, Sterile, Low Retention, racked	10x96 pcs
ECTD00020	Primo® filter tips 2 - 20 µl, Sterile, Low Retention, racked	10x96 pcs
ECTD00100	Primo® filter tips 2 - 100 µl, Sterile, Low Retention, racked	10x96 pcs
ECTD00200	Primo® filter tips 2 - 200 µl, Sterile, Low Retention, racked	10x96 pcs
ECTD00300	Primo® filter tips 2 - 300 µl, Sterile, Low Retention, racked	10x96 pcs
ECTD01005	Primo® filter tips 100 - 1000 µl, Sterile, Low Retention, racked	8x96 pcs

Primo® tips:

- Without filter
- Not sterile
- Autoclavable
- Available in: Bulk, Rack and Refill Kit

Code	Description	Format
Bag Not Sterile		
ECTD10010	Primo® tips 0.1 - 10 µl, clear, bag	1000 pcs
ECTD10011	Primo® tips 0.1 - 10 µl, long, clear, bag	1000 pcs
ECTD10200	Primo® tips 2 - 200 µl, clear, bag	1000 pcs
ECTD10300	Primo® tips 2 - 300 µl, clear, bag	1000 pcs
ECTD51000	Primo® tips 100 - 1000 µl, clear, bag	1000 pcs
Rack Not Sterile		
ECTD50010RN	Primo® tips 0.1 - 10 µl, Racked	10x96 pcs
ECTD50011RN	Primo® tips 0.1 - 10 µl, Long, Racked	10x96 pcs
ECTD50200RN	Primo® tips 2 - 200 µl, Clear, Racked	10x96 pcs
ECTD50300RN	Primo® tips 2 - 300 µl, Clear, Racked	10x96 pcs
ECTD51000RN	Primo® tips 100 - 1000 µl, Clear, Racked	8x96 pcs
Refill Not Sterile		
ECTD50010RL	Primo® tips 0.1 - 10 µl, Clear, Refill Kit	10x96 pcs
ECTD50011RL	Primo® tips 0.1 - 10 µl, Long, Clear, Refill Kit	10x96 pcs
ECTD50200RL	Primo® tips 2 - 200 µl, Clear, Refill Kit	10x96 pcs
ECTD51000RL	Primo® tips 100 - 1000 µl, Clear, Refill Kit	5x96 pcs

Primo® rack for Refill Kit:

Code	Description	Format
ECTD50010E	Primo® rack for 10 µl reload tips	50 pcs
ECTD50011E	Primo® rack for 10 µl long reload tips	50 pcs
ECTD50200E	Primo® rack for 200 µl reload tips	50 pcs
ECTD51000E	Primo® rack for 1000 µl reload tips	32 pcs

