

**NEW  
PRODUCT**



# TUBES, STRIPS, FOILS AND PLATES

## Introducing Our All-New Range of PCR and qPCR Consumables

### MAIN FEATURES OF PRIMO® TUBES, STRIPS AND PLATES

- Made from high-quality polypropylene and manufactured in an ISO 9001:2015 certified production facility.
- Free from DNA, RNA, DNases, RNase, pyrogens and PCR inhibitors to ensure accurate and consistent performance.
- Available in multiple skirt options.
- Compatible with the most used devices on in the market.



### ADDITIONAL FEATURES OF PRIMO® PLATES

- Wells with alphanumeric coding for easy sample tracking.
- Raised well rims for a perfect sealing.
- Frosted and white plates designed to enhance signal-to-noise ratios.

### MAIN FEATURES OF PRIMO® FOIL

- Adaptable to a broad range of temperatures from -40°C to 121°C without loss of integrity.
- Perfect for withstanding fluctuation of PCR temperature without peeling off.
- For colder storage in the range from -40°C to -80°C, Primo® Aluminium Foil is the best choice.
- Perfect for cold storage samples.

## ORDERING INFORMATION

CODE	DESCRIPTION	FORMAT
EPMT02SF	Primo® PCR tubes 0,2ml, with flat cap, bag, sterile	500
 EPMS0208F	Primo® PCR 8-strip 0,2 ml, <b>with single attached flat cap, bag</b>	12x10
EPMS0208O	Primo® PCR 8-strip 0,2 ml, with separate optical cap, bag	12x10
EPMT01LPF	Primo® PCR tubes 0,1 ml, with flat cap, low profile, bag	1000
EPMS0108LPO	Primo® PCR 8-strip 0,1 ml, with separate optical cap, low profile, bag	12x10
EPMF0500	Primo® Foil, PCR plate adhesive sealing	10x10
EPMF0560A	Primo® Foil, qPCR plate, optical adhesive sealing (similar to ABI)	10x10
EPMF0560B	Primo® Foil, qPCR plate, optical adhesive sealing (similar to Bio-Rad)	10x10
EPMF0550	Primo® Foil, PCR plate, aluminium adhesive sealing	10x10
EPMP0296NSK	Primo® 96 well plate non-skirted, 0,2ml, clear wells	5x10
EPMP0296SSK	Primo® 96 well plate semi skirted, 0,2ml, clear wells	5x10
EPMP0196NSK	Primo® 96 well plate non skirted, 0,1ml, clear wells	5x10
 EPMP0196WSK	Primo® 96 well plate semi skirted, 0,1ml, <b>white well</b>	5x10
EPMP0196FSSK	Primo® 96 well plate frosted, semi skirted, 0,1ml.	5x10
EPMP0196WFSK	Primo® 96 well plate fully skirted, 0,1ml, clear well	5x10

FOR MORE INFORMATION ABOUT THIS PRODUCT,  
**PLEASE CONTACT YOUR EUROCLONE REPRESENTATIVE**