

**NEW
PRODUCT**

PRIMO[®] *ecologik* PIPETTOR

The pipettor Primo[®]Eco-Logik, thanks to its unique features, is an essential tool for all laboratories.

The main features are: LED display that shows the charge status and set flow rate, adjustable aspiration and dispensing speed, electronic braking system (EBS[™]), and long-lasting rechargeable lithium battery, which allows for 8 hours of continuous operation on a full charge.

The pipettor is equipped with accessories that enhance its handness and can be used with plastic or glass pipettes from 0.1 to 100 ml.

Code	Description	Format
ECE2000	Primo [®] Eco-Logik Pipettor 0.1-100ml	1 pcs
ECE2000BAT	Primo [®] Eco-Logik Pipettor replacement lithium battery	1 pcs
ECE2000FT02	Primo [®] Eco-Logik Pipettor replacement filter 0.2um	1 pcs
ECE2000FT045	Primo [®] Eco-Logik Pipettor replacement filter 0.45um	1 pcs
ECE2000STAND	Primo [®] Eco-Logik Pipettor charging stand	1 pcs



FOR MORE INFORMATION ABOUT THIS PRODUCT,
PLEASE CONTACT YOUR EUROCLONE REPRESENTATIVE



<https://www.euroclonogroup.it/categories/1/13252/Manuali-digitali>