


**NEW
 PRODUCT**



**TRY RNAGarde TODAY
 AND ELEVATE YOUR RESEARCH!**
 Our Solution for Super RNA Stability

**TIRED OF STRUGGLING WITH
 RNA DEGRADATION IN YOUR
 EXPERIMENTS?**


No more worries about RNA integrity, with **RNAGarde** our stabilization solution that is designed to preserve your samples with high efficacy, maintaining their suitability for several downstream applications such as gene expression analysis, sequencing, and RT-qPCR.

**EFFECTIVE AT
 DIFFERENT TEMPERATURES**


TEMPERATURE	TIME
Frozen (-20°C / -80°C)	Indefinite
Refrigerated (-20°C / -80°C)	Up to 30 days
Room Temperature	Up to 7 days
Elevated (37°C)	24 hours

**EFFECTIVE WITH DIFFERENT
 BIOLOGICAL SAMPLES**


**ANIMAL AND
 PLANT TISSUES**




**CELL
 LINES**



**WHITE
 BLOOD CELLS**



**BACTERIA
 YEAST**



- **IMMEDIATE RNase INACTIVATION**

With a rapid tissue permeation offers maximum RNA stabilization and protection.

- **SIMPLIFIED STORAGE & HANDLING**

Eliminates the need for immediate freezing in liquid nitrogen or storing in lab freezers. Samples can be stored at room temperature, allowing for convenient sample collection at different time points without immediate processing.

- **SUITABLE FOR FIELD COLLECTION**

Specially tested for on-site tissue collection, perfect for locations like hospitals where immediate RNA processing is not feasible.

- **SAFE & NON-TOXIC**

An aqueous solution that poses no harm, allowing seamless downstream tissue processing.

- **VERSATILITY IN RNA ISOLATION & HIGH-QUALITY**

Compatible with a broad range of RNA isolation techniques, high quality RNA (mRNA or total)

READY TO IMPROVE YOUR RNA EXPERIMENTS?

Ensure the integrity, reliability, and reproducibility of your RNA samples with Euroclone **RNA*Garde*** solution for preservation and protection.



CODE - EMR605100
DESCRIPTION - RNA*Garde*
FORMAT - 100 ml

FOR MORE INFORMATION
ABOUT THIS PRODUCT,
**PLEASE CONTACT YOUR
EUROCLONE REPRESENTATIVE**



<https://www.euroclonogroup.it/categories/1/13252/Manuali-digitali>

Euroclone SpA Società a Socio Unico

Via Figino, 20/22 - 20016 Pero (MI) - +39 02.381951 - +39 02.38101465 - info@euroclone.it - www.euroclone.it

Quality Management System Certified to ISO 9001 and ISO 13485 international standard