



Euroclone has the exclusive distribution of Enzo products in Italy.

Today at a special price.



## Enzo Life Sciences


Enzo Life Sciences offers the broadest (and expanding) palette of IHC stain colors to choose from, giving you the ultimate flexibility in planning your **multiplex Immunohistochemistry protocol**.

- Widest panel of color choices
- Unique yellow, black and blue HRP chromogen
- Increased flexibility with HRP & AP options
- High definition color development
- Compatible with any automation system and with manual use

VALID UNTIL 30/06/2020

**20%**

PROMO CODE **PRPP20-1**



**DISCOUNT ON THE LARGEST COLLECTION  
OF HIGH DEFINITION CHROMOGENS FOR IHC AND ISH**

**Condition of sale:** Minimum order value 150.00 € | Cost contribution of 30.00 € for order less than 400.00 € | Cost contribution of 15.00 € for shipment with dry ice | Promo code PRPP20-1 | Offer not combinable with other discount and/or current offer.

## HIGH DEFINITION CHROMOGENS FOR IHC AND ISH

VALID UNTIL 30/06/2020

**20%**

PROMO CODE **PRPP20-1**

**DISCOUNT ON THE LARGEST  
COLLECTION OF HIGH DEFINITION  
CHROMOGENS FOR IHC AND ISH**

Codice	Descrizione	Format
<b>HRP Chromogens</b>		
EN9501510030	HIGHDEF® Blue IHC chromogen (HRP) #ADI-950-151-0030	30 mL
EN9501700030	HIGHDEF® Yellow IHC chromogen (HRP) #ADI-950-170-0030	30 mL
EN9501710030	HIGHDEF® Black IHC chromogen (HRP) #ADI-950-171-0030	30 mL
EN9502100030	HIGHDEF® Red IHC Chromogen (HRP) #ADI-950-210-0030	30 mL
ENACC1050200	HIGHDEF® DAB Chromogen/Substrate Set #ENZ-ACC105-0200	200 mL
EN9501820015	HIGHDEF® IHC enhancer (DAB, plus) #ADI-950-182-0015	15 mL
<b>AP Chromogens</b>		
EN9501400030	HIGHDEF® Red IHC chromogen (AP) #ADI-950-140-0030	30 mL
EN9501410030	HIGHDEF® Red IHC chromogen (AP, plus) #ADI-950-141-0030	30 mL
EN9501500030	HIGHDEF® Blue IHC chromogen (AP) #ADI-950-150-0030	30 mL
ENACC1300030	HIGHDEF® Green AP Chromogen/Substrate #ENZ-ACC130-0030	30 mL

## THE COMPLETE SOLUTION FOR YOUR IMMUNOHISTOCHEMISTRY

Codice	Descrizione	Format
<b>Antigen Retrieval Reagents</b>		
EN9502710500	IHC antigen retrieval reagent (citrate, pH 7.0) #ADI-950-271-0500	500 mL
EN9502710500	IHC antigen retrieval reagent (citrate, pH 7.0) #ADI-950-271-0500	500 mL
EN9502730500	IHC antigen retrieval reagent (EDTA, pH 8.0) #ADI-950-273-0500	500 mL
EN9502800015	IHC enzyme antigen retrieval reagent #ADI-950-280-0015	15 mL
ENACC1120100	Antigen Retrieval Reagent, pH 6 (10X) #ENZ-ACC112-0100	100 mL
ENACC1130100	Antigen Retrieval Reagent, pH 9 (10X) #ENZ-ACC113-0100	100 mL
<b>Detection*</b>		
ENACC1030150	POLYVIEW PLUS HRP (Anti-Rabbit) reagent #ENZ-ACC103-0150	150 tests
ENACC1040150	POLYVIEW PLUS HRP (Anti-Mouse) reagent #ENZ-ACC104-0150	150 tests
ENACC1100150	POLYVIEW PLUS AP (Anti-Rabbit) reagent #ENZ-ACC110-0150	150 tests
ENACC1140150	POLYVIEW PLUS AP (Anti-Mouse) reagent #ENZ-ACC114-0150	150 tests
ENKIT1590150	POLYVIEW PLUS HRP-DAB (Anti-Rabbit) Kit #ENZ-KIT159-0150	150 tests
ENKIT1600150	POLYVIEW PLUS HRP-DAB (Anti-Mouse) Kit #ENZ-KIT160-0150	150 tests
EN9501000001	MULTIVIEW (Mouse-HRP/Rabbit-AP) IHC kit #ADI-950-100-0001	100 tests
ENKIT1810150	MULTIVIEW PLUS (Mouse-HRP/Rabbit-AP) IHC Kit (Brown/Green) #ENZ-KIT181-0150	150 tests
ENACC1020150	SAVIEW PLUS HRP reagent #ENZ-ACC102-0150	150 tests
ENACC1110150	SAVIEW PLUS AP reagent #ENZ-ACC111-0150	150 tests
<b>Counterstain</b>		
ENACC1060100	HIGHDEF® Hematoxylin #ENZ-ACC106-0100	100 mL
<b>Mounting Medium</b>		
EN9502600025	HIGHDEF® IHC fluoromount #ADI-950-260-0025	25 mL
EN9502610030	HIGHDEF® IHC mount #ADI-950-261-0030	30 mL
<b>Controls</b>		
EN9500200025	MULTIVIEW Universal Tissue Microarray #ADI-950-020-0025	25 slides
EN9502310025	IHC universal negative control reagent #ADI-950-231-0025	25 mL
<b>Accessories</b>		
EN9502320001	IHC PAP (mini) pen #ADI-950-232-0001	1 each
EN9502330001	IHC PAP pen #ADI-950-233-0001	1 each
EN9502340100	IHC tissue primer #ADI-950-234-0100	100 mL
EN9502351000	IHC Wash Buffer (10X) #ADI-950-235-1000	1000 mL
EN9502400100	IHC PBS tablets #ADI-950-240-0100	100 each
EN9502440250	IHC diluent (primary Ab) #ADI-950-244-0250	250 mL
ENACC1090100	Streptavidin Blocker/Diluent #ENZ-ACC109-0100	100 mL
ENACC1080100	Antibody Blocker/Diluent #ENZ-ACC108-0100	100 mL

\*POLYVIEW® are Biotin-free, nanopolymer detection reagents

\*MULTIVIEW® IHC kits allow efficient detection of multiple antigens in a single tissue sample

PR 2038\_Promo Power\_Enzo\_0220

**Euroclone SpA Società a Socio Unico:**

Via Figino, 20/22 - 20016 Pero (MI) - +39 02.381951 - +39 02.38101465 - info@euroclone.it - www.euroclone.it

Quality Management System Certified to ISO 9001 and ISO 13485 international standard