



OFFERTA SPECIALE

PROMO


POWER

AFFRETTATI!

OFFERTA VALIDA
FINO AL 30/06/2019

CHIEDI AL TUO
AGENTE DI ZONA
L'OFFERTA A TE
DEDICATA PER LA
PROMO POWER!

CST:
EPILOPE TAG & LOADING CONTROLS ANTIBODIES



ANTICORPI

Per un periodo limitato di tempo una selezione di Anticorpi...ad un prezzo IRRIPETIBILE

Epitope Tag Antibodies

Gli epitopi Tag semplificano l'identificazione di proteine ricombinanti. Utilizzando anticorpi tag-specifici è possibile monitorare l'espressione di queste proteine e la visualizzazione in cellule o tessuti.

Loading Controls Antibodies

Questi anticorpi sono usati per confermare che la quantità di proteine caricate sia la stessa in tutte le lanes di un gel. I loading control sampler kit sono ideali per determinare qual è il controllo più adatto.

Epitope Tag Antibodies

CODICE	DESCRIZIONE	APPLICAZIONE
BK14682S	Digoxigenin (D8Q9J) Rabbit mAb	W
BK14681S	Dinitrophenol (D1D6) Rabbit mAb	W
BK13419S	E-Tag Antibody	W, IP,
BK8146S	DYKDDDDK Tag (9A3) Mouse mAb (Binds to same epitope as Sigma's Anti-FLAG® M2 Antibody)	IHC-P, W, F, IP, IF-IC
BK2368S	DYKDDDDK Tag Antibody (Binds to same epitope as Sigma's Anti-FLAG® M2 Antibody)	W, F, IP,
BK14793S	DYKDDDDK Tag (D6W5B) Rabbit mAb (Binds to same epitope as Sigma's Anti-FLAG® M2 Antibody)	IHC-P, W, F, IP, IF-IC, ChIP
BK2956S	GFP (D5.1) XP® Rabbit mAb	IHC-P, W, F, IF-IC
BK2955S	GFP (4B10) Mouse mAb	IHC-P, W, F,
BK2555S	GFP Antibody	IHC-P, W,
BK2448S	Glu-Glu Tag Antibody	W, IP, IF-IC
BK9325S	GSK-3β-Tag Antibody	W
BK2625S	GST (91G1) Rabbit mAb	W, IP,
BK2622S	GST Antibody	W, IP,
BK2624S	GST (26H1) Mouse mAb	W, IP, IF-IC
BK3724S	HA-Tag (C29F4) Rabbit mAb	IHC-P, W, F, IP, IF-IC, ChIP
BK2367S	HA-Tag (6E2) Mouse mAb	IHC-P, W, F, IF-IC
BK12698S	His-Tag (D311O) XP® Rabbit mAb	W, F, IP, IF-IC ChIP
BK2365S	His-Tag Antibody	W, IP
BK2366S	His-Tag (27E8) Mouse mAb	W, IP
BK2372S	β-Gal (14B7) Mouse mAb	W, F, IP, IF-IC
BK2396S	MBP Tag (8G1) Mouse mAb	W, IP, IF-IC
BK43590S	mCherry (E5D8F) Rabbit mAb	W, IF-IC, IF-F
BK2278S	Myc-Tag (71D10) Rabbit mAb	W, F, IP, IF-IC
BK2276S	Myc-Tag (9B11) Mouse mAb	IHC-P, W, F, IP, IF-IC, ChIP
BK2272S	Myc-Tag Antibody	W, F, IF-IC
BK92013S	Myc-Tag (E7F9B) Mouse mAb	W, IF-IC
BK12774S	S-Tag (D2K2V) XP® Rabbit mAb	W, F, IP, IF-IC
BK13246S	T7-Tag (D9E1X) XP® Rabbit mAb	W, IP, IF-IC
BK13202S	V5-Tag (D3H8Q) Rabbit mAb	W, F, IP, IF-IC
BK53061S	mNeonGreen Tag Antibody	W
BK68083S	HPC4-Tag Antibody	W, IP

Gli anticorpi sono disponibili anche coniugati a biglie, fluorofori, HRP e biotina. Sul sito cst-science.com/comp-tags la gamma completa di Epitope Tag Antibodies

Loading Controls Antibodies

CODICE	DESCRIZIONE	APPLICAZIONE
BK9774T	Loading Control Antibody Sampler Kit (Mouse)	
BK4670T	Loading Control Antibody Sampler Kit (HRP Conjugate)	
BK5142T	Loading Control Antibody Sampler Kit	W
BK12748S	Pan-Actin (D18C11) Rabbit mAb (HRP Conjugate)	W
BK4970T	β-Actin (13E5) Rabbit mAb	W, IHC-P, IHC-F, IF-IC, F
BK5125S	β-Actin (13E5) Rabbit mAb (HRP Conjugate)	W
BK8457S	β-Actin (D6A8) Rabbit mAb	W, IF-IC, F
BK12262S	β-Actin (8H10D10) Mouse mAb (HRP Conjugate)	W
BK3700S	β-Actin (8H10D10) Mouse mAb	W, IHC-P, IF-IC, F
BK4967S	β-Actin Antibody	W
BK5247S	COX IV (3E11) Rabbit mAb (HRP Conjugate)	W
BK4850S	COX IV (3E11) Rabbit mAb	W, IP, IHC-P, IHC-F, IF-IC, F
BK4844S	COX IV Antibody	W, IP, IHC-P, F
BK8884S	GAPDH (D16H11) XP® Rabbit mAb (HRP Conjugate)	W
BK5174S	GAPDH (D16H11) XP® Rabbit mAb	W, IHC-P, IF-IC
BK2118S	GAPDH (14C10) Rabbit mAb	W, IHC-P, IF-IC, F
BK5014S	GAPDH (14C10) Rabbit mAb (Biotinylated)	W, IF-IC, F
BK3683S	GAPDH (14C10) Rabbit mAb (HRP Conjugate)	W
BK4499S	Histone H3 (D1H2) XP® Rabbit mAb	W, IHC-P, IF-IC
BK5192S	Histone H3 (3H1) Rabbit mAb (HRP Conjugate)	W
BK9717S	Histone H3 (3H1) Rabbit mAb	W
BK3638S	Histone H3 (96C10) Mouse mAb	W
BK9715S	Histone H3 Antibody	W
BK5346S	β-Tubulin (9F3) Rabbit mAb (HRP Conjugate)	W
BK2128S	β-Tubulin (9F3) Rabbit mAb	W, IHC-P, IF-IC, F
BK2146S	β-Tubulin Antibody	W, IHC-P, IF-IC, F
BK2148S	α/β-Tubulin Antibody	W, IHC-P, IF-IC, F
BK9099S	α-Tubulin (11H10) Rabbit mAb (HRP Conjugate)	W
BK2125S	α-Tubulin (11H10) Rabbit mAb	W, IHC-P, IF-IC, F
BK3873S	α-Tubulin (DM1A) Mouse mAb	W, IHC-P, IF-IC
BK12351S	α-Tubulin (DM1A) Mouse mAb (HRP Conjugate)	W
BK2144S	α-Tubulin Antibody	W, IHC-P, IF-IC, F
BK13901S	Vinculin (E1E9V) XP® Rabbit mAb #13901	W, IHC, F

Sul sito cst-science.com/comp-loading-controls la gamma completa dei loading control antibodies disponibili.